
Invacare® Etude Plus® HC

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This manual **MUST** be given to the user of the product.
BEFORE using this product, read this manual and save for future reference.

Yes, you can.

© 2015 Invacare® Corporation

All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

This bed has been tested and approved according to IEC 60601-2-52

This bed has undergone a risk analysis according to EN ISO 14971.

Invacare is certified according to ISO 9001 and ISO 13485

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Contents

This manual **MUST** be given to the user of the product.
BEFORE using this product, read this manual and save for future reference.

I General	4
I.1 About this manual	4
I.2 Intended Use	4
I.3 Service life	4
I.4 Warranty	4
2 Safety	5
2.1 General warnings	5
2.2 Labels and symbols on the product	5
3 Setup	7
3.1 Receiving the bed	7
3.2 Assembly of the bed	7
3.3 Wiring	7
3.4 Dismantle the bed	8
4 Operating the bed	9
4.1 Hand control	9
4.2 Castors and brakes	9
4.2.1 Castor brake	9
4.3 Adjusting the leg section	9
4.4 Wooden side panel	10
4.5 Support handles / Swivel handles	10
4.6 Operating side rails	11
4.6.1 Steel side rails	11
4.6.2 Verso HC side rail	11
4.7 Lifting pole	12
4.7.1 Mounting the lifting pole	12
4.7.2 Changing the handle height	12
4.8 Emergency release of the backrest or thigh/leg section	12
4.9 Mattress support extension	13
5 Accessories	14
5.1 Accessories / Components	14
5.2 Mattresses	14
5.3 Transport brackets	14
6 Maintenance	16
6.1 Maintenance	16
6.1.1 Checklist maintenance	16
6.2 Inspection after relocation - Prepare for new user	16
6.2.1 Checklist - After relocation	16
6.3 Cleaning and disinfection	16
6.3.1 Cleaning methods	16
6.4 Lubrication	17
7 After Use	18
7.1 Waste disposal	18
8 Troubleshooting	19
8.1 Troubleshooting electrical system	19
9 Technical data	20
9.1 Dimensions	20
9.2 Mattress dimensions	20
9.3 Weights	20
9.4 Environmental conditions	20
9.5 Electrical system	20
9.6 Electromagnetic compliance (EMC)	21

I General

I.1 About this manual

This user manual contains important information about the handling of the product. In order to ensure safety when using the product, read the user manual carefully and follow the safety instructions.

All indications of right and left are based on a user lying on his back in the bed.

Please note that there may be sections in this user manual, which are not relevant to your product, since this manual applies to all existing modules (on the date of printing).

If a problem should arise in connection with the delivered product, please contact your *Invacare®* supplier. An address list is shown on the back side at this manual.

Symbols in this manual

In this user manual, warnings are indicated by symbols. The warning symbols are accompanied by a heading that indicates the severity of the danger.

WARNING

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

CAUTION

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.

IMPORTANT

Indicates a hazardous situation that could result in damage to property if it is not avoided.

Tips and Recommendations

Gives useful tips, recommendations and information for efficient, trouble-free use.

I.2 Intended Use

- The bed has been developed for domestic care and long term environments.
- The bed is intended for indoor use only.
- The bed has been developed for users over 12 years of age weighing 45-165 kg, taller than 150 cm.
- Safe Working Load (SWL): 200 kg (patient/user + accessories).
- Max User Weight: 165 kg.
- The bed is not intended for psychiatric patients.
- The bed is not intended for transportation of users. It is mobile within a room with a patient in it. Castors are lockable.

The bed is defined to be used according to:

- *Application environment 3*; Long-term care in a medical area where medical supervision is required and monitoring is provided if necessary and medical electrical equipment used in medical procedures may be provided to help maintain or improve the condition of the user.
- *Application environment 4*; Care provided in a domestic area where medical electrical equipment is used to alleviate or compensate for an injury, disability or disease.

This bed is height adjustable between 27 and 67 cm to enable safe use without side rails and to allow easier access in and out of bed.

CAUTION!

Risk of injury

Minimum working height of the bed is 27 cm.

Using an improper/incorrect posture while working at low working heights can lead to injuries for the carer.

- Be aware of your body posture when attending to the user.

WARNING!

Incorrect use could lead to hazardous situations.

Invacare accepts no liability for any use, change or assembly of the product, other than stated in this user manual.

I.3 Service life

The expected service life of this product is seven years when used daily and in accordance with the safety instructions, maintenance intervals and correct use, stated in this manual. The effective service life can vary according to frequency and intensity of use.

I.4 Warranty

The warranty covers all material and production defects for two years from the date of delivery, provided it can be demonstrated that such defects were present before delivery. All manufacturing faults or defects must be promptly reported.

Invacare® may repair the fault or replace the component. The warranty provided by *Invacare®* does not cover additional costs (transport, packaging, labour, sundry expenses, etc. are for the customer's account).

The warranty does not cover:

- Damage caused during transport that is not directly reported to the forwarder at the moment of delivery.
- Repairs performed by unauthorized centers and personnel.
- Parts subject to normal wear.
- Malicious damages or damages caused by improper use of the bed.

2 Safety

2.1 General warnings

WARNING!
Risk of injury

There is a risk of entrapment of fingers in the bed moving parts.

- Pay attention to your fingers.

WARNING!
Risk of entrapment / suffocation

There's a risk of entrapment / suffocation between mattress support, side rail and bed end or between moving parts and objects placed nearby the bed.

- The bed must not be used by persons under 12 years of age, or by persons with a body size equivalent to an average 12 years old or smaller.
- The bed, in combination with side rails must not be used for persons under 45 kg or under 150 cm height.
- The hand control must not be used by children.
- Due to mattress compression, an increased risk may occur over time. Periodically monitor gaps between the bed, mattress and/or side rail. Replace mattress if the gaps may lead to entrapment.
- Make sure that there is nothing under, over, or near the bed that can limit the movement of the bed or the mattress support; such as other medical devices (lifters), furniture, window frames and storage boxes.

WARNING!
Risk due to electromagnetic interference

Electromagnetic interference between the bed and other electrical products can occur.

- To reduce or eliminate such electromagnetic interference, increase the distance between the bed and the products or switch them off.

This medical bed can be used together with medical electrical equipment connected to the heart (intracardially) or blood vessels (intravascularly) provided that following points are respected:

- The bed should be equipped with means for potential equalization connection marked out by a symbol shown in the back of this manual.
- Medical electrical equipment should not be fixed on the bed's metallic accessories such as side rails, lifting pole, drip rod, bed ends, etc.
- The medical electrical equipment power supply cord should be kept clear of the accessories or any other moving parts of the bed.

WARNING!
Risk of slipping through the openings

The bed fulfils all requirements regarding maximum distances. However, it is possible that persons with small body dimensions slip through the openings between the side rails or through the opening between the side rail and the mattress support.

- Pay special attention, if the bed is used for the care of persons with small body dimensions.

CAUTION!

- Always lower the bed to the lowest position before leaving the user in the bed unattended.
- When entering or exiting the bed, the height adjustment should be used. Furthermore, the backrest can be used as support. However, the thigh- and leg section must be horizontal, otherwise, there is risk of overloading of the mattress support.

DANGER!
Risk of injury or damage to property

- Do not roll the castors over the main power cord.
- Do not bring main power cord into moving parts.
- Disconnect the plug from the mains before moving the bed.
- Make sure that no wires (mains or from other equipment) are jammed or damaged, when the bed is used.
- Keep bed components and accessories at least 30 cm away from a heated surface and not in direct sunlight.


DANGER!
Risk of fire

Use of wrong type of oxygen administering equipment can result in a fire hazard.

- Use only masked or nasal type oxygen administering equipment in conjunction with the electric bed.

2.2 Labels and symbols on the product

The product label is placed on the frame of the bed:

Product label

A	Type number of the bed
B	Serial number
C	Date of production
D	Manufacturer; Name and Address
E	Electrical information. See Chapter Technical Data.
F	Degree of protection. See also Technical data.
G	Symbols. See explanations in section; Symbols on the product.

Abbreviations

lin	Current in	Max	Maximum
Uin	Voltage in	min	minutes
Int.	Intermittency	kg	kilogram
AC	Alternating current	A	Ampere

Symbols

All symbols might not be present on all products.


	Max. patient weight
	Max. safe working load
	CLASS II equipment
	Type B Applied Part.

	Potential equalization
	Refer to the user manual.
	Waste disposal and recycling information.

3 Setup

3.1 Receiving the bed

When you receive the bed, check the packaging. If the bed shows any signs of damage upon delivery, please follow terms of delivery

Scope of delivery:

- Ⓐ Backrest (mattress support, upper half)
- Ⓑ Mattress retainers
- Ⓒ Hand control
- Ⓓ Thigh/Leg section (mattress support, lower half)
- Ⓔ Motor, bed end
- Ⓕ Castors
- Ⓖ Motor, backrest
- Ⓗ Wooden side panel
- Ⓘ Motor, thigh section
- ⓰ Bed end

Optional

- Support handle
- Swivel support
- Side rail
- Lifting pole
- Mattress support extension
- Back-up battery


3.2 Assembly of the bed


WARNING!

- After each assembly, check that all fittings are tightened correctly and that the function of all parts is given.

Mattress support

1. Place the two inserts in the upper half of the mattress support. The inserts must be mounted in such a way that one extends further out of the side tube than the other.
2. Loosely screw in the two thumb screws.
3. Push the lower half of the mattress support onto the two inserts.
4. Tighten with two thumb screws.
5. Retighten the two thumb screws at the head end of the top frame.
6. Turn up the four mattress support brackets.
7. Turn the handle of the leg section up.
8. Press the bracket Ⓐ on the locking ring* and turn the locking ring on the bed ends to the position "open".

9. Latch the mattress support to the bed ends and press firmly into position.
10. Make sure that the bed ends are in the same level. It is possible to level the bed in high position (1) and low position (2).
11. Turn the locking ring to the position "lock".

* The locking rings ensure that the mattress support is not accidentally lift clear off the bed ends.

Control Box

The control box is latched to the backrest motor. The control box is provided with a label with symbols showing where to connect the motor plugs:


- Backrest motor
- Thigh section's motor
- Bed end motor
- Hand control

3.3 Wiring

IMPORTANT!

- The cables must be mounted in such a way that they are kept clear off the floor and do not block the castors. Mount the mains cable on the hook for this purpose, see image.


In order to prevent the cables from being damaged when activating the motors, follow the instructions below.

1. Connect the cables of the bed end motors (head and foot end) directly to the control box.
2. The plug of the thigh motor must be directed through the opening of the supporter for the backrest motor and then connected to the control box.


3. Connect the main cable plug to the power socket.
4. Run the motors of the bed ends to their top position.
5. Place the motor cable of the foot end on the hook at the foot end.

6. Run the backrest to its top position.

7. Place the head end motor cables on the four hooks at the head end.

8. Attach the locking cam over the plugs in the control box.


3.4 Dismantle the bed

1. Dismantle side rails or wooden side panel and lifting pole.
2. Bring the bed to its lowest position and adjust all mattress support sections to horizontal position.
3. Disconnect main power supply. Roll the cable onto the hook of the head end of the bed.
4. Dismantle the locking cam from the control box with a tool, e.g. a flat headed screwdriver.
5. Disconnect all cables from the control box.
6. Separate the mattress support from the bed ends.
7. Divide the upper and lower halves of the mattress support.

4 Operating the bed

4.1 Hand control

-  The locking function is optional. If your control has no key  it is not possible to lock functions.

- ▲ Up
▼ Down

- Ⓐ Sitting position
Ⓑ Thight section
Ⓒ Height adjustment of the bed
Ⓓ Tilt function

- Ⓔ Key for locking and unlocking the electrical functions (optional)
Insert the key below each button to respectively lock and unlock the hand control functions.

1. Use the buttons  to  on the hand control to operate the electrical functions of the bed.
2. Insert the key  below each button and turn it to respectively lock and unlock the hand control functions (hand control unit with locking function only).

-  If the bed is used by restless or confused persons or persons with spasms, the bed should be equipped with a lockable hand control unit to block the hand control functions.

Tilt function

The bed can be equipped with a hand control (HL84) that activates the tilt of the mattress support. The tilt function is for user's comfort when sitting in the bed.

CAUTION! Risk of injury

- Only medically trained persons must operate the tilt function.
- The tilt function is NOT a Trendelenburger function and must not be used for medical treatment.
- When the tilt function is not needed, lock the tilt function with the key  on the hand control HL84.

4.2 Castors and brakes

CAUTION!

Risk of trapping/squeezing

- All brakes are foot-operated.
– Do not release the brake with the fingers.

CAUTION!

Risk of injury


- User can fall, while entering or exiting the bed, if brakes are not locked.
– Always lock the brakes before the user is moving in or out of bed or when attending to the user.
– At least one castor at the head end and one castor in the foot end must be locked.

Castors might make marks on different types of absorbing floor surfaces - such as untreated or badly treated floors. To prevent marks, *Invacare®* recommends to place a suitable kind of protection between the castors and the floor.


4.2.1 Castor brake

Locking the brake

Step on the outer pedal .


Unlocking the brake

Push up the release button .


4.3 Adjusting the leg section

Operate the leg section by lifting the mattress handle:

1. **Up:** Lift the mattress handle on the leg section.
2. **Down:** Lift the mattress handle on the leg section right up and then lower the leg section.

4.4 Wooden side panel

Press the brackets ① down, onto the top frame.


4.5 Support handles / Swivel handles

One or more support handles can be mounted on the top frame to support the user when getting in or out of bed.


WARNING! **Risk of serious injury**

The support handles are not to be used as side rails; to keep the user from falling out of the bed
– Only use the support handles as a support when getting in and out of bed.


Mounting the Support handles / Swivel handles

1.


- Loosen the finger screws ①.
- Rotate the metal plates ② out to the side.

2.

Press the brackets onto the top frame.


WARNING! **Risk of squeezing, getting trapped or suffocated**

Pay attention to the distance between support/swivel handle and bed end.


Pay attention to the distance between two support handles.

- The safety distance ① must be bigger than 390 mm.
- The safety distance ② must be either 0–60 mm or 318 mm or bigger.
- Never use the support handles in combination with side rails.

The Swivel handle has to be mounted with the movable part towards the center of the bed.

3.


- Rotate the metal plates ② back to close the brackets around the frame.
- Re-tighten the finger screws ①.

Operating the Swivel handles


Open and close the handle

1.


Pull the swivelling handle Ⓐ up to disengage it.

2.

Rotate the handle inwards or outwards to/from the bed.

3.

Push the handle down to secure it. Make sure it is properly locked before using it as support.

4.6 Operating side rails

WARNING!

Risk of entrapment or suffocation

There's a risk of entrapment or suffocation between mattress support, side rail and bed end.
– Always ensure correct fitting of the side rails.

4.6.1 Steel side rails


WARNING!

Risk of entrapment

There's a risk of entrapment or suffocation between mattress support, side rail and bed end.
– When using steel side rails, always make sure that the distance between the bed end and the side rail's handle upper edge is less than 6 cm in the head end and more than 32 cm in the foot end of the bed.

A sticker is placed on the side rail to indicate the minimum space between side rail and the head end of the bed:

CAUTION!

There's a risk of entrapment of fingers on assembly and operation of the side rail.
– Pay attention to your hands when using side rails.

Mounting, dismounting and operating

The side rail must be mounted according to the separate mounting instruction for the chosen side rail.

CAUTION!


If the side rail is not locked properly it can fall down.
– Pull/push the top bar of the side rail to ensure that the locking mechanism is properly engaged.

4.6.2 Verso HC side rail

Mounting the side rail


1. Unfold the side rail.

2.


Open the brackets by turning the finger screws Ⓐ anti clockwise (1). Press the fingers screws down and fold them backwards completely out of the notches (2).

3.

Press the U-formed side rail brackets onto the top frame.


4.

Fold the finger screws forwards (under the top frame) and up into the notches (1). Close the brackets by turning the finger screws clockwise (2).

5. Make sure the finger screws are tightened and the side rail properly attached (not loose).

Operating the side rail

Folding down the side rail

1. Hold the top bar of the side rail with one hand and push the release button (A) with the other hand.
2. Fold down the side rail sideways to its lowest position.

Folding up the side rail

1. Hold the top bar of the side rail and pull up until it locks with a clicking sound.

4.7 Lifting pole

WARNING! Risk of injury

The bed can tip if the handle is used, while the lifting pole is turned away from the bed.


- The lifting pole always has to be positioned with the handle hanging over the bed area.
- Do not exceed the maximum load of the lifting pole; 80 kg.

4.7.1 Mounting the lifting pole

1. Remove the plug from the lifting pole tube at the head end of the bed where the lifting pole is to be placed.
2. Insert the lifting pole into the lifting pole tube and fix with the thumb screw.


4.7.2 Changing the handle height

1.

Hold the cord out while sliding the plastic cord lock (A) upwards or downwards until the handle reaches the preferred height.

2.

Lock the handle by pressing back the cord in the lock (A) and pull the handle downwards.

IMPORTANT!

After adjustment of the lifting pole handle

- Check that the two cords above the cord lock, are parallel and inside the cord lock.
- Make sure the cord is properly locked by pulling the handle hard.

4.8 Emergency release of the backrest or thigh/leg section

An emergency release of the backrest or thigh/leg section could be necessary in the case of power- or motor failure.


CAUTION!

Risk of injury

- A minimum of two persons are required to release a mattress section.
- An emergency release of the height adjustment is NOT possible.
- Remove the plug from the mains socket before emergency release of the mattress support.

1. Both persons hold the mattress section.
- 2.

One of them pulls out the cotter pin (from the motor in question).

3. Both persons slowly lower the mattress section until it is completely down.


4.9 Mattress support extension

For users, taller than two meters it is recommended to use a mattress support extension.

1. Disassemble the bed and remove the inserts.
2. Mount the mattress support extension between both parts of the mattress support and reassemble the bed.

Fitting the mattress support extension

5 Accessories

5.1 Accessories / Components

For accessories, please see more information in the instructions delivered with the accessory.

Important

- Use only original accessories and spare parts. Spare parts lists and extra user manuals can be ordered from Invacare.

The tools necessary for mounting/dismounting the accessories of the bed are:

- Allen key
- adjustable spanner

5.2 Mattresses

WARNING!

Safety aspects regarding combination of side rails and mattresses:

To get the highest possible safety level, when using side rails on the bed, the minimum and maximum measures for mattresses, must be respected.

- For correct mattress measures see *mattress table: 9 Technical data, page 20*


WARNING!

Risk of entrapment and/or suffocation

- The user could get trapped and/or suffocate, if the horizontal space, between the mattress side and the inside of the side rail, is too big. Follow the minimum width (and length) of mattresses in combination with a given side rail, as stated in the *mattress table: 9 Technical data, page 20*

- Be aware that using very thick or soft mattresses (low density), or a combination of these, increases the risk.

WARNING!

Risk of falling


The user can fall over the edge and get seriously injured, if the vertical distance, between the top of the mattress and the edge of the side rail/bed end, is too short.

- Always keep a minimum distance of 22 cm. The maximum mattress height in combination with a given


side rail is stated in the *mattress table: 9 Technical data, page 20*

5.3 Transport brackets


1.


2.

3.

4.

6 Maintenance

6.1 Maintenance

A service contract can be made in the countries, where *Invacare®* has its own sales company. In certain countries *Invacare®* offers courses in service and maintenance of the bed. Spare parts lists and additional user manuals are available from *Invacare®*.

With normal daily use, service must be carried out according to the check list after two years use and thereafter every second year.

IMPORTANT!

- The mattress support must be supported during service inspections, to prevent accidental lowering.
- Only personnel who have received the necessary instructions or training may perform service and maintenance of the bed.
- After reconditioning the bed, or if bed functions change, service must be carried out according to the check list.

Before Use

- Ensure that all manual and electrical parts functions correctly and are in a secure state.

After three months

- Ensure that all manual and electrical parts are functioning, and tighten bolts, screws, nuts, etc.

Every year

- We recommend a safety test comprising the motors' performance and mechanical state.

Every second year

- Service should be performed according to the check list.

Motors, hand control and control units are serviced by exchanging the faulty part.

6.1.1 Checklist maintenance

Checkpoints

- ☐ Visual inspection of all parts of the bed (no deformations)
- ☐ Circlips, cotter pins and plastic fixing ring - properly locked and intact.
- ☐ Screws - tightened.
- ☐ Weldings - intact (no wear and tear)
- ☐ Castors (rolling smoothly)
- ☐ Castor fittings - tightened.
- ☐ Castor brakes - locking properly.
- ☐ All motors - running properly (with regular speed and at low noise).
- ☐ Wires - correctly wired and undamaged.
- ☐ Electric plugs - undamaged.
- ☐ Side rail locking and moving system - properly locking and running smoothly.
- ☐ Accessories - correctly assembled and correct function.
- ☐ Damaged coating - repaired.
- ☐ Leakage of grease - check for grease.

CAUTION!

Risk of injury or damage to property

Grease leaks can lead to accidental falls or fire.
– Contact service provider if leaks of grease or other residue is detected.

6.2 Inspection after relocation - Prepare for new user

IMPORTANT!

When the bed has been relocated; before given to a new user it has to be thoroughly inspected.
– Inspection must be done by a trained professional.
– For regular maintenance see maintenance chart.

6.2.1 Checklist - After relocation

Checkpoints

- ☐ Check that the inserts connecting the two mattress support halves are fully inserted and locked.
- ☐ Check that the mattress support is correctly mounted to the bed ends and that the locking ring is engaged.
- ☐ Check the locking of the motors (pipe pins correctly mounted).
- ☐ Check the electronic wiring for the motors (wires not squeezed).
- ☐ Check that the enclosing of the motors are intact (no cracks that enables fluids to pierce).
- ☐ Check that the locking cam is properly mounted over the contacts in the control box.
- ☐ Check the sections of the mattress support by using the hand control to activate all functions of the moving parts.
- ☐ Check the function of the brakes.
- ☐ Check the locking function(s) of the side rails.
- ☐ Check that all bed slats are intact and not loose.

6.3 Cleaning and disinfection

IMPORTANT!

Wrong fluids or methods can harm or damage you product.
– Follow instructions carefully for either Non-washable or Washable components.
– Never use corrosive fluids (alkalines, acid, cellulose thinner, acetone etc)
– Never use a solvent that changes the structure of the plastic or dissolves the attached labels.
– Always make sure that the bed is carefully dried before taking it into use again.

6.3.1 Cleaning methods

Electrical components

Method: Wipe off with a wet cloth or brush

Max temp: 40 °C

Solvent/chemicals: Water

Metal components

Method: Wipe off with a wet cloth or brush. Water may be pressurized, but not high pressure or steam.

Max temp: 40 °C

Solvent/chemicals: Household detergent or soap and water, 6-8 pH

Wood (and textile bands, if included in the design)

Method: Wipe off with a wet cloth or brush.

Max temp: 40 °C

Solvent/chemicals: Household detergent or soap and water, 6-8 pH

Textiles (including upholstery and mattresses)

Method: See attached label on each product.

6.4 Lubrication

Lubrication plan

We recommend lubricating the bed according to the following instructions:

- Points of rotation in mattress support and base frame - lubricate with medically clean oil.

- Motor attachment points to mattress platform - lubricate with medically clean oil.
- All of the motors' tension rod - lubricate with medically clean oil.
- All motor bearings - lubricate with medically clean oil.

The wooden side rails gliding system must not be lubricated with oil - otherwise the wooden bars will move sluggishly.

7 After Use

7.1 Waste disposal

Waste disposal/recycle must comply with the laws and regulations for waste handling in each country.

Invacare® is continuously working towards ensuring that the company's impact on environment, locally and globally, is reduced to a minimum.

We comply with the current environment legislation (e.g. WEEE and RoHS directive).

We only use REACH compliant materials and components.

- All wooden parts must be dismantled and sent for incineration.
- All electric parts must be dismantled and be disposed of as electric components.
- Plastic parts must be sent for incineration or recycling.
- Steel parts and castors must be disposed of as waste metals.

8 Troubleshooting

8.1 Troubleshooting electrical system

Symptom	Possible cause	Remedy
Mains indicator does not light up	Mains are not connected	Connect mains
	Fuse in the control unit is blown	* Replace the control unit
	Control unit is defective	* Replace the control unit
Mains indicator lights up, but the motor is not running. The relay in the control unit makes a clicking noise.	Motor plug is not fully inserted into the control unit.	Insert the motor plug properly into the control unit
	The motor is defective.	* Replace the motor
	Motor cable is damaged.	* Replace the cable
	Control unit is defective	* Replace the control unit
Mains indicator lights up, but the motor is not running. No relay sound is heard from the control unit.	Control unit is defective	* Replace the control unit
	Hand control is defective	* Replace the hand control
Control unit is in order except for one direction on one channel.	Control unit is defective	* Replace the control unit
	Hand control is defective	* Replace the hand control
Motor is running, but the piston rod does not move.	Motor is damaged	* Replace the motor
The motor cannot lift full load.		
Motor noise, but no movement of piston rod.		
Piston rod operates inwards and not outwards.		

* Service and maintenance of the bed must only be performed by personnel who have received the necessary instruction or training.

Risk of personal injury and damage to the product.

– The bed must be unplugged from the main power source before opening or repairing electrical parts.

WARNING!

Risk of personal injury and damage to the product.


– The bed must be unplugged from the main power source before opening or repairing electrical parts.

9 Technical data

9.1 Dimensions

All measurements are stated in cm. All angles are stated in degrees. All measurements and angles are stated without tolerances.


Invacare® reserves the right to change the stated measurements and angles.

	Etude Plus HC	
	Low	High
A	219	219
B	102	102
C	20-60	27-67
D	138-178	145-185

The angle between the lower leg section and horizontal is adjustable from 0° to 15°.

Bed with tilt function

9.2 Mattress dimensions

	Mattress measures (in cm):			
	Min. height	Max. height	Min. width	Min. length
Side rail: Verso HC	14	16	90	190

Min. density: 38kg/m³

9.3 Weights

	Max. patient weight (provided that the weight of the mattress and the accessories do not exceed 35 kg)	165 kg
	Max. safe working load (patient + accessories)	200 kg

Etude Plus standard bed, complete, excluding accessories	77.3 kg
Etude Plus auto-regression, complete, excluding accessories	82.7 kg
Etude Plus bed end with extra plate – 1 piece	20.6 kg

Etude Plus bed end without extra plate – 1 piece	17.9 kg
Mattress support, standard head end	24.1 kg
Mattress support with auto-regression, head end	27.9 kg
Mattress support with auto-regression, head end, without actuator and control box	24.1 kg
Mattress support, foot end	17.4 kg
Verso II side rail – 1 piece	7.4 kg
Lifting pole	4.2 kg
Mattress support extension (10 cm)	3.2 kg
Mattress support extension (20 cm)	5.5 kg
Support handle 40 x 30 cm	1,6 kg
Support handle 40 x 40 cm	2,0 kg
Swivel handle	1,8 kg

9.4 Environmental conditions

	Storage and transportation	Operation
Temperature	-10°C to +50°C	+5°C to +40°C
Relative humidity	20% to 75%	
Atmospheric pressure	800 hPa to 1060 hPa	

Be aware that when a bed has been stored under low temperatures, it must be adjusted to operating conditions before use.

9.5 Electrical system

Voltage supply: U _{in} = 120 Volt, AC, 50/60 Hz (AC = Alternating current)
Maximum current input: I _{in} max. = 1A
Intermittent (periodic motor operation): I _{nt} = max. 10%, 2 min ON / 18 min OFF
Degree of protection: IPx6* or IPx4** (depending on configuration)
The control unit, external power supply and motors are protected according to IPx6. When the bed is equipped with a HL80 hand control, it is protected according to IPx4
Insulation class: Class II
Type B Applied Part
Applied Part complying with the specified requirements for protection against electrical shock according to IEC60601-1.
(An applied parts is a part of the medical equipment which is designed to come into physical contact with the user or parts that are likely to be brought into contact with the user.)
Sound level: 45 – 50 dB (A)

* IPX6 classification means that the electrical system is protected against water projected in powerful jets from any direction.

** IPX4 classification means that the electrical system is protected against water splashed against the component from any direction.

The bed has no isolator (main switch). If the bed needs to be electrically disconnected, it has to be unplugged from the mains socket.

9.6 Electromagnetic compliance (EMC)

Guidance and manufacturer's declaration – electromagnetic emission

The medical bed is intended for use in the electromagnetic environment specified below. The customer or the user of the bed should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11(partly)	Group I	The medical bed uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11(partly)	Class B	The medical bed is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

The medical bed is intended for use in the electromagnetic environment specified below. The customer or the user of the bed should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (>95% dip in U_T) for 5 sec	< 5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the medical bed requires continued operation during power mains interruptions, it is recommended that the medical bed be powered from an un-interruptible power supply or a battery. U_T is the a. c. mains voltage prior to application of the test level.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Conducted RF IEC 61000-4-6	3 V	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the medical bed including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m	3 V/m	$d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). ^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the medical bed is used exceeds the applicable RF compliance level above, the medical bed should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the medical bed.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than $[V1]$ V/m.

At 80 MHz and 800 MHz, the higher frequency range applies.

Recommended separation distances between portable and mobile RF communications equipment and the medical bed

The medical bed is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the medical bed can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the medical bed as recommended below, according to the maximum output power of the communications equipment

Rated maximum output of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$	800 MHz to 2,5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above the recommended separation, distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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