

# INVACARE® MATTRESSES



Now you can offer your customers the comfort, protection and durability of superior **Homecare Mattresses** by Invacare, along with **Invacare® Mattress Covers**. Choose from innerspring or foam mattresses. Mattress covers will keep your Invacare® mattress clean and fresh, and protect it from stains, odors, dust and incontinence. Covers are made of non-allergenic waterproof vinyl and are available in zippered and contoured styles.

Invacare beds, mattresses and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit [www.invacare.com/bedsafety](http://www.invacare.com/bedsafety) for additional safety and entrapment information.

### Invacare® Innerspring Mattress

Model # 5185 , 5185XL

High quality .312 innerspring design

Durable 11-ounce vinyl cover is anti-bacterial, anti-static, acid-resistant and waterproof for easy use and care

5185XL (not shown) is available for use with extender kit

### Invacare® Economy Foam Mattress

Model # 5180 , 5184

High quality 1.5 lb. density foam provides comfort and durability in a one-piece construction

### Invacare® Bariatric Foam Mattress

Model # BARMATT(39/42/48/EXP)

Constructed of latex-free materials

Durable 11-ounce vinyl cover is anti-bacterial and waterproof for easy use and care

Solid inner core

Mattress is reversible and flippable

Constructed of high resilience foam for extended life and service for the homecare provider

Easy to clean

HCPCS Code: E0271

HCPCS Code: E01272



5185



BARMATT39

1. Invacare® MC0195-1 zippered mattress cover comes individually wrapped; 12 per case
2. Invacare® MC0065-1 contoured mattress cover comes individually wrapped; 12 per case

**Invacare Corporation**  
www.invacare.com

#### USA


One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

#### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2012

Invacare Corporation.  
All rights reserved.  
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.  
Form No. 11-219  
111225



Model #	Dimensions LxWxD	Gauge	Regulation Compliance	Limited Warranty
MC0165-1, MC0195-1	80"x36"x6"	4"	N/A	30 Days
5180	80"x36"x5.5"	N/A	16CFR1632/33	1 Year
5184	84"x36"x5.5"	N/A	16CFR1632/33	1 Year
5185	80"x16"x6"	N/A	16CFR1632/33	1 Year
5185XL	84"x36"x6"	N/A	16CFR1632/33	1 Year
BARMATT39	80"x39"x6.5"	N/A	16CFR1632/33	1 Year
BARMATT42	80"x42"x6.5"	N/A	16CFR1632/33	1 Year
BARMATT48	80"x48"x6.5"	N/A	16CFR1632/33	1 Year
BARMATTEXP	80"x39"-48"x6.5"	N/A	16CFR1632/33	1 Year

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

