

INVACARE® MATRX®

MINI BACK



The **Invacare Matrix Mini Back** features Mini Set hardware, a specially designed miniature version of popular the Easy Set hardware. Every element of this back support has been carefully crafted and constructed, keeping pediatric users at the forefront of the design. The Mini Back is contoured to fit children of all shapes and sizes, combining exceptional support in an appealing low profile package.

Call to order

1.800.333.6900
invacare.com

1.888.433.6818
motionconcepts.com



matrx
SEATING SERIES

Features and Functional Benefits

Cover – Foam and cover overlap the shell for comfort and protection. Water-resistant and breathable. Available in a choice of Pink Lollipop, Blizzard Blue or Kid*ab*ra print.

Hardware – “Mini Set” hardware installs and adjusts in minutes. Small attachment point for minimal interference on back canes. Easily adjusts for growth – depth (4”), angle (40°) and even rotation. Adjustable mounting pins allow fit to chairs same or 1” wider than back. Quick release for convenient portability.

Back Shell – Lightweight, durable aluminum shell, designed for easy attachment of optional lateral supports and headrest.

Adjustable Lumbo-Sacral Support – Included.

Headrest Mount– The Mini Back Shell is pre-drilled to accept Motion Concepts, Stealth and Otto Bock™ headrest mounts.

Swing-Away Lateral Supports – 2.5” of vertical adjustment on back pan. Mounts inside back shell without interfering with depth/angle adjustment. Easy push-button release. Available in black, blue & pink.

Lateral pads – XS: 4.25” H x 4.25” L

HCPCS Code: E2615



1. Swing-Away
2. 3” Contour Depth



Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com
info@motionconcepts.com

USA
700 Ensminger Rd., Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905)-695-0134
Fax: 905-695-0138

© 2016
Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Matrix is a registered trademark of Motion Concepts. Otto Bock is a trademark of Otto Bock HealthCare. Specifications are subject to change without notification.
Form No. 15-034
Rev. 09/16
160959



Model #	Height	Width	Depth* Adjustment	Angle Adjustment	Product Weight	Weight Capacity	Limited Warranty
MINI	7", 10", 12" & 14"	10" -14"	4"	40° range +/- 20°	0.9 lb. and up	150 lb.	2 Year

*Limitation of depth/angle adjustment if used on chairs narrower than model width.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program and should consult an attorney or other advisor to discuss specific situations in further detail.

