

INVACARE® MATRX®  
**ELITE™ TR**  
BACK



FREEDOM DESIGNS® P.R.O. CG™ Tilt-In-Space Wheelchair with Elite™ TR Back\*

The **Invacare Matrx Elite TR Back** is specially designed for use in manual or power tilt-only or tilt/recline applications. It includes an extra layer of SuperSoft HR foam for improved immersion and comfort, and incorporates full support across the shoulder area. Bariatric Elite TR Back available with 500 lb. capacity.

## Call to order

1.800.333.6900  
invacare.com

1.888.433.6818  
motionconcepts.com



**matrx**  
SEATING SERIES

## Features and Functional Benefits

**Cover** – Oversized cover featuring SuperSoft HR foam layered over firm HR foam for comfort and protection. Attractive Meshtex cover is water-resistant and breathable.

**Hardware** – “Easy Set” hardware installs and adjusts in minutes – all adjustments can be made with user seated in the wheelchair. Other hardware styles may be substituted for unique mounting applications; Split, Fixed Space Saver, Extended Range, as well as many versions specific to power wheelchair applications including ESR Mount.

**Back Shell** – Stylish, durable aluminum shell includes contoured support in shoulder area, slotted for lateral support hardware attachment.

**Adjustable Lumbo-Sacral Support** – Included.

**Headrest Adapter** – Required for easy, durable and safe mounting of headrest hardware and to provide offset for clearance of foam overlap.

## Positioning Options

**Lateral Supports** – Full range of fixed, offset and swing-away laterals available for order separately.

**Fixed** – 2.5" of vertical adjustment on back pan. Mounts inside back shell without interfering with depth/angle adjustment. Available in standard 5" H x 7.5" L and pediatric 3" H x 5.5" L sizes.

**Offset Fixed** – Provides an additional 3" of vertical position (higher or lower) relative to mounting slots. Available 5" H x 7.5" L.

**Swing-Away** – 2.5" of vertical adjustment on back pan. Mounts inside back shell without interfering with depth/angle adjustment. Easy push-button release, four pad sizes (see below).

**Offset Swing-Away** – Provides an additional 3" of vertical position (higher or lower). Mounts inside back shell without interfering with depth/angle adjustment. Easy push-button release, four pad sizes (see below).

**Lateral Pads for Swing-Away Laterals** – XS: 3.25" H x 4.25" L, S: 3.5" H x 5.25" L, M: 3.75" H x 5.75" L, L: 4.25" H x 6.5" L

## HCPCS Code: E2620



1. Fixed Lateral Pads
2. Offset Fixed Lateral Pads
3. Swing-Away Laterals
4. Offset, Swing-Away Laterals
5. 3" Contour Depth
6. Privacy Flap



**Invacare Corporation**  
www.invacare.com

**USA**  
One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

**Canada**  
570 Matheson Blvd. E.  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

**Motion Concepts**  
www.motionconcepts.com  
info@motionconcepts.com

**USA**  
700 Ensinger Rd., Suite 112  
Tonawanda, NY 14150  
(888) 433-6818  
Fax: (888) 433-6834

**Canada**  
84 Citation Dr., Unit 1-7  
Concord, ON L4K 3C1  
(866) 748-7943  
(905) 695-0134  
Fax: 905-695-0138

Model #	Height	Width	Depth** Adjustment	Angle Adjustment	Product Weight	Weight Capacity	Limited Warranty
PBE TR	16", 18", 20"	15" - 20"	4"	40° range +/- 20°	3.5 lb. and up	300 lb.	2 Year
PBE TRHD	20"	21" - 24"	4"	40° range +/- 20°	8.5 lb. and up	500 lb.	2 Year

\*FREEDOM DESIGNS® P.R.O. CG™ Tilt-In-Space Wheelchair weight capacity limit is 250 lbs.

\*\*Limitation of depth/angle adjustment if used on chairs narrower than model width

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program and should consult an attorney or other advisor to discuss specific situations in further detail.



© 2016  
Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols ™ and  
©. All trademarks are  
owned by or licensed to  
Invacare Corporation unless  
otherwise noted. Matrix  
is a registered trademark  
of Motion Concepts.  
Specifications are subject  
to change without  
notification.  
Form No. 15-035  
Rev. 09/16  
160958