



Yes, you can.®

INVACARE®

# AQUATEC® COMMODES

OCEAN ERGO / OCEAN ERGO SP

---

Ocean Ergo



Ocean Ergo SP



---

Family of shower chair commodes have been designed with user comfort, independence and dignity in mind. Based on the well-established Aquatec® Ocean range, the Ocean Ergo family comes packed with a range of new features designed to make personal care activities easier for both the user and caregiver

A new and more ergonomic design encourages an upright and stable seating position. Ergo shower commode chairs are manufactured from stainless steel and are rust proof. The improved caster design eliminates the collection of debris in the axle area. The ergonomically designed collection pan secures close to the seat ensuring peace of mind for the user and caregiver. The Ocean Ergo components are easily removed for cleaning in an autoclave or the entire shower commode can be used in a chair washer. The Ocean Ergo SP comes complete with 24" rear wheels with stainless steel quick release axles and wheel locks promotes independence and privacy.

## FEATURES

- Ergonomic seat plate with a 5° posterior angle to support natural toileting position
- Tool-free seat height adjustment range of 5" to fit over most toilets Latex free
- Enhanced frame design for greater over toilet positioning and forward stability
- Four 5" locking caster with one directional lock to enable straight tracking with ease
- Ergonomically shaped flip up armrest flip up for transfer
- 24" Rear Wheel Kit with Quick Release Axles and Wheel Lock standard on Ergo SP
- XL models with 22" backrest and 398lb weight capacity available in Standard and SP

## ACCESSORIES



1. Soft seat with Hygiene recess
2. Universal soft seat - Ocean Ergo Series
3. Variable soft seat - Ocean/VIP/Ergo Series
4. Hygiene Backrest Support
5. Splash guard - Ocean Ergo Series

## SPECIFICATIONS

MODEL	ERGO	ERGO SP	ERGO XL	ERGO SP XL
Length	35.3"	41"	35"	41"
Width	22"	26"	25"	25"
Height	40.5" - 45.5"	40.5" - 45.5"	40.5" - 45.5"	40.5" - 45.5"
Product Weight	37.5 lbs	50.7 lbs	39.7 lbs	50.7 lbs
Width Between Arms/Handles	18"	18"	21.75"	21.75"
Seat Width	19"	19"	19"	19"
Back Width	19"	19"	22"	22"
Seat Depth	19"	19"	19"	19"
Seat Angle	5°	5°	5°	5°
Seat to Footrest	11.5" - 22.5"	11.5" - 22.5"	11.5" - 22.5"	11.5" - 22.5"
Height Seat to Floor	18.5" - 23"	18.5" - 23"	18.5" - 23"	18.5" - 23"
Capacity	330 lbs	330 lbs	398 lbs	398 lbs

## WARRANTY

Aquatec® OceanErgo Series shower commodes are covered under a two-year limited warranty against manufacturing defects. **Manufacturer's warranty does not cover normal wear.**

All Aquatec® OceanErgo shower commode accessories are covered for a period of one year.