



NAME \_\_\_\_\_ PHONE NO \_\_\_\_\_

ACCOUNT # \_\_\_\_\_ P.O. \_\_\_\_\_ DATE \_\_\_\_\_

ADDRESS \_\_\_\_\_

<b>Pre-Configured Stock Models</b>		
<input type="checkbox"/> 9SL_PTO_34745	9SL/WD66/88/BH1517	\$1279
9000SL 16" x 16" Adult Frame with Fixed Height Space - Saver Desk Arm		
<input type="checkbox"/> 9SL_PTO_34749	9SL/WD66/78/BH1517	\$1279
9000SL 18" x 16" Adult Frame with Fixed Height Space - Saver Full Arm		
<input type="checkbox"/> 9SL_PTO_34751	9SL/WD66/88/BH1517	\$1279
9000SL 18" x 16" Adult Frame with Fixed Height Space - Saver Desk Arm		
<input type="checkbox"/> 9SL_PTO_34752	9SL/WD66/89/BH1517	\$1329
9000SL 18" x 16" Adult Frame with Adjustable Height Space - Saver Desk Arm		
<input type="checkbox"/> 9SL_PTO_47874	9SL/WD06/88/BH1517/34P/ADULT	\$1393
9000SL 20" x 16" Adult Frame with Fixed Height Space - Saver Desk Arm		
<b>Stock Front Rigging</b>		
<input type="checkbox"/> T93HAP	Hemi Footrests, Aluminum Footplates with heel loops	\$159
<input type="checkbox"/> T93HCP	Hemi Footrests, Composite Footplates with heel loops	\$159
<input type="checkbox"/> T94HAP	Hemi Elevating Legrests; Alum. Footplates & Padded Calf Pads	\$309
<input type="checkbox"/> T94HCP	Hemi Elevating Legrests; Comp. Footplates & Padded Calf Pads	\$309
<input type="checkbox"/> AHL4A	Hemi Smart Artic. Legrests; Alum. plates & Padded Calf Pads	\$399
<input type="checkbox"/> T93HEP	Economy Hemi Footrests; Comp. Footplates without Heel Loops	\$127
<input type="checkbox"/> T94HEP	Economy Hemi Elev. Legrests; Comp. Footplates & Calf Pads	\$268
<b>Stock Options</b>		
<input type="checkbox"/> 1497	Oxygen Holder Pkg Assembly	\$198
<input type="checkbox"/> 1496	IV/Oxygen Holder Pkg Assembly	\$198
<input type="checkbox"/> 1360	Adj. Rear Anti-Tippers (8" Casters)	\$109
<input type="checkbox"/> 1360A	Adj. Rear Anti-Tippers (6" Casters)	\$123
<input type="checkbox"/> A1576	Wheel Lock Extensions	\$63
<i>Elastic cord wheel lock extension. Hemi STF's come standard with non-cord wheel lock extension.</i>		

CHAIR FEATURES AND OPTIONS
<ul style="list-style-type: none"> <li>• Chrome steel frame</li> <li>• Steel crossbraces</li> <li>• Dual axle position</li> <li>• Urethane 24" rear wheels and 8" casters</li> <li>• Seat-to-Floor height; Adult at 19½" and Hemi at 17½" with 24" rear wheels and 8" casters</li> <li>• Sectional height back with 10° bend, 15" - 19" in 1" increments</li> <li>• Space-saver flip back non-removable arms</li> <li>• Adult frame chair with non-removable flip back desk length arms weighs less than 34 lbs.</li> <li>• Lifetime warranty on frame and crossbraces</li> </ul>

WEIGHT LIMITATION 250 LBS.
<p>Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.</p>

## FRAME & SEAT

### FRAME TYPE

<input checked="" type="checkbox"/> 9SL	9000SL	\$1,215
---	--------	---------

### FRAME WIDTH & DEPTH

<input type="checkbox"/> WD66	16" X 16"	No Charge
<input type="checkbox"/> WD68	16" X 18"	No Charge
<input type="checkbox"/> WD86	18" X 16"	No Charge
<input type="checkbox"/> WD88	18" X 18"	No Charge
<input type="checkbox"/> WD06	20" X 16"	No Charge
<input type="checkbox"/> WD08	20" X 18"	No Charge

### SEAT TO FLOOR HEIGHT

*Based on 24" wheel and 8" caster, see seat-to-floor matrix. See last page for Matrices.*

<input type="checkbox"/> Adult	(19.5" STF)	No Charge
<input type="checkbox"/> Hemi	(17.5" STF)	No Charge

## ARM & BACK

### ARM TYPE

*Fixed height arm 10", adj height arm 10"-14"*

<input type="checkbox"/> 19	Full Length Adj.Height Conventional(10-14")	\$372
<i>Flip back conventional removable arm style. Available with Adj. Angle Back only.</i>		
<input type="checkbox"/> 29	Desk Length Adj.Height Conventional(10-14")	\$372
<i>Flip back conventional removable arm style. Available with Adj. Angle Back only.</i>		
<input type="checkbox"/> 78	Full Length Fixed Height Space-Saver(10")	\$261
<i>Flip back non-removable arm style. Not available with adj.angle back type.</i>		
<input type="checkbox"/> 79	Full Length Adj.Height Space-Saver(10-14")	\$372
<i>Flip back non-removable arm style. Not available with adj.angle back type.</i>		
<input type="checkbox"/> 88	Desk Length Fixed Height Space-Saver(10")	\$261
<i>Flip back non-removable arm style. Not available with adj.angle back type.</i>		
<input type="checkbox"/> 89	Desk Length Adj.Height Space-Saver(10-14")	\$358
<i>Flip back non-removable arm style. Not available with adj.angle back type.</i>		

### BACK TYPE

<input type="checkbox"/> ANGLE	Adj Angle Back w/ 10° Bend	\$172
<i>Not available with space-saver arms. Standard with anti-tippers. Adjustable height 15-19", factory set at 16".</i>		

### BACK HEIGHT

*Available with Adj. Angle Back only. BH1517 is factory set at 16" High, BH1719 is factory set at 18" high*

<input type="checkbox"/> BH1517	Adjustable Back Height 15"-17"	<b>Standard</b>
<input type="checkbox"/> BH1719	Adjustable Back Height 17"-19"	No Charge

**COLOURS****FRAME FINISH**

<input type="checkbox"/> 08P	Deep Purple	No Charge
<input type="checkbox"/> 24P	Wet Black	No Charge
<input type="checkbox"/> 34P	Satin (Flat) Black	<b>Standard</b>
<input type="checkbox"/> 71P	Silver Vein	No Charge

**UPHOLSTERY COLOURS**

<input type="checkbox"/> U240	Black Nylon Upholstery	<b>Standard</b>
<input type="checkbox"/> U67	Black Vinyl	No Charge

**CASTERS****CASTER OPTIONS**

<input type="checkbox"/> AT902	6" Composite/Urethane 3 Spoke	\$120
<input type="checkbox"/> AT903	8" X 1-1/4" Composite/Urethane	<b>Standard</b>
<input type="checkbox"/> 1299	8" X 1" Solid Rubber Casters	No Charge
<input type="checkbox"/> 1292	8" X 1-1/4" Pneumatic Casters	No Charge
<input type="checkbox"/> 1492	8" X 1-1/4" Pneu w/Flat Free Insert	\$128

**WHEEL OPTIONS**

<input type="checkbox"/> U1763BK	20" Composite/Urethane Tires	No Charge
<input type="checkbox"/> U1764BK	22" Composite/Urethane Tires	No Charge
<input type="checkbox"/> 2244BK	24" Composite/Pneumatic Tires	\$231
<input type="checkbox"/> 4244BK	24" Composite/Pneumatic-Flat Free Inserts	\$353
<input type="checkbox"/> U2244BK	24" Composite/Urethane Treaded Tires	No Charge
<input type="checkbox"/> U2222C	24" Composite/Urethane Corded Tires	<b>Standard</b>

**WHEEL LOCKS**

<input type="checkbox"/> 9695	Hill Holder	\$149
<i>Not available with urethane rear tires.</i>		
<input type="checkbox"/> 9690B	Pull to Lock Wheel Locks (Standard on Hemi)	No Charge
<input type="checkbox"/> 9691B	Push to Lock Wheel Locks (Standard on Adult)	No Charge

**HANDRIM OPTIONS**

<input type="checkbox"/> AHDRM	Aluminum Handrims	\$43
<i>Not available with standard U2222 rear wheels.</i>		
<input type="checkbox"/> COM	Composite Handrim	Standard
<i>Only Available with 24" Rear Wheels.</i>		
<input type="checkbox"/> L8702	Omit Handrims	(Less Net \$15)

**POSITIONING STRAPS****SEAT POSITIONING STRAPS**

<input type="checkbox"/> 1311BK	Airline Buckle Seat Positioning Straps	\$84
<input type="checkbox"/> 1515	Auto Style Buckle Seat Positioning Straps	\$84
<input type="checkbox"/> 1321BK	Hook & Loop Seat Positioning Straps	\$66

**ACCESSORIES****BACK ACCESSORIES**

<input type="checkbox"/> L8601	Omit Back Upholstery	(Less net \$20)
<input type="checkbox"/> 1160	Carrying Pocket	\$83
<input type="checkbox"/> STC	Stencil Back Upholstery	\$59

**SEAT ACCESSORIES**

<input type="checkbox"/> L8651	Omit Seat Upholstery	(Less Net \$15)
<input type="checkbox"/> 9120	Solid Seat Insert	\$144

**ARM ACCESSORIES**

<input type="checkbox"/> NRAS	Non-Removable Arm with Screws	No Charge
<input type="checkbox"/> L8405	Omit Arm Pads	(Less Net \$15)

**MISCELLANEOUS ACCESSORIES**

<input type="checkbox"/> 9357L	Amputee Attachment	\$131
<input type="checkbox"/> 1577QR	Quick Release Axles	\$142
<input type="checkbox"/> 2218	IV Rod, Telescopic	\$185
<input type="checkbox"/> 2226	Front Anti-Tippers	\$158

**STOCK OPTIONS**

<input type="checkbox"/> 1497	Oxygen Holder Package Assembly	\$198
<input type="checkbox"/> 1496	IV/Oxygen Holder Package Assembly	\$198
<input type="checkbox"/> 1360	Adjustable Rear Anti-Tipper (8" Casters)	\$109
<i>Standard with angle adjustable back.</i>		
<input type="checkbox"/> 1360A	Adjustable Rear Anti-Tipper (6" Casters)	\$123
<input type="checkbox"/> A1576	Wheel Lock Extensions (8")	\$63

*Elastic cord wheel lock extension. Hemi STF's come standard with non-cord wheel lock extension.*

**FRONT RIGGING OPTIONS**

**FRONT RIGGING**

- AHL4A Hemi Smartleg Artic. Legrest 16" wide Alum. Ftpts \$399
- AHL4C Hemi Smartleg Artic. Legrest 16" wide Comp. Ftpts \$392
- T93HA Hemi Swingaway Footrest Alum. Ftpts \$159
- T93HC Hemi Swingaway Footrests Comp. Ftpts \$159
- T94HA Hemi Sw Elevating Legrests Alum. Ftpts \$309
- T94HC Hemi Sw Elevating Legrests Comp. Ftpts \$302
- T93HAB Hemi Bolt-Thru Footrest, Alum. Ftpts \$158
- T93HCB Hemi Bolt-Thru Footrest, Comp. Ftpts \$159
- T94HAB Hemi Bolt-Thru Legrest, Alum. Ftpts \$302
- T94HCB Hemi Bolt-Thru Legrest, Comp. Ftpts \$309
- T93HE Hemi Footrest w/o Loop Composite Ftpts \$119
- T94HE Hemi Legrest Comp. w/ Blk Non-Removable Calfpad \$239

**FOOTPLATE OPTIONS**

- 1350 Extra Large Aluminum Footplates \$106

**FOOTREST / LEGREST ACCESSORIES**

- 1820 Non-Removable Device \$41
- Not available with AHL4A or AHL4C Articulating legrests*
- 1346 Heel Loops \$63
- Heel Loops are standard on Swingaway Footrests.*
- L8515 Omit Calf Pads (Less net \$10)
- L8689 Omit Heel Loops (Less net \$10)

Seat-to-Floor Matrix 6" Casters			
Seat-to-Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
15.5"	Top	20"	Top
16.5"	Middle	22"	Top
17.5"	Bottom	24"	Top
17.5"	Bottom	20"	Bottom

14 Not available with AHL4A or AHL4C Articulating legrests

Seat-To-Floor Matrix 8" Casters			
Seat-to-Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
17.5"	Top	20"	Bottom
17.5"	Top	24"	Top
18.5"	Middle	22"	Bottom
19.5"	Bottom	24"	Bottom

**Invacare Canada reserves the right to change prices and specifications without notice.**

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324 (Canada) Fax: (905) 890-5244 (800) 668-5478 (Canada)  
 Invacare®, and Yes, you can., are trademarks of Invacare Corporation.  
 WWW.INVACARE.CA