



NAME \_\_\_\_\_ PHONE NO \_\_\_\_\_

ACCOUNT# \_\_\_\_\_ P.O. \_\_\_\_\_ DATE \_\_\_\_\_

ADDRESS \_\_\_\_\_

**Pre-Configured Stock Models**

**Stock Chairs**

<input type="checkbox"/> 9XT_PTO_29153	9XT/WD86/88 - 9000XT 18" x 16" Adult, Desk Length Fixed Ht Space Saver Arm	\$1649
<input type="checkbox"/> 9XT_PTO_29154	9XT/WD66/88 - 9000XT 18" x 16" Adult, Desk Length Fixed Ht Space Saver Arm	\$1649
<input type="checkbox"/> 9XT_PTO_29155	9XT/WD86/78 9000XT 18" x 16" Adult, Full Length Fixed Ht Space Saver Arm	\$1649
<input type="checkbox"/> 9XT_PTO_33131	9XT/WD06/88 9000XT 20" x 16" Adult, Desk Length Fixed Ht Space Saver Arm	\$1649

**Stock Front Rigging**

<input type="checkbox"/> T93HAP	Hemi Footrests, Aluminum Footplates with heel loops	\$159
<input type="checkbox"/> T93HCP	Hemi Footrests, Composite Footplates with heel loops	\$159
<input type="checkbox"/> T94HAP	Hemi Elevating Legrests; Alum. Footplates & Padded Calf Pads	\$309
<input type="checkbox"/> T94HCP	Hemi Elevating Legrests; Comp. Footplates & Padded Calf Pads	\$309
<input type="checkbox"/> AHL4A	Hemi Smart Artic. Legrests; Alum. Footplates & Padded Calf Pads	\$399
<input type="checkbox"/> T93HEP	Economy Hemi Footrests; Comp. Footplates without Heel Loops	\$127
<input type="checkbox"/> T94HEP	Economy Hemi Elev. Legrests; Comp. Footplates & Calf Pads	\$239

**Stock Options**

<input type="checkbox"/> 1497	Oxygen Holder Pkg Assembly	\$198
<input type="checkbox"/> 1496	IV/Oxygen Holder Pkg Assembly	\$198
<input type="checkbox"/> 1360	Adj. Rear Anti-Tippers (8" Casters)	\$109
<input type="checkbox"/> 1360A	Adj. Rear Anti-Tippers (6" Casters)	\$123
<input type="checkbox"/> A1576	Wheel Lock Extensions	\$63

*Elastic cord wheel lock extension.*

**FRAME & SEAT**

**FRAME TYPE**

<input type="checkbox"/> 9XT	9000XT	\$1,809
------------------------------	--------	---------

**FRAME WIDTH & DEPTH**

<input type="checkbox"/> WD46	14" X 16"	No Charge
<input type="checkbox"/> WD48	14" X 18"	No Charge
<input type="checkbox"/> WD56	15" X 16"	No Charge
<input type="checkbox"/> WD58	15" X 18"	No Charge
<input type="checkbox"/> WD66	16" X 16"	No Charge
<input type="checkbox"/> WD67	16" X 17"	No Charge
<input type="checkbox"/> WD68	16" X 18"	No Charge
<input type="checkbox"/> WD76	17" X 16"	No Charge
<input type="checkbox"/> WD77	17" X 17"	No Charge
<input type="checkbox"/> WD78	17" X 18"	No Charge
<input type="checkbox"/> WD86	18" X 16"	No Charge
<input type="checkbox"/> WD87	18" X 17"	No Charge
<input type="checkbox"/> WD88	18" X 18"	No Charge
1 <input type="checkbox"/> WD96	19" X 16"	No Charge
1 <input type="checkbox"/> WD97	19" X 17"	No Charge
1 <input type="checkbox"/> WD98	19" X 18"	No Charge
1 <input type="checkbox"/> WD06	20" X 16"	No Charge
1 <input type="checkbox"/> WD07	20" X 17"	No Charge
1 <input type="checkbox"/> WD08	20" X 18"	No Charge
1,2 <input type="checkbox"/> WD26	22" X 16"	\$206
1,2 <input type="checkbox"/> WD27	22" X 17"	\$330
1,2 <input type="checkbox"/> WD28	22" X 18"	\$330

**SEAT TO FLOOR HEIGHT**

3 <input type="checkbox"/> Adult	Adult (19-½" STF)	No Charge
3 <input type="checkbox"/> H	Hemi (17-½" STF)	No Charge
4.5 <input type="checkbox"/> SHE	Super Hemi (15-½" STF)	\$63
6 <input type="checkbox"/> SH	Super Low (14 ½" STF)	No Charge

**WEIGHT LIMITATION 250 LBS.**

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional knowledgeable of the particular limitations of the product.

**FOOTNOTES**

- 1 Not available with Super Low seat-to-floor.
- 2 Not available with adjustable angle back type.
- 3 Based on 24" wheels and 8" casters. See seat-to-floor matrix.
- 4 Based on 20" wheels and 6" casters. See seat-to-floor matrix.
- 5 Requires of AT902 (6" casters), U1763BK (20" composite wheels, urethane tires) AHDRM (aluminum handrims) and A1576 (wheel lock extensions). Seat-to-floor height is approx. 15"
- 6 Not available on ADJC and ADJCPB, upcharge will apply if required.

## ARM & BACK

### 7 ARM TYPE

<input type="checkbox"/> 18	Full Length Fixed Height Conventional(10")	\$251
<input type="checkbox"/> 19	Full Length Adj.Height Conventional(10-14")	\$358
<input type="checkbox"/> 28	Desk Length Fixed Height Conventional(10")	\$251
<input type="checkbox"/> 29	Desk Length Adj.Height Conventional(10-14")	\$358
2 <input type="checkbox"/> 78	Full Length Fixed Height Space-Saver(10")	\$251
2 <input type="checkbox"/> 79	Full Length Adj.Height Space-Saver(10-14")	\$358
2 <input type="checkbox"/> 88	Desk Length Fixed Height Space-Saver(10")	\$251
2 <input type="checkbox"/> 89	Desk Length Adj.Height Space-Saver(10-14")	\$358
2,8 <input type="checkbox"/> 79JY	Full Length Adjustable Ht Space-Saver(8-10")	\$358
2,8 <input type="checkbox"/> 89JY	Desk Length Adjustable Ht Space-Saver(8-10")	\$358

### BACK TYPE

<input type="checkbox"/> ADJ	Adjustable Back w/10° Bend (15"-19")	<b>Standard</b>
<input type="checkbox"/> S1418	Straight Back w/Push Handles(14"-18")	No Charge
1,9,10 <input type="checkbox"/> ANGLE	Adj Angle Back w/10° Bend	\$165
1,9,10 <input type="checkbox"/> ANGLE-S	Adj Angle Back-Straight	\$165

## COLOURS

### FRAME FINISH

<input type="checkbox"/> 08P	Deep Purple	No Charge
<input type="checkbox"/> 17P	Cobalt Blue	No Charge
<input type="checkbox"/> 24P	Wet Black	No Charge
<input type="checkbox"/> 34P	Satin (Flat) Black	<b>Standard</b>
<input type="checkbox"/> 71P	Silver Vein	No Charge
<input type="checkbox"/> 88P	Wet Black w/Twilight Sparkle	No Charge
<input type="checkbox"/> 105P	Midnight Blue	No Charge

### CANDY COLOURS

<input type="checkbox"/> 61P	Electric Red	No Charge
<input type="checkbox"/> 62P	Electric Blue	No Charge
<input type="checkbox"/> 66P	Electric Purple	No Charge
<input type="checkbox"/> 104P	Electric Teal	No Charge

### UPHOLSTERY COLOURS

<input type="checkbox"/> U240	Black Nylon Upholstery	No Charge
12 <input type="checkbox"/> U245	Black Sail Cloth Adjust. Tension Seat	<b>Standard</b>
13 <input type="checkbox"/> U67	BlackVinyl	No Charge

### NAUGAHYDE

13 <input type="checkbox"/> WRAP	Wrap Around Vinyl Upholstery	No Charge
13,14 <input type="checkbox"/> SOUPH	Screwed On Back Upholstery	\$140

## CASTERS

### CASTER OPTIONS

<input type="checkbox"/> AT902	6" Composite/Urethane 3 Spoke	No Charge
1 <input type="checkbox"/> 1228	8" X 1 3/4" Semi-Pneumatic Casters Lt. Grey	\$150
1 <input type="checkbox"/> AT903	8" x 1 1/4" Composite/Urethane	<b>Standard</b>
1 <input type="checkbox"/> 1299	8" X 1" Solid Rubber Casters	No Charge
1 <input type="checkbox"/> 1292	8" X 1-1/4" Pneumatic Casters	No Charge
1 <input type="checkbox"/> 1492	8" X 1-1/4" Pneu w/Flat Free Insert	\$128

### 1,15 CASTER JOURNAL OPTION

16 <input type="checkbox"/> ADJC	Adjustable Caster Journal
16 <input type="checkbox"/> ADJCPB	Adj. Precision Bearing Caster Journal

### FOOTNOTES

- 7 Removable arm style, flipback.
- 8 Only Available with 20" Rear Wheels
- 9 Not available with adjustable tension back upholstery.
- 10 Standard with rear anti-tippers.
- 11 Only available with Super Low seat-to-floor.
- 12 U245 sail cloth adjustable tension seat upholstery available w/U240 black nylon back upholstery only.
- 13 Vinyl upholstery must choose WRAP or SOUPH.
- 14 Not available with Adjustable Angle or Adjustable Angle-S back type.

### See Additional Charts on Page 6

- \$63 Seat to Floor matrix w/6" Casters & Offset Axle Bracket
- No Charge Seat to Floor matrix w/8" Casters & Offset Axle Bracket

**WHEELS**

**WHEEL OPTIONS**

<input type="checkbox"/>	1763BK	20" Composite/Pneumatic Tires	\$231	
<input type="checkbox"/>	4763BK	20" Composite/Pneumatic-Flat Free Inserts	\$349	
<input type="checkbox"/>	U1763BK	20" Composite/Urethane Tires	No Charge	
<input type="checkbox"/>	1764BK	22" Composite/Pneumatic Tires	\$231	
<input type="checkbox"/>	4764BK	22" Composite/Pneumatic-Flat Free Inserts	\$349	
<input type="checkbox"/>	U1764BK	22" Composite/Urethane Tires	No Charge	
1	<input type="checkbox"/>	2244BK	24" Composite/Pneumatic Tires	<b>Standard</b>
1	<input type="checkbox"/>	4244BK	24" Composite/Pneumatic-Flat Free Inserts	\$128
1	<input type="checkbox"/>	5244BK	24" Heavy Duty Comp/Pneu.-Flat Free Inserts	\$420
1,17	<input type="checkbox"/>	U2244BK	24" Composite/Urethane Treaded Tires	No Charge
1,18	<input type="checkbox"/>	U2222C	24" Composite/Urethane Corderd Tires	No Charge

**19 WHEEL OPTIONS**

**ONE ARM DRIVE (OAD)**

20,17	<input type="checkbox"/>	9456L	One Arm Drive w/U2222C LH	\$672
20,17	<input type="checkbox"/>	9456R	One Arm Drive w/U2222C RH	\$672

**CLD ONE ARM DRIVE**

1,19,20,21	<input type="checkbox"/>	LDR96RI	CLD One Arm Drive Right, 6" Caster	\$1,349
1,19,20,21	<input type="checkbox"/>	LDR96LI	CLD One Arm Drive Left, 6" Caster	\$1,349
1,19,20,21	<input type="checkbox"/>	LDR98RI	CLD One Arm Drive Right, 8" Caster	\$1,349
1,19,20,21	<input type="checkbox"/>	LDR98LI	CLD One Arm Drive Left, 8" Caster	\$1,349

**HANDRIM OPTIONS**

22	<input type="checkbox"/>	AHDM	Aluminum Handrims	<b>Standard</b>
23	<input type="checkbox"/>	COM	Composite Handrims	No Charge
	<input type="checkbox"/>	1512	Plastic Coated Aluminum Handrims	\$114
	<input type="checkbox"/>	BHDM	Black Handrims (Alum for U2222)	No Charge
17,23	<input type="checkbox"/>	1224-8	8 Alum Oblique Proj,hdrm,plstic Coat	\$307
17,23	<input type="checkbox"/>	1225-8	8 Alum Vert Proj,hdrm,plstic Coat	\$307

**REAR AXLE**

1,16	<input type="checkbox"/>	OFFAXL	Offset Axle	No Charge
------	--------------------------	--------	-------------	-----------

**16 WHEEL LOCKS**

24	<input type="checkbox"/>	9695	Hill Holder	\$149
	<input type="checkbox"/>	A1567	Pull To Lock Wheel Lock-Clamp On	No Charge
	<input type="checkbox"/>	1721	High Push To Lock Wheel Lock-Clamp On	No Charge
	<input type="checkbox"/>	9690B	Pull to Lock Wheel Lock (Bolt Thru)(STD on Hemi)	No Charge
	<input type="checkbox"/>	9691B	Push to Lock Wheel Lock (Bolt Thru)(STD on Adult)	No Charge
	<input type="checkbox"/>	9690	Pull To Lock Wheel Locks (Disk Style)	No Charge
	<input type="checkbox"/>	9691	Push To Lock Wheel Locks (Disk Style)	No Charge

**POSITIONING STRAPS**

**SEAT POSITIONING STRAPS**

<input type="checkbox"/>	1311BK	Airline Buckle Seat Positioning Strap	\$84
<input type="checkbox"/>	1515	Auto Style Buckle Seat Positioning Strap	\$84
<input type="checkbox"/>	1321BK	Hook and Loop Seat Positioning Strap	\$66

**CHEST POSITIONING STRAPS**

<input type="checkbox"/>	7311BK	Airline Buckle Chest Positioning Strap	\$84
<input type="checkbox"/>	9515	Auto Style Chest Positioning Strap	\$84
<input type="checkbox"/>	7321BK	Hook and Loop Chest Positioning Strap	\$66

**FOOTNOTES**

- 15 Fixed caster journals are standard.
- 16 Not available with Adjustable Angle or Adjustable Angle-S back type.
- 17 High profile treaded urethane tires.
- 18 Low profile non-treaded Urethane
- 19 The Cyclical Lever Drive requires treaded rear wheels and is only available on 16, 17, and 18" deep. Hemi STF must use Space Saver arms. Weight limit is250 lbs.
- 20 Not available w/amputee attachment or offset axle
- 21 Not available with adjustable caster journal..
- 22 Not available w/U2222 rear wheels.
- 23 Not available on 20" and 22" composite rear wheels.
- 24 Not available with urethane tires or 20" rear wheels.
- 25 Only available w/Vinyl upholstery. 3/4" block letters.

**ACCESSORIES & STOCK OPTIONS**

**BACK ACCESSORIES**

20	<input type="checkbox"/>	ADJTENS	Adjustable Tension Back Upholstery	\$142
	<input type="checkbox"/>	1160	Carrying Pocket	\$83
25	<input type="checkbox"/>	STC	Stencil Back Upholstery	\$59
	<input type="checkbox"/>	RB	Dual Back Liner	\$91
	<input type="checkbox"/>	L8601	Omit Back Upholstery	Less net (20.00)
	<input type="checkbox"/>	BP02	Back Pack	\$105

**SEAT ACCESSORIES**

26	<input type="checkbox"/> 9120	Solid Seat Insert	\$138
	<input type="checkbox"/> RS	Dual Seat Liner	\$87
	<input type="checkbox"/> L8651	Omit Seat Upholstery	Less net (15.00)

**ARM ACCESSORIES**

	<input type="checkbox"/> NRAS	Non-Removable Arm With Screws	No Charge
27	<input type="checkbox"/> 1300	Utility Tray	\$268

**MISCELLANEOUS ACCESSORIES**

28,29	<input type="checkbox"/> 9357L	Amputee Attachment	\$131
30	<input type="checkbox"/> 1305L	Crutch & Cane Carrier	\$66
28,29	<input type="checkbox"/> 1577QR	Quick Release Axles	\$142
30	<input type="checkbox"/> 2218	IV Rod, Telescopic	\$185
	<input type="checkbox"/> 2226	Front Anti-Tippers	\$158
11	<input type="checkbox"/> 9357SH	Super Low Amputee Attachment	\$131
11	<input type="checkbox"/> 1360SH	Super Low Adjustable Rear Anti-Tipper	\$109

**STOCK OPTIONS**

28,30,31	<input type="checkbox"/> 1497	Oxygen Holder Package Assembly	\$198
28,30,32	<input type="checkbox"/> 1496	IV/Oxygen Holder Package Assembly	\$198
29,30,33,34	<input type="checkbox"/> 1360	Adjustable Rear Anti-Tipper (8" Casters)	\$109
29,30,35	<input type="checkbox"/> 1360A	Adjustable Rear Anti-Tipper (6" Casters)	\$123
36	<input type="checkbox"/> A1576	Wheel Lock Extensions	\$63

**FRONT RIGGING OPTIONS**

**FRONT RIGGING**

11,37	<input type="checkbox"/> SH904M	Super Low Elevating Legrest 11"-17"	\$356
11,37	<input type="checkbox"/> SH904MS	Super Low Elevating Legrests 14"wide	\$356
11	<input type="checkbox"/> SH93E	Super Low Ext Sw Footrests 11"-17"	\$230
11	<input type="checkbox"/> SH93ES	Super Low Ext Sw Footrests 14"wide	\$230
37	<input type="checkbox"/> AHL4A	Hemi Smartleg Artic. Legrests 16" wide Alum. Ftpts.	\$399
37	<input type="checkbox"/> AHL4C	Hemi Smartleg Artic. Legrests 16" wide Comp. Ftpts.	\$399
37	<input type="checkbox"/> AHS4A	Hemi Smartleg Articulating Legrests 14" wide	\$399
	<input type="checkbox"/> T93HA	Hemi Sw Footrests Alum. Ftpts.	\$159
	<input type="checkbox"/> T93HC	Hemi Sw Footrests Comp. Ftpts.	\$159
37	<input type="checkbox"/> T94HS	Hemi Sw Elevating Legrests Alum Ftpt for 14" wide	\$302
37	<input type="checkbox"/> T94HA	Hemi Sw Elevating Legrests Alum. Ftpts.	\$309
37	<input type="checkbox"/> T94HC	Hemi Sw Elevating Legrests Comp. Ftpts.	\$302
	<input type="checkbox"/> T93HS	Hemi Sw Footrests Alum. Ftpts. 14"wide	\$158
	<input type="checkbox"/> T963A	60 Deg Footrest Aluminum Footplates	\$158
	<input type="checkbox"/> T963C	60 Deg Footrest Composite Footplates	\$158
	<input type="checkbox"/> T963S	60 Deg Footrest Composite Footplate 14"wide	\$158
	<input type="checkbox"/> T93HAB	Hemi Bolt-Thru Footrests, Aluminum Ftpts	\$158
	<input type="checkbox"/> T93HCB	Hemi Bolt-Thru Footrests, Composite Ftpts	\$158
37	<input type="checkbox"/> T94HAB	Hemi Bolt-Thru Legrests, Aluminum Ftpts	\$302
37	<input type="checkbox"/> T94HCB	Hemi Bolt-Thru Legrests, Composite Ftpts	\$302

**FOOTPLATE OPTIONS**

	<input type="checkbox"/> 1350	Extra Large Aluminum Footplates	\$106
	<input type="checkbox"/> AT5543	Adjustable Angle Flip Up Footplates	\$132

**FOOTREST / LEGREST ACCESSORIES**

	<input type="checkbox"/> 1820	Non-Removable Device	\$41
	<input type="checkbox"/> 1337	Calf Strap	\$63
	<input type="checkbox"/> IMPACTG	Impact Guards	No Charge
39	<input type="checkbox"/> 1346	Heel Loops	\$63

**FOOTNOTES**

- 26 Not available on 22" wide chairs.
- 27 Not available with 24" rear wheel set in Hemi position.
- 28 Not available with offset axle.
- 29 Not available with any one arm drive.
- 30 Not available with traditional one arm drive, but is available with CLD one arm drive.
- 31 For 0° and 6° seat angles. See service manual for details.
- 32 Not available with amputee attachment.
- 34 Standard w/Adjustable Angle Back.
- 35 Only available with 6" casters.
- 36 Elastic cord wheel lock extensions. Hemi seat to floors come standard with non-cord wheel lock extensions.
- 37 Elevating legrests available in chrome or black finish.
- 38 Only available on T944 elevating legrests.
- 39 Standard on footrests.
- 40 Only available w/H93M front riggings.

Seat-To-Floor w 6" Casters and Offset Axle Bracket					
Seat-To Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position	Axle Plate Mounting On Frame	
				Top Lug	Bottom Lug
15.5"	Top Hole	22"	#5 Hole	#4 Hole	#1 Hole
15.5"	Top Hole	20"	#3 Hole	#4 Hole	#1 Hole
16.5"	Middle Hole	24"	#5 Hole	#4 Hole	#1 Hole
16.5"	Middle Hole	22"	#3 Hole	#4 Hole	#1 Hole
16.5"	Middle Hole	20"	#2 Hole	#4 Hole	#1 Hole
17.5"	Bottom Hole	24"	#3 Hole	#4 Hole	#1 Hole
17.5"	Bottom Hole	22"	#2 Hole	#4 Hole	#1 Hole
17.5"	Bottom Hole	20"	#3 Hole	#4 Hole	#1 Hole

For a complete list of seat-to-floor heights and angles that are configurable by the dealer refer to the service manual


Seat-To-Floor Matrix w 8" Casters and offset Axle Bracket-					
Seat-To Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position	Axle Plate Mounting On Frame	
				Top Lug	Bottom Lug
17.5"	Top Hole	20"	# 3 Hole	#1 Hole	#4 Hole
17.5"	Top Hole	22"	# 2 Hole	#4 Hole	#1 Hole
17.5"	Top Hole	24"	# 3 Hole	#4 Hole	#1 Hole
18.5"	Middle Hole	22"	# 3 Hole	#1 Hole	#4 Hole
18.5"	Middle Hole	24"	# 2 Hole	#4 Hole	#1 Hole
19.5"	Bottom Hole	20"	# 3 Hole	#1 Hole	#4 Hole

For a complete of seat-to-floor heights and angles that are configurable by the dealer refer to the service manual.

SEAT-TO-FLOOR MATRIX 6" CASTERS			
Seat to Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
15 1/2"	Top	20"	Top
16 1/2"	Middle	22"	Top
17 1/2"	Bottom	24"	Top
17 1/2"	Bottom	20"	Bottom

SEAT-TO-FLOOR MATRIX 8" CASTERS			
Seat to Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
17 1/2"	Top	20"	Bottom
17 1/2"	Top	24"	Top
18 1/2"	Middle	22"	Bottom
19 1/2"	Bottom	24"	Bottom

SEAT-TO-FLOOR MATRIX 6" CASTERS & SUPER LOW FRAME			
Seat to Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
13 1/2"	Top	20"	Only One
14 1/2"	Bottom	22"	Only One



Invacare Canada reserves the right to change prices and specifications without notice.  
 HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324 (Canada) Fax: (905) 890-5244 (800) 668-5478 (Canada)  
 Invacare®, and Yes, you can., are trademarks of Invacare Corporation.  
 WWW.INVACARE.CA