



NAME \_\_\_\_\_ PHONE NO \_\_\_\_\_

ACCOUNT# \_\_\_\_\_ P.O. \_\_\_\_\_ DATE \_\_\_\_\_

ADDRESS \_\_\_\_\_

**Pre-Configured Stock Models**

**Stock Chairs**

<input type="checkbox"/> 9XT_PTO_29153	9XT/WD86/88 - 9000XT 18" x 16" Adult, Desk Length Fixed Ht Space Saver Arm	\$1649
<input type="checkbox"/> 9XT_PTO_29154	9XT/WD66/88 - 9000XT 18" x 16" Adult, Desk Length Fixed Ht Space Saver Arm	\$1649
<input type="checkbox"/> 9XT_PTO_29155	9XT/WD86/78 9000XT 18" x 16" Adult, Full Length Fixed Ht Space Saver Arm	\$1649
<input type="checkbox"/> 9XT_PTO_33131	9XT/WD06/88 9000XT 20" x 16" Adult, Desk Length Fixed Ht Space Saver Arm	\$1649

**Stock Front Rigging**

<input type="checkbox"/> T93HAP	Hemi Footrests, Aluminum Footplates with heel loops	\$159
<input type="checkbox"/> T93HCP	Hemi Footrests, Composite Footplates with heel loops	\$159
<input type="checkbox"/> T94HAP	Hemi Elevating Legrests; Alum. Footplates & Padded Calf Pads	\$309
<input type="checkbox"/> T94HCP	Hemi Elevating Legrests; Comp. Footplates & Padded Calf Pads	\$309
<input type="checkbox"/> AHL4A	Hemi Smart Artic. Legrests; Alum. Footplates & Padded Calf Pads	\$399
<input type="checkbox"/> T93HEP	Economy Hemi Footrests; Comp. Footplates without Heel Loops	\$127
<input type="checkbox"/> T94HEP	Economy Hemi Elev. Legrests; Comp. Footplates & Calf Pads	\$239

**Stock Options**

<input type="checkbox"/> 1497	Oxygen Holder Pkg Assembly	\$198
<input type="checkbox"/> 1496	IV/Oxygen Holder Pkg Assembly	\$198
<input type="checkbox"/> 1360	Adj. Rear Anti-Tippers (8" Casters)	\$109
<input type="checkbox"/> 1360A	Adj. Rear Anti-Tippers (6" Casters)	\$123
<input type="checkbox"/> A1576	Wheel Lock Extensions	\$63

*Elastic cord wheel lock extension.*

**FRAME & SEAT**

**FRAME TYPE**

<input type="checkbox"/> 9XT	9000XT	\$1,809
------------------------------	--------	---------

**FRAME WIDTH & DEPTH**

<input type="checkbox"/> WD46	14" X 16"	No Charge
<input type="checkbox"/> WD48	14" X 18"	No Charge
<input type="checkbox"/> WD56	15" X 16"	No Charge
<input type="checkbox"/> WD58	15" X 18"	No Charge
<input type="checkbox"/> WD66	16" X 16"	No Charge
<input type="checkbox"/> WD67	16" X 17"	No Charge
<input type="checkbox"/> WD68	16" X 18"	No Charge
<input type="checkbox"/> WD76	17" X 16"	No Charge
<input type="checkbox"/> WD77	17" X 17"	No Charge
<input type="checkbox"/> WD78	17" X 18"	No Charge
<input type="checkbox"/> WD86	18" X 16"	No Charge
<input type="checkbox"/> WD87	18" X 17"	No Charge
<input type="checkbox"/> WD88	18" X 18"	No Charge
1 <input type="checkbox"/> WD96	19" X 16"	No Charge
1 <input type="checkbox"/> WD97	19" X 17"	No Charge
1 <input type="checkbox"/> WD98	19" X 18"	No Charge
1 <input type="checkbox"/> WD06	20" X 16"	No Charge
1 <input type="checkbox"/> WD07	20" X 17"	No Charge
1 <input type="checkbox"/> WD08	20" X 18"	No Charge
1,2 <input type="checkbox"/> WD26	22" X 16"	\$206
1,2 <input type="checkbox"/> WD27	22" X 17"	\$330
1,2 <input type="checkbox"/> WD28	22" X 18"	\$330

**SEAT TO FLOOR HEIGHT**

3 <input type="checkbox"/> Adult	Adult (19-½" STF)	No Charge
3 <input type="checkbox"/> H	Hemi (17-½" STF)	No Charge
4.5 <input type="checkbox"/> SHE	Super Hemi (15-½" STF)	\$63
6 <input type="checkbox"/> SH	Super Low (14 ½" STF)	No Charge

**WEIGHT LIMITATION 250 LBS.**

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional knowledgeable of the particular limitations of the product.

**FOOTNOTES**

- 1 Not available with Super Low seat-to-floor.
- 2 Not available with adjustable angle back type.
- 3 Based on 24" wheels and 8" casters. See seat-to-floor matrix.
- 4 Based on 20" wheels and 6" casters. See seat-to-floor matrix.
- 5 Requires of AT902 (6" casters), U1763BK (20" composite wheels, urethane tires) AHDRM (aluminum handrims) and A1576 (wheel lock extensions). Seat-to-floor height is approx. 15"
- 6 Not available on ADJC and ADJCPB, upcharge will apply if required.

**Invacare® 9000 XT Wheelchair Accessories & Options**

**NOTES**

## ARM & BACK

### 7 ARM TYPE

<input type="checkbox"/>	18	Full Length Fixed Height Conventional(10")	\$251	
<input type="checkbox"/>	19	Full Length Adj.Height Conventional(10-14")	\$358	
<input type="checkbox"/>	28	Desk Length Fixed Height Conventional(10")	\$251	
<input type="checkbox"/>	29	Desk Length Adj.Height Conventional(10-14")	\$358	
2	<input type="checkbox"/>	78	Full Length Fixed Height Space-Saver(10")	\$251
2	<input type="checkbox"/>	79	Full Length Adj.Height Space-Saver(10-14")	\$358
2	<input type="checkbox"/>	88	Desk Length Fixed Height Space-Saver(10")	\$251
2	<input type="checkbox"/>	89	Desk Length Adj.Height Space-Saver(10-14")	\$358
2,8	<input type="checkbox"/>	79JY	Full Length Adjustable Ht Space-Saver(8-10")	\$358
2,8	<input type="checkbox"/>	89JY	Desk Length Adjustable Ht Space-Saver(8-10")	\$358

### BACK TYPE

<input type="checkbox"/>	ADJ	Adjustable Back w/10° Bend (15"-19")	<b>Standard</b>	
<input type="checkbox"/>	S1418	Straight Back w/Push Handles(14"-18")	No Charge	
1,9,10	<input type="checkbox"/>	ANGLE	Adj Angle Back w/10° Bend	\$165
1,9,10	<input type="checkbox"/>	ANGLE-S	Adj Angle Back-Straight	\$165

## COLOURS

### FRAME FINISH

<input type="checkbox"/>	08P	Deep Purple	No Charge
<input type="checkbox"/>	17P	Cobalt Blue	No Charge
<input type="checkbox"/>	24P	Wet Black	No Charge
<input type="checkbox"/>	34P	Satin (Flat) Black	<b>Standard</b>
<input type="checkbox"/>	71P	Silver Vein	No Charge
<input type="checkbox"/>	88P	Wet Black w/Twilight Sparkle	No Charge
<input type="checkbox"/>	105P	Midnight Blue	No Charge

### CANDY COLOURS

<input type="checkbox"/>	61P	Electric Red	No Charge
<input type="checkbox"/>	62P	Electric Blue	No Charge
<input type="checkbox"/>	66P	Electric Purple	No Charge
<input type="checkbox"/>	104P	Electric Teal	No Charge

### UPHOLSTERY COLOURS

<input type="checkbox"/>	U240	Black Nylon Upholstery	No Charge	
12	<input type="checkbox"/>	U245	Black Sail Cloth Adjust. Tension Seat	<b>Standard</b>
13	<input type="checkbox"/>	U67	BlackVinyl	No Charge

### NAUGAHYDE

13	<input type="checkbox"/>	WRAP	Wrap Around Vinyl Upholstery	No Charge
13,14	<input type="checkbox"/>	SOUPH	Screwed On Back Upholstery	\$140

## CASTERS

### CASTER OPTIONS

<input type="checkbox"/>	AT902	6" Composite/Urethane 3 Spoke	No Charge	
1	<input type="checkbox"/>	1228	8" X 1 3/4" Semi-Pneumatic Casters Lt. Grey	\$150
1	<input type="checkbox"/>	AT903	8" x 1 1/4" Composite/Urethane	<b>Standard</b>
1	<input type="checkbox"/>	1299	8" X 1" Solid Rubber Casters	No Charge
1	<input type="checkbox"/>	1292	8" X 1-1/4" Pneumatic Casters	No Charge
1	<input type="checkbox"/>	1492	8" X 1-1/4" Pneu w/Flat Free Insert	\$128

### 1,15 CASTER JOURNAL OPTION

16	<input type="checkbox"/>	ADJC	Adjustable Caster Journal	No Charge
16	<input type="checkbox"/>	ADJCPB	Adj. Precision Bearing Caster Journal	No Charge

### FOOTNOTES

7 Removable arm style, flipback.

8 Only Available with 20" Rear Wheels

9 Not available with adjustable tension back upholstery.

10 Standard with rear anti-tippers.

\$150 11 Only available with Super Low seat-to-floor.

**Standard** 12 U245 sail cloth adjustable tension seat upholstery available w/U240 black nylon back upholstery only.

No Charge 13 Vinyl upholstery must choose WRAP or SOUPH.

\$128 14 Not available with Adjustable Angle or Adjustable Angle-S back type.

### See Additional Charts on Page 6

\$63 Seat to Floor matrix w/6" Casters & Offset Axle Bracket

No Charge Seat to Floor matrix w/8" Casters & Offset Axle Bracket

# WHEELS

## WHEEL OPTIONS

<input type="checkbox"/>	1763BK	20" Composite/Pneumatic Tires	\$231
<input type="checkbox"/>	4763BK	20" Composite/Pneumatic-Flat Free Inserts	\$349
<input type="checkbox"/>	U1763BK	20" Composite/Urethane Tires	No Charge
<input type="checkbox"/>	1764BK	22" Composite/Pneumatic Tires	\$231
<input type="checkbox"/>	4764BK	22" Composite/Pneumatic-Flat Free Inserts	\$349
<input type="checkbox"/>	U1764BK	22" Composite/Urethane Tires	No Charge
<input type="checkbox"/>	2244BK	24" Composite/Pneumatic Tires	<b>Standard</b>
<input type="checkbox"/>	4244BK	24" Composite/Pneumatic-Flat Free Inserts	\$128
<input type="checkbox"/>	5244BK	24" Heavy Duty Comp/Pneu.-Flat Free Inserts	\$420
<input type="checkbox"/>	U2244BK	24" Composite/Urethane Treaded Tires	No Charge
<input type="checkbox"/>	U2222C	24" Composite/Urethane Corderd Tires	No Charge

## 19 WHEEL OPTIONS

### ONE ARM DRIVE (OAD)

<input type="checkbox"/>	9456L	One Arm Drive w/U2222C LH	\$672
<input type="checkbox"/>	9456R	One Arm Drive w/U2222C RH	\$672

### CLD ONE ARM DRIVE

<input type="checkbox"/>	LDR96RI	CLD One Arm Drive Right, 6" Caster	\$1,349
<input type="checkbox"/>	LDR96LI	CLD One Arm Drive Left, 6" Caster	\$1,349
<input type="checkbox"/>	LDR98RI	CLD One Arm Drive Right, 8" Caster	\$1,349
<input type="checkbox"/>	LDR98LI	CLD One Arm Drive Left, 8" Caster	\$1,349

## HANDRIM OPTIONS

<input type="checkbox"/>	AHDM	Aluminum Handrims	<b>Standard</b>
<input type="checkbox"/>	COM	Composite Handrims	No Charge
<input type="checkbox"/>	1512	Plastic Coated Aluminum Handrims	\$114
<input type="checkbox"/>	BHDM	Black Handrims (Alum for U2222)	No Charge
<input type="checkbox"/>	1224-8	8 Alum Oblique Proj,hdrm,plstic Coat	\$307
<input type="checkbox"/>	1225-8	8 Alum Vert Proj,hdrm,plstic Coat	\$307

## REAR AXLE

<input type="checkbox"/>	OFFAXL	Offset Axle	No Charge
--------------------------	--------	-------------	-----------

## 16 WHEEL LOCKS

<input type="checkbox"/>	9695	Hill Holder	\$149
<input type="checkbox"/>	A1567	Pull To Lock Wheel Lock-Clamp On	No Charge
<input type="checkbox"/>	1721	High Push To Lock Wheel Lock-Clamp On	No Charge
<input type="checkbox"/>	9690B	Pull to Lock Wheel Lock (Bolt Thru)(STD on Hemi)	No Charge
<input type="checkbox"/>	9691B	Push to Lock Wheel Lock (Bolt Thru)(STD on Adult)	No Charge
<input type="checkbox"/>	9690	Pull To Lock Wheel Locks (Disk Style)	No Charge
<input type="checkbox"/>	9691	Push To Lock Wheel Locks (Disk Style)	No Charge

# POSITIONING STRAPS

## SEAT POSITIONING STRAPS

<input type="checkbox"/>	1515	Auto Style Buckle Seat Positioning Strap	\$84
<input type="checkbox"/>	1321BK	Hook and Loop Seat Positioning Strap	\$66

## CHEST POSITIONING STRAPS

<input type="checkbox"/>	7311BK	Airline Buckle Chest Positioning Strap	\$84
<input type="checkbox"/>	9515	Auto Style Chest Positioning Strap	\$84
<input type="checkbox"/>	7321BK	Hook and Loop Chest Positioning Strap	\$66

## FOOTNOTES

- 15 Fixed caster journals are standard.
- 16 Not available with Adjustable Angle or Adjustable Angle-S back type.
- 17 High profile treaded urethane tires.
- 18 Low profile non-treaded Urethane
- 19 The Cyclical Lever Drive requires treaded rear wheels and is only available on 16, 17, and 18" deep. Hemi STF must use Space Saver arms. Weight limit is 250 lbs.
- 20 Not available w/amputee attachment or offset axle
- 21 Not available with adjustable caster journal..
- 22 Not available w/U2222 rear wheels.
- 23 Not available on 20" and 22" composite rear wheels.
- 24 Not available with urethane tires or 20" rear wheels.
- 25 Only available w/Vinyl upholstery. 3/4" block letters.

# ACCESSORIES & STOCK OPTIONS

## BACK ACCESSORIES

<input type="checkbox"/>	ADJTENS	Adjustable Tension Back Upholstery	\$142
<input type="checkbox"/>	1160	Carrying Pocket	\$83
<input type="checkbox"/>	STC	Stencil Back Upholstery	\$59
<input type="checkbox"/>	RB	Dual Back Liner	\$91
<input type="checkbox"/>	L8601	Omit Back Upholstery	Less net (20.00)
<input type="checkbox"/>	BP02	Back Pack	\$105

# Invacare® 9000 XT Wheelchair Accessories & Options

# NOTES

## SEAT ACCESSORIES

<input type="checkbox"/>	9120	Solid Seat Insert	\$138
<input type="checkbox"/>	RS	Dual Seat Liner	\$87
<input type="checkbox"/>	L8651	Omit Seat Upholstery	Less net (15.00)

NOTE: All specifications and dimensions are approximate.  
For the most current pricing information visit [www.invacare.ca](http://www.invacare.ca)

## ARM ACCESSORIES

<input type="checkbox"/> NRAS	Non-Removable Arm With Screws	No Charge
27 <input type="checkbox"/> 1300	Utility Tray	\$268

## MISCELLANEOUS ACCESSORIES

28,29 <input type="checkbox"/> 9357L	Amputee Attachment	\$131
30 <input type="checkbox"/> 1305L	Crutch & Cane Carrier	\$66
28,29 <input type="checkbox"/> 1577QR	Quick Release Axles	\$142
30 <input type="checkbox"/> 2218	IV Rod, Telescopic	\$185
<input type="checkbox"/> 2226	Front Anti-Tippers	\$158
11 <input type="checkbox"/> 9357SH	Super Low Amputee Attachment	\$131
11 <input type="checkbox"/> 1360SH	Super Low Adjustable Rear Anti-Tipper	\$109

## STOCK OPTIONS

28,30,31 <input type="checkbox"/> 1497	Oxygen Holder Package Assembly	\$198
28,30,32 <input type="checkbox"/> 1496	IV/Oxygen Holder Package Assembly	\$198
29,30,33,34 <input type="checkbox"/> 1360	Adjustable Rear Anti-Tipper (8" Casters)	\$109
29,30,35 <input type="checkbox"/> 1360A	Adjustable Rear Anti-Tipper (6" Casters)	\$123
36 <input type="checkbox"/> A1576	Wheel Lock Extensions	\$63

## FRONT RIGGING OPTIONS

### FRONT RIGGING

11,37 <input type="checkbox"/> SH904M	Super Low Elevating Legrest 11"-17"	\$356
11,37 <input type="checkbox"/> SH904MS	Super Low Elevating Legrests 14"wide	\$356
11 <input type="checkbox"/> SH93E	Super Low Ext Sw Footrests 11"-17"	\$230
11 <input type="checkbox"/> SH93ES	Super Low Ext Sw Footrests 14"wide	\$230
37 <input type="checkbox"/> AHL4A	Hemi Smartleg Artic. Legrests 16" wide Alum. Ftpts.	\$399
37 <input type="checkbox"/> AHL4C	Hemi Smartleg Artic. Legrests 16" wide Comp. Ftpts.	\$399
37 <input type="checkbox"/> AHS4A	Hemi Smartleg Articulating Legrests 14" wide	\$399
<input type="checkbox"/> T93HA	Hemi Sw Footrests Alum. Ftpts.	\$159
<input type="checkbox"/> T93HC	Hemi Sw Footrests Comp. Ftpts.	\$159
37 <input type="checkbox"/> T94HS	Hemi Sw Elevating Legrests Alum Ftpt for 14" wide	\$302
37 <input type="checkbox"/> T94HA	Hemi Sw Elevating Legrests Alum. Ftpts.	\$309
37 <input type="checkbox"/> T94HC	Hemi Sw Elevating Legrests Comp. Ftpts.	\$302
<input type="checkbox"/> T93HS	Hemi Sw Footrests Alum. Ftpts. 14"wide	\$158
<input type="checkbox"/> T963A	60 Deg Footrest Aluminum Footplates	\$158
<input type="checkbox"/> T963C	60 Deg Footrest Composite Footplates	\$158
<input type="checkbox"/> T963S	60 Deg Footrest Composite Footplate 14"wide	\$158
<input type="checkbox"/> T93HAB	Hemi Bolt-Thru Footrests, Aluminum Ftpts	\$158
<input type="checkbox"/> T93HCB	Hemi Bolt-Thru Footrests, Composite Ftpts	\$158
37 <input type="checkbox"/> T94HAB	Hemi Bolt-Thru Legrests, Aluminum Ftpts	\$302
37 <input type="checkbox"/> T94HCB	Hemi Bolt-Thru Legrests, Composite Ftpts	\$302

### FOOTNOTES

- 26 Not available on 22" wide chairs.  
 27 Not available with 24" rear wheel set in Hemi position.  
 28 Not available with offset axle.  
 29 Not available with any one arm drive.  
 30 Not available with traditional one arm drive, but is available with CLD one arm drive.  
 31 For 0° and 6° seat angles. See service manual for details.  
 32 Not available with amputee attachment.  
 34 Standard w/Adjustable Angle Back.  
 35 Only available with 6" casters.  
 36 Elastic cord wheel lock extensions. Hemi seat to floors come standard with non-cord wheel lock extensions.  
 37 Elevating legrests available in chrome or black finish.  
 38 Only available on T944 elevating legrests.  
 39 Standard on footrests.  
 40 Only available w/H93M front riggings.

### FOOTPLATE OPTIONS

<input type="checkbox"/> 1350	Extra Large Aluminum Footplates	\$106
<input type="checkbox"/> AT5543	Adjustable Angle Flip Up Footplates	\$132

### FOOTREST / LEGREST ACCESSORIES

<input type="checkbox"/> 1820	Non-Removable Device	\$41
<input type="checkbox"/> 1337	Calf Strap	\$63
<input type="checkbox"/> IMPACTG	Impact Guards	No Charge
39 <input type="checkbox"/> 1346	Heel Loops	\$63

Seat-To-Floor w 6" Casters and Offset Axle Bracket					
Seat-To Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position	Axle Plate Mounting On Frame	
				Top Lug	Bottom Lug
15.5"	Top Hole	22"	#5 Hole	#4 Hole	#1 Hole
15.5"	Top Hole	20"	#3 Hole	#4 Hole	#1 Hole
16.5"	Middle Hole	24"	#5 Hole	#4 Hole	#1 Hole
16.5"	Middle Hole	22"	#3 Hole	#4 Hole	#1 Hole
16.5"	Middle Hole	20"	#2 Hole	#4 Hole	#1 Hole
17.5"	Bottom Hole	24"	#3 Hole	#4 Hole	#1 Hole
17.5"	Bottom Hole	22"	#2 Hole	#4 Hole	#1 Hole
17.5"	Bottom Hole	20"	#3 Hole	#4 Hole	#1 Hole

For a complete list of seat-to-floor heights and angles that are configurable by the dealer refer to the service manual


Seat-To-Floor Matrix w 8" Casters and offset Axle Bracket-					
Seat-To Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position	Axle Plate Mounting On Frame	
				Top Lug	Bottom Lug
17.5"	Top Hole	20"	# 3 Hole	#1 Hole	#4 Hole
17.5"	Top Hole	22"	# 2 Hole	#4 Hole	#1 Hole
17.5"	Top Hole	24"	# 3 Hole	#4 Hole	#1 Hole
18.5"	Middle Hole	22"	# 3 Hole	#1 Hole	#4 Hole
18.5"	Middle Hole	24"	# 2 Hole	#4 Hole	#1 Hole
19.5"	Bottom Hole	20"	# 3 Hole	#1 Hole	#4 Hole

For a complete of seat-to-floor heights and angles that are configurable by the dealer refer to the service manual.

SEAT-TO-FLOOR MATRIX 6" CASTERS			
Seat to Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
15 1/2"	Top	20"	Top
16 1/2"	Middle	22"	Top
17 1/2"	Bottom	24"	Top
17 1/2"	Bottom	20"	Bottom

SEAT-TO-FLOOR MATRIX 8" CASTERS			
Seat to Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
17 1/2"	Top	20"	Bottom
17 1/2"	Top	24"	Top
18 1/2"	Middle	22"	Bottom
19 1/2"	Bottom	24"	Bottom

SEAT-TO-FLOOR MATRIX 6" CASTERS & SUPER LOW FRAME			
Seat to Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
13 1/2"	Top	20"	Only One
14 1/2"	Bottom	22"	Only One



Invacare Canada reserves the right to change prices and specifications without notice.  
 HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324 (Canada) Fax: (905) 890-5244 (800) 668-5478 (Canada)  
 Invacare®, and Yes, you can., are trademarks of Invacare Corporation.  
 WWW.INVACARE.CA