



NAME _____ PHONE NO _____
 ACCOUNT # _____ P.O. _____ DATE _____
 ADDRESS _____

FRAME & SEAT

CHAIR TYPE

9RC 9000XT Recliner \$2,618

FRAME WIDTH & DEPTH

<input type="checkbox"/> WD46	14" Wide X 16" Deep	No Charge
<input type="checkbox"/> WD47	14" Wide X 17" Deep	\$145
<input type="checkbox"/> WD48	14" Wide X 18" Deep	\$145
<input type="checkbox"/> WD66	16" Wide X 16" Deep	No Charge
<input type="checkbox"/> WD67	16" Wide X 17" Deep	\$145
<input type="checkbox"/> WD68	16" Wide X 18" Deep	\$145
<input type="checkbox"/> WD86	18" Wide X 16" Deep	No Charge
<input type="checkbox"/> WD87	18" Wide X 17" Deep	\$145
<input type="checkbox"/> WD88	18" Wide X 18" Deep	\$145
<input type="checkbox"/> WD96	19" Wide X 16" Deep	No Charge
<input type="checkbox"/> WD97	19" Wide X 17" Deep	\$145
<input type="checkbox"/> WD98	19" Wide X 18" Deep	\$145
<input type="checkbox"/> WD06	20" Wide X 16" Deep	No Charge
<input type="checkbox"/> WD07	20" Wide X 17" Deep	\$145
<input type="checkbox"/> WD08	20" Wide X 18" Deep	\$145
<input type="checkbox"/> WD26	22" Wide X 16" Deep	\$206
<input type="checkbox"/> WD27	22" Wide X 17" Deep	\$333
<input type="checkbox"/> WD28	22" Wide X 18" Deep	\$333

FRAME TYPE

R Recliner \$774

This is an automatic pick in the configurator.

ARM & BACK

ARM TYPE

Fixed Height Arms 9", Adjustable Height Arms 9"-13", removable style.

<input type="checkbox"/> 18	Full Length Fixed Height Conventional	\$261
<input type="checkbox"/> 28	Desk Length Fixed Height Conventional	\$261
<input type="checkbox"/> 19	Full Length Adj. Height Conventional	\$372
<input type="checkbox"/> 29	Desk Length Adj. Height Conventional	\$372
<input type="checkbox"/> 19JY	Full length Adj. Height Conventional (8"-10")	\$372
<input type="checkbox"/> 29JY	Desk length Adj. Height Conventional (8"-10")	\$372

BACK TYPE

Add additional 11" for the headrest extension.

<input type="checkbox"/> 2R	Reclining Back 23" high	No Charge
<input type="checkbox"/> L	Reclining Back 23" high w/Lift Package	\$130

COLOURS

FRAME FINISH

<input type="checkbox"/> 17P	Cobalt Blue	No Charge
<input type="checkbox"/> 61P	Electric Red	No Charge
<input type="checkbox"/> 24P	Wet Black	No Charge
<input type="checkbox"/> 34P	Satin (Flat) Black	Standard
<input type="checkbox"/> 71P	Silver Vein	No Charge

CHAIR FEATURES AND OPTIONS

- ◆ Chrome plated steel frame
- ◆ Choice of frame colours available
- ◆ Conventional adjustable and fixed-height arm styles
- ◆ 24" composite/Pneumatic tire
- ◆ Degree of recline from 90° to 168°
- ◆ 22" high back and 10" headrest extension
- ◆ {removable} with pillow standard
- ◆ Non-swing back removable arm
- ◆ Attendant operated wheel locks
- ◆ Lifetime warranty on frame and crossbraces
- ◆ Standard with rear anti-tippers
- ◆ STF = 19.5" with 24" wheels
- ◆ STF = 18.5" with 22" wheels
- ◆ STF = 17.5" with 20" wheels

WEIGHT LIMITATION 250 LBS.

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

UPHOLSTERY COLOURS

All back upholstery is screwed on.

- U240 Black Nylon Upholstery **Standard**
- U245 Black Sail Cloth Adjust. Tension Seat No Charge
- U245 sail cloth adjustable tension seat upholstery available w/U240 black nylon back upholstery only.*
- U550 Midnight Blue Vinyl Upholstery No Charge
- U67 Black Vinyl No Charge
- U75 Dove Grey Vinyl No Charge
- U86 Ruby Vinyl No Charge
- U88 Camel Vinyl No Charge

CASTERS

CASTER OPTIONS

- 1228 8" X 1 3/4" Semi-Pneu Caster Lt. Grey \$150
- AT903 8" X 1-1/4" Composite/Urethane **Standard**
- 1290 8" X 1" Solid Rubber Casters Charcoal No Charge
- 1299 8" X 1" Solid Rubber Casters Dark Grey No Charge
- 1292 8" X 1-1/4" Pneumatic Casters No Charge
- 1492 8" X 1-1/4" Pneu w/Flat Free Insert \$128

Seat-to-Floor Matrix 8" Casters			
Seat to Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
17.50"	Top	20"	1 Position Only
18.50"	Middle	22"	1 Position Only
19.50"	Bottom	24"	1 Position Only

WHEELS

WHEEL OPTIONS

- 1763BK 20" Composite/Pneumatic Tires \$231
- 4763BK 20" Composite/Pneumatic-Flat Free Inserts \$349
- U1763BK 20" Composite/Urethane Tires No Charge
- 1764BK 22" Composite/Pneumatic Tires \$231
- 4764BK 22" Composite/Pneumatic-Flat Free Inserts \$349
- U1764BK 22" Composite/Urethane Tires No Charge
- 2244BK 24" Composite/Pneumatic Tires **Standard**
- 4244BK 24" Composite/Pneumatic-Flat Free Inserts \$128
- U2244BK 24" Composite/Urethane Treaded Tires No Charge
- U2222C 24" Composite/Urethane Corded Tires No Charge

Not available with AHDRM or 1225-12.

HANDRIM OPTIONS

- AHDRM Aluminum Handrims **Standard**
- Not available with U2222C rear wheels.*
- COM Composite Handrims No Charge
- Available with 24" rear wheels only.*
- 1512 Plastic Coated Aluminum Handrims \$114
- 1224-8 8 Alum Oblique Proj. Hdrm. Plastic Coated \$307
- Available with 24" rear wheels only.*
- 1225-8 8 Alum Vert. Proj. Hdrm., Plastic Coated \$307
- Available with 24" rear wheels only.*

WHEEL LOCKS

- Attendant Operated **Standard**
- 1329BK Toggle Wheel Lock Extension (8") \$63
- 9696R Push to Lock Wheel Locks (Patient Operated) \$48
- Attendant operated wheel locks are standard and are not available with patient operated wheel locks.*
- 1245 Pull to Lock Wheel Locks (Patient Operated) No Charge
- Attendant operated wheel locks are standard and are not available with patient operated wheel locks.*

POSITIONING STRAPS

SEAT POSITIONING STRAPS

- 1515 Auto Style Buckle Seat Positioning Strap \$84
- 1321BK Hook & Loop Seat Positioning Strap \$66

CHEST POSITIONING STRAPS

- 9515 Auto Style Chest Positioning Strap \$84

ACCESSORIES & OPTIONS

BACK ACCESSORIES

- 1160 Carrying Pocket \$83
- 3040 Torso Support for Recliner \$175

<input type="checkbox"/>	STC	Stencil Back Upholstery	\$59
		<i>3/4" block letters. Only available with vinyl upholstery.</i>	
<input type="checkbox"/>	3046	Head Support for Recliner	\$163
<input type="checkbox"/>	RB	Dual Back Liner	\$91
<input type="checkbox"/>	L8601	Omit Back Upholstery	Less Net (\$20)

SEAT ACCESSORIES

<input type="checkbox"/>	9120	Solid Seat Insert	\$144
<input type="checkbox"/>	RS	Dual Seat Liner	\$90

ARM ACCESSORIES

<input type="checkbox"/>	NRAS	Non-Removable Arm With Screws	No Charge
<input type="checkbox"/>	1300	Utility Tray	\$268

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/>	2219	IV Rod, Telescopic (Recliner)	\$163
<input type="checkbox"/>	2226	Front Anti-Tippers	\$158
<input type="checkbox"/>	1100	Tool Kit	No Charge

FRONT RIGGING OPTIONS

FRONT RIGGING

<input type="checkbox"/>	AHL4A	Hemi Smartleg Artic. Legrests 16" wide Alum. Footplates.	\$399
		<i>Legrest available in black finish only.</i>	
<input type="checkbox"/>	AHL4C	Hemi Smartleg Artic. Legrests 16" wide Comp. Footplates.	\$399
		<i>Legrest available in black finish only.</i>	
<input type="checkbox"/>	AHS4A	Hemi Smartleg Articulating Legrests 14"-15" Wide	\$392
		<i>Legrest available in black finish only. Only available on 14" wide.</i>	
<input type="checkbox"/>	T93HA	Hemi Sw Footrests Aluminum Footplates.	\$159
<input type="checkbox"/>	T93HC	Hemi Sw Footrests Composite Footplates.	\$159
<input type="checkbox"/>	T94HA	Hemi Sw Elevating Legrests Aluminum Footplates.	\$309
		<i>Legrest available in black finish only.</i>	
<input type="checkbox"/>	T94HC	Hemi Sw Elevating Legrests Composite Footplates.	\$302
		<i>Legrest available in black finish only.</i>	
<input type="checkbox"/>	T94HS	Hemi Sw Elevating Legrests Aluminum Footplate 14" wide	\$302
		<i>Legrest available in black finish only. Only available on 14" wide.</i>	
<input type="checkbox"/>	T93HS	Hemi Sw Footrests Aluminum Footplates. 14" wide	\$158
		<i>Only available on 14" wide.</i>	
<input type="checkbox"/>	T93HAB	Hemi Bolt-Thru First, Aluminum Footplates	\$158
<input type="checkbox"/>	T93HCB	Hemi Bolt-Thru First, Composite Footplates	\$158
<input type="checkbox"/>	T94HAB	Hemi Bolt-Thru Legrest, Aluminum Footplates	\$302
<input type="checkbox"/>	T94HCB	Hemi Bolt-Thru Legrest, Composite Footplates	\$302
<input type="checkbox"/>	T93HSB	Hemi Bolt-Thru First, Aluminum Footplates 14" wide	\$158
		<i>Only available on 14" wide.</i>	
<input type="checkbox"/>	T94HSB	Hemi Bolt-Thru Legrest, Aluminum Footplates 14" wide	\$158
		<i>Only available on 14" wide.</i>	

FOOTPLATE OPTIONS

<input type="checkbox"/>	1350	Extra Large Aluminum Footplates	\$106
<input type="checkbox"/>	AT5543	Adjustable Angle Flip-up Footplates	\$132

FOOTREST / LEGREST ACCESSORIES

<input type="checkbox"/>	1820	Non-Removable Device	\$41
		<i>Not Available with AHL4A or AHL4C elevating legrests.</i>	
<input type="checkbox"/>	1337	Calf Strap	\$63
<input type="checkbox"/>	1600BK	Heel Loops With Ankle Straps	\$84
<input type="checkbox"/>	1348	Legrest Panel	\$93
		<i>Not available with Hemi minimum Swingaway footrests.</i>	
<input type="checkbox"/>	1346	Heel Loops	\$63
		<i>Standard on footrests.</i>	
<input type="checkbox"/>	L8689	Omit Heel Loops	Less Net (\$10)
<input type="checkbox"/>	1319	Black Non-Padded Calf Pads	Less Net (\$10)
		<i>Not available with AHL4A or AHL4C elevating legrests.</i>	

Invacare Canada reserves the right to change prices and specifications without notice.

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324 (Canada) Fax: (905) 890-5244 (800) 668-5478 (Canada)
 Invacare®, and Yes, you can., are trademarks of Invacare Corporation.
 WWW.INVACARE.CA