



NAME _____ PHONE NO _____

ACCOUNT # _____ P.O. _____ DATE _____

ADDRESS _____

FRAME & SEAT

FRAME TYPE

☒ 9JYLT 9000 Jymni Light Weight Chair \$1,539

FRAME WIDTH & DEPTH

☐ WD00 10" X 10" No Charge
Not available with 17" seat-to-floor height, unless chair is a recliner.
☐ WD20 12" X 10" No Charge
Not available with 17" seat-to-floor height, unless chair is a recliner.
☐ WD22 12" X 12" No Charge
Not available with 17" seat-to-floor height, unless chair is a recliner.
☐ WD42 14" X 12" No Charge
Not available with 17" seat-to-floor height, unless chair is a recliner.
☐ WD44 14" X 14" No Charge
☐ WD64 16" X 14" No Charge

SEAT TO FLOOR HEIGHT

☐ 15STF 15" Seat To Floor Height No Charge
Only available with 20" rear wheels
☐ 16STF 16" Seat To Floor Height **Standard**
Only available with 20" rear wheels
☐ 17STF 17" Seat To Floor Height No Charge
Only available with 22" rear wheels.

COLOUR

FRAME FINISH

☐ 08P Deep Purple No Charge
☐ 24P Wet Black No Charge
☐ 34P Satin (Flat) Black No Charge
☐ 71P Silver Vein No Charge
☐ 61P Electric Red No Charge
☐ 62P Electric Blue No Charge
☐ 66P Electric Purple No Charge
☐ 104P Electric Teal No Charge
☐ 105P Midnight Blue No Charge

UPHOLSTERY COLOURS

☐ U240 Black Nylon Upholstery **Standard**
Recliner is only available in U240 black nylon upholstery. Recliner back height 16" plus 11" headrest extension.
☐ U550 Midnight Blue Vinyl Upholstery No Charge
☐ U67 Black Vinyl Upholstery No Charge

ARM & BACK

ARM TYPE

All above selections come standard with fixed height arms 7 1/2", adjustable height arms 8" - 10" in 1" increments. Removable flip-back arms.

☐ 18JY Full Length Fixed Height Conventional \$261
ANGLE and Recliner back must use conventional arms. Arms do not flip back on the recliner.
☐ 19JY Full Length Adjustable Height Conventional \$372
ANGLE and Recliner back must use conventional arms. Arms do not flip back on the recliner.
☐ 28JY Desk Length Fixed Height Conventional \$261
ANGLE and Recliner back must use conventional arms. Arms do not flip back on the recliner.
☐ 29JY Desk Length Adj Height Conventional \$372
ANGLE and Recliner back must use conventional arms. Arms do not flip back on the recliner.
☐ 78JY Full Length Fixed Height Space-Saver \$261
☐ 79JY Full Length Adjustable Height Space-Saver \$372
☐ 88JY Desk Length Fixed Height Space-Saver \$261
☐ 89JY Desk Length Adjustable Ht Space-Saver \$372

BACK TYPE

All below selections come standard with fixed height arms 7 1/2", adjustable height arms. 8"H - 10"H in 1" increments. Removable flip-back arms

☐ ADJBH Adjustable Ht Back with 10° Bend **Standard**
Adjustable back height 14"-15". Factory Set at 14".
☐ ANGLE Adjustable Angle Back with 10° Bend \$172
RECLINER and angle back must use conventional arms. Arms do not flip back on RECLINER.
This selection comes standard with rear anti-tippers.
☐ R Recliner Frame \$774
Recliner is only available in U240 black nylon upholstery. Recliner back height 16" plus 11" headrest extension.
Recliner back must use conventional arms. Comes standard with rear anti-tippers. Arms do not flip back on the recliner.
This selection must be made when WD42 & 17STF are selected.

POSITIONING STRAPS**SEAT POSITIONING STRAPS**

<input type="checkbox"/> 1311BK	Airline Buckle Seat Positioning Strap	\$84
<input type="checkbox"/> 1515	Push Button Style Seat Positioning Strap	\$84
<input type="checkbox"/> 1321BK	Hook & Loop Seat Positioning Strap	\$63

CHEST POSITIONING STRAPS

<input type="checkbox"/> 9515	Push Button Style Chest Positioning Strap	\$84
-------------------------------	---	------

CASTER AND WHEEL OPTIONS**CASTER OPTIONS**

<input type="checkbox"/> 1687	6" X 1" Dark Grey Comp. Urethane Casters
<input type="checkbox"/> 1686	6" X 1" Composite Urethane Casters
<input type="checkbox"/> AT902	6" Composite Urethane. 3 Spoke

Standard

No Charge

No Charge

WHEEL OPTIONS

<input type="checkbox"/> 1763BK	20" Composite Pneumatic Tires
<input type="checkbox"/> 4763BK	20" Composite Pneumatic with Flat Free Inserts
<input type="checkbox"/> U1763BK	20" Composite Urethane Tires
<input type="checkbox"/> 1764BK	22" Composite Pneumatic Tires
<input type="checkbox"/> 4764BK	22" Composite Pneumatic with Flat Free Inserts
<input type="checkbox"/> U1764BK	22" Composite Urethane Tires

\$231

\$349

Standard

\$231

\$349

No Charge

HANDRIM OPTIONS

<input type="checkbox"/> BHDRM	Black Handrims	Standard
--------------------------------	----------------	----------

ACCESSORIES**BACK ACCESSORIES**

<input type="checkbox"/> 1483JY	Stroller Handle	\$156
<i>Not available on Recliner or ANGLE back type.</i>		
<input type="checkbox"/> 9130	Solid Back Insert	\$142
<i>Only available in 14" back height.</i>		
<input type="checkbox"/> L8601	Omit Back Upholstery	Less net (20.00)
<input type="checkbox"/> 1483R	Headrest with Push Handles Extension	\$165

SEAT ACCESSORIES

<input type="checkbox"/> 9120	Solid Seat Insert	\$144
<input type="checkbox"/> L8651	Omit Seat Upholstery	Less net (20.00)

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/> 9357JY	1" Forward Positioner	\$104
<i>Comes standard with rear anti-tippers. Not available with Recliner, Angle back type, or 10" and 12" deep chairs.</i>		
<input type="checkbox"/> A1576	Wheel Lock Extensions	\$63
<input type="checkbox"/> 9691	Push To Lock Wheel Locks	No Charge
<i>Pull-to-lock wheel locks are standard.</i>		
<input type="checkbox"/> 1577QR	Quick Release Axles	\$137
<input type="checkbox"/> 1487	IV/Oxygen Holder Package Assembly	\$198
<i>Not available on Recliner or ANGLE back type.</i>		
<input type="checkbox"/> 1488	Oxygen Holder Package Assembly	\$198
<i>Not available on Recliner or ANGLE back type.</i>		
<input type="checkbox"/> 1360	Adjustable Rear Anti-Tippers	\$109
<i>Anti-tippers for 22" rear wheels.</i>		
<input type="checkbox"/> 1360B	Adj Rear Anti-Tippers 20" Rear Wheels	\$109
<i>Anti-tippers for 20" rear wheels.</i>		

FRONT RIGGING OPTIONS**FRONT RIGGING**

<input type="checkbox"/> J932A	70° Footrests with Aluminum Footplates	\$158
<input type="checkbox"/> J932A-14	80° Footrests with Alum. Footplates-12 & 14" deep	\$158
<i>Not available on 10" deep frames.</i>		
<input type="checkbox"/> J933A	90° Footrests with Aluminum Footplates	\$158
<input type="checkbox"/> J934A	Elevating Legrests with Aluminum Footplates	\$302
<i>Knee to Heel Range is Required with J934A.</i>		
<input type="checkbox"/> J934A-14	Elevating Legrests with Alum Footplates-14" deep	\$309
<i>Not available on 10" & 12" deep frames.</i>		

FOOTREST / LEGREST ACCESSORIES

<input type="checkbox"/> 1820	Non-Removable Device	\$41
<input type="checkbox"/> 1346	Heel Loops	\$63
<i>Standard on swingaway footrests.</i>		

KNEEL TO HEEL RANGE (only for J934A Elevating Legrest choice)*Knee-to-Heel range must be selected if J934A is chosen*

<input type="checkbox"/> SHEXT	7"- 10" Kneel to Heel (Short Extended)	No Charge
<input type="checkbox"/> MDEXT	9"-12" Kneel to Heel (Medium Extended)	No Charge
<input type="checkbox"/> LGEXT	11"-14" Kneel to Heel (Long Extended)	No Charge

Seat to Floor Matrix 6" Casters			
Seat to Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
15"	Top	20"	Top
16"	Middle	20"	Bottom
17"	Bottom	22"	Bottom

Weight Limitation 125 Pounds	
Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.	

Invacare Canada reserves the right to change prices and specifications without notice.

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324 (Canada) Fax: (905) 890-5244 (800) 668-5478 (Canada)

Invacare®, and Yes, you can., are trademarks of Invacare Corporation.

WWW.INVACARE.CA

NOTE: All specifications and dimensions are approximate.

For the most current pricing information visit www.invacare.ca