



**Invacare® 9000 TOPAZ
Wheelchair
AADL Customers Only
AADL # W689
Effective April 1, 2022**

NAME _____ PHONE NO _____

ACCOUNT # _____ P.O. _____ DATE _____

ADDRESS _____

FRAME & SEAT

FRAME TYPE

9TPZ 9000 Topaz \$3,299

CUSTOM PACKAGES

HDPG Heavy Duty Package \$456
Reinforced frame weight capacity of 1000 pounds. Only available in 28" and 30" Frame widths.

FRAME WIDTH

W20 20" Seat Width. No Charge
Not available with HDPG option.
 W22 22" Seat Width. No Charge
Not available with HDPG option.
 W24 24" Seat Width. No Charge
Not available with HDPG option.
 W26 26" Seat Width. No Charge
Not available with HDPG option.
 W28 28" Seat Width. No Charge
 W30 30" Seat Width. No Charge

FRAME DEPTH

D18 18" Seat Depth No Charge
 D20 20" Seat Depth \$333

SEAT TO FLOOR HEIGHT

ADULT Adult 19.5" No Charge
Based on 24" wheels and 8" casters.
 H Hemi 17.5" No Charge
Based on 24" wheels and 8" casters.
 9TSH Super Hemi 15.5" No Charge
Not available with HDPG option. Only available with 6" casters and 20" rear wheels

ARM & BACK SELECTIONS

ARM TYPE

Fixed height arms 9.5", adjustable height arms 9"-12" removable arm style.

T18 Full Length Fixed Height Conventional \$261
 T19 Full Length Adj Height Conventional \$372
 T28 Desk Length Fixed Height Conventional \$261
 T29 Desk Length Adj Height Conventional \$372

BACK TYPE

ADJH Adjustable Height Back w/10° Bend **Standard**
Adjustable back height 16" - 18". Factory set at 16".

CASTER & WHEEL SELECTIONS

CASTER OPTIONS

1259 6" X 2" Semi-Pneumatic Caster \$132
Not available with HDPG option. Only available with Super Hemi seat-to-floor.
 1228BK 8" X 1 3/4" Semi-Pneumatic Casters Lt. Grey **Standard**
Not available with HDPG option. Only available with Adult and Hemi STF and AHDRM handrims
 1228HD 8" X 2" Semi-Pneumatic Casters Black No Charge
Only available with HDPG option.

WHEEL OPTION

5763BK 20" HD Pneu Tires w/Flat Free Inserts \$420
Not available with HDPG option. Only available with Super Hemi seat-to-floor. *Check configurator for 700 lb weight capacity.*
 U5222SP 24" Spoked HD Urethane Treaded Tires 1000 Lb Capacity No Charge
Only available with Adult and Hemi STF and AHDRM handrim. Only available with HDPG option.
 U5700SP Spoked HD Urethane Treaded Tires 700 Lb Capacity No Charge

HANDRIM OPTIONS

AHDRM Aluminum Handrims No Charge

Weight Limit 700 Pounds
Invacare recommends all wheelchairs to be fitted to the consumer under the guidance of a health care professional knowledgeable of the particular limitations of the product.

FRAME FINISH & UPHOLSTERY

FRAME FINISH

- 24P Wet Black No Charge
- 71P Silver Vein **Standard**

UPHOLSTERY COLOURS

- U240 Black Nylon Upholstery **Standard**
- U550 Midnight Blue Vinyl Upholstery No Charge
- U67 Black Vinyl No Charge

FRONT RIGGING OPTIONS

FRONT RIGGING

- TZ33A Swingaway Footrests w/Extra Large Footplates \$434
- TZ34A Elevating Legrests w/Extra Large Footplates \$395

FOOTREST / LEGREST ACCESSORIES

- 1337 Calf Strap \$63
 - 1346 Heel Loops \$63
- Standard on swingaway footrests.*

SEAT ACCESSORIES

SEAT POSITIONING STRAPS

- T1515 Auto Style Buckle Seat Positioning Strap \$84

BACK ACCESSORIES

BACK ACCESSORIES

- T1483 Push Bar \$156
- This accessory shipped in a separate carton.*

SEAT ACCESSORIES

- T9120 Solid Seat Insert \$140
- This accessory shipped in a separate carton. Manufactured in vinyl upholstery only.*

ACCESSORIES

These accessories shipped in a separate carton. Manufactured in vinyl upholstery only.

ACCESSORIES

- T1496 IV/Oxygen Holder Package Assembly \$264
- T1497 Oxygen Holder Package Assembly \$198

MISCELLANEOUS ACCESSORIES

- T1358 Rear Anti-Tippers \$109
 - A1567 Pull to Lock Clamp on Wheel locks (Std on Hemi). No Charge
 - 1721 High Push to Lock Clamp on Wheel locks (Std on Adult) No Charge
- Only available with HDPG option.*
- 1305L Crutch & Cane Carrier \$66
 - A1576 Wheel Lock Extensions \$63
 - 1335XL Non-Folding Device - Topaz \$63

Seat-to-Floor 6" Casters			
Seat-to-Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
15 1/2"	1 Position Only	20"	Top

Seat-to-Floor 8" Casters			
Seat-to-Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
17 1/2"	Top	24"	Top
19 1/2"	Bottom	24"	Bottom

Invacare Canada reserves the right to change prices and specifications without notice.

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324 (Canada) Fax: (905) 890-5244 (800) 668-5478 (Canada)
 Invacare®, and Yes, you can., are trademarks of Invacare Corporation.
 WWW.INVACARE.CA