

# INVACARE® CANADA DOLOMITE® ROLLATORS



Yes, you can.®



**Yes, you can.®**

Your mobility is your life and Dolomite Rollators are designed to help you live life to the fullest.

Dolomite Rollators improve the mobility of people with many different conditions, reducing their risk of falling and increasing their independence and quality of life.

If you are living with any of the following conditions, Invacare recommends the appropriate rollator.

## **Fractured Hip or Hip Replacement**

---

Recommended Rollator:

- All models, depending on user's size and daily mobility requirements

## **Stroke Patients**

---

Recommended Rollator:

- Alpha with one-hand brake, Symphony, with one-hand brake and HP Handle

## **Respiratory Conditions**

---

Recommended Rollator:

- All models (except Soprano and Jazz) fitted with the Oxygen Tank Holder

## **Arthritis**

---

Recommended Rollator:

- Alpha Basic/Alpha Advanced with Forearm Supports, & models that can be fitted with Anatomical Handles: Symphony, Legacy, Maxi+, Jazz

## **Osteoporosis**

---

Recommended Rollator:

- All models, depending on user's size and daily mobility requirements

## **Parkinson's**

---

Recommended Rollator:

- All models that can be fitted with Slow-Down Brakes: Symphony, Legacy, Maxi+, Alpha

# Dolomite® Jazz Rollator



reddot design award  
winner 2013



Model no. JAZZLOW

## Dolomite Jazz Rollator

### Basket, anatomical grips, curb climbers and backrest standard.

Model no. JAZZSTD            JAZZ Standard  
Model no. JAZZLOW         JAZZ Low  
Model no. JAZZLOWHH     JAZZ Low with High Handles  
Colour: Anthracite Grey

The Dolomite Jazz is the epitome of function and style. Jazz provides superior ergonomics and control, making it the perfect companion to an active lifestyle. The rigid and durable frame can be folded effortlessly through the unique gate hinge and locked in place for storage or travel. Carefully crafted to provide the best in safety and mobility, Jazz boasts several patents and was awarded the red dot award: product design by an international jury of experts.

## Features

- Simple to operate, one-touch height adjustment - Patented release mechanism is contoured and provides effortless height adjustments requiring no tools
- Memory locking handle adjustment - Intelligent handle height memory function is a Dolomite-patented design ensuring the user's handle height is always returned to the proper setting
- Hybrid non-linear braking - Delivers an even braking force remaining the same at both the beginning and at the end of the braking cycle. Hybrid system provides a clean look and simple adjustments
- Comfortable, waterproof seat - Provides additional comfort while using Jazz in the seated position. Seat pads are easily removed for cleaning
- Convenient folding basket - Provides space for daily necessities and holds up to 15 lbs. Built-in carry handles allow the basket to be removed and used as a tote

## Specifications

	Low	LowHH	Standard
Distance between rear wheels	20"	20"	20"
Minimum handle height	26 ½"	30.3"	29 ½"
Maximum handle height	31.5"	35.4"	39"
Seat height	20"	20"	24"
Seat dimensions	18" x 11"	18" x 11"	18" x 11"
Width folded	8 ½"	8 ½"	8 ½"
Overall width	23 ½"	23 ½"	23 ½"
Overall length	27"	27"	27"
Width between handles	18"	18"	18"
Rollator weight (lbs.)	16	16.8	17
Weight capacity (lbs.)	331	331	331

# Dolomite® Legacy Rollators



Legacy Deluxe with padded seat and backrest

Model no. 12050-38-86

## Dolomite Legacy Rollator

### Basket and adjustable backrest standard.

Model no. 12050-38-86	LEGACY Standard Candy Apple Red
Model no. 12050-37-86	LEGACY Standard Navy Blue
Model no. 12052-38-86	LEGACY Low Candy Apple Red
Model no. 12052-37-86	LEGACY Low Navy Blue
Model no. 12054-38-86	LEGACY Super Low Candy Apple Red
Model no. 12054-37-86	LEGACY Super Low Navy Blue

## Dolomite Legacy Deluxe Rollator

### Padded backrest, padded seat and basket standard.

Model no. 12050-38-86D	LEGACY Deluxe Standard Candy Apple Red
Model no. 12050-37-86D	LEGACY Deluxe Standard Navy Blue
Model no. 12052-38-86D	LEGACY Deluxe Low Candy Apple Red
Model no. 12052-37-86D	LEGACY Deluxe Low Navy Blue
Model no. 12054-38-86DLX	LEGACY Deluxe Super Low Candy Apple Red
Model no. 12054-37-86D	LEGACY Deluxe Super Low Navy Blue

The Dolomite Legacy is more than just a name, it is a Legacy. For over 10 years the Legacy has brought daily mobility back to hundreds of thousands throughout 80 countries worldwide. Founded on Dolomite's patented flip-up seat technology, Legacy provides proven safety and support, allowing you to regain your independence and enjoy all of life's daily activities with confidence.

## Features

- Flip-up seat ensures proper support inside walking frame and reach access
- Unique forward-facing handles promote proper walking posture
- Patented brake system with adjustable tension for precise control
- Large 8" SLII wheels with superior bearings for increased maneuverability
- Optional tall handles increase handle height up to 39"

## Specifications

	Super Low	Low	Standard
Distance between rear wheels	19"	19"	19"
Minimum handle height	28"	30 ½"	33 ½"
Maximum handle height	31"	33 ½"	36 ½"
Min. & max. tall handle heights	N/A	34 - 37"	35 - 39"
Seat height	19"	21 ½"	24"
Overall width	23"	23"	23"
Overall length	23"	26"	26"
Width between handles	17"	17"	17"
Rollator weight (lbs.)	15	15	15
Weight capacity (lbs.)	275	275	275

# Dolomite® Symphony Rollator



Model no. 12160

## Dolomite Symphony Rollator

### Basket and back support standard.

Model no. 12160	SYMPHONY Standard
Model no. 12161	SYMPHONY Low
Model no. 12162	SYMPHONY Super Low
Colour: Silver	

Made of ultra-light aluminum, the Dolomite Symphony represents the ultimate in safety, comfort and ergonomic design. Finished in an elegant silver colour, the Dolomite Symphony combines style and sophistication with the weight advantage of aluminum.

### Features

- Lightweight aluminum frame
- Dolomite Advanced Brake System (DABS II)
- Flip-up seat for greater reach access and gait training
- Adjustable back support provides comfort and safety
- Infinite handle height and angle adjustments
- Anatomical handles help prevent pain and numbness in the hands

### Specifications

	Super Low	Low	Standard
Distance between rear wheels	18 ½"	18 ½"	18 ½"
Minimum handle height	25"	27 ¾"	30"
Maximum handle height	30"	33 ¾"	36"
Min. & max. tall handle heights	N/A	N/A	N/A
Seat height	18"	21 ½"	24 ½"
Overall width	23"	23"	23"
Overall length	22"	24 ½"	27 ½"
Width between handles	17"	17"	17"
Rollator weight (lbs.)	14	14	14 ½
Weight capacity (lbs.)	275	275	275

# Dolomite® Maxi+ Rollator



Model no. 12120-37-746

## Dolomite Maxi+ Rollator

### Basket and back support standard.

Model no. 12120-37-746      MAXI+ Standard  
Model no. 12122-37-746      MAXI+ Low  
Colour: Navy Blue

Dolomite's Maxi+ rollator offers high quality and function for users requiring additional seating space and increased weight capacity. The Maxi+ has a reinforced frame with a maximum weight capacity of 440 lbs.

## Features

- Forward swept handles to promote better posture
- Flip-up seat for greater reach access and gait training
- Anatomical handles
- Patented braking system
- Extra wide seating area
- Large 8" SLII wheels with superior bearings for increased maneuverability

## Specifications

	Low	Standard
Distance between rear wheels	22"	22"
Minimum handle height	29"	32"
Maximum handle height	35"	38 ½"
Seat height	22"	25 ½"
Overall width	27 ¼"	27 ¼"
Overall length	25 ½"	29"
Width between handles	21 ½"	21 ½"
Rollator weight (lbs.)	22	22
Weight capacity (lbs.)	440	440

# Dolomite® Soprano Rollator



## Dolomite Soprano Rollator

### Basket and back support standard.

Model no. 12030-44-86      SOPRANO Standard

Colour: Champagne

The Dolomite Soprano is a heavy-duty walker with large 10" wheels, providing excellent performance on rough, uneven ground, making it ideal for active outdoor users. With a cable-free braking system, a wide, sturdy frame and reinforced seat, the Dolomite Soprano allows for an increased weight capacity of 355 lbs.

## Features

- Cable-free braking system
- Frame – Sturdy heavy-duty frame provides excellent stability
- Adjustable handles – Infinite adjustments to the required height and adjustable inward angle up to a maximum of 15°
- Hand grips – Gently wedge-shaped hand grips provide a comfortable hand position for all hand sizes
- Flip-up seat – Extra large seat flips up for greater reach access
- Folding – Stands when folded for easy storage and transport
- 10" wheels – Large solid wheels provide excellent traction and a comfortable ride

## Specifications

	Soprano
Distance between rear wheels	22"
Minimum handle height	31"
Maximum handle height	39"
Seat height	24"
Overall width	25 ½"
Overall length	30"
Width between handles	18"
Folded in upright position	19"
Folded in transport position	13"
Wheel diameter	10"
Rollator weight (lbs.)	21
Weight capacity (lbs.)	355

# Dolomite® Alpha Rollators



## Dolomite Alpha Advanced/Basic Rollators

### Basket and food & beverage tray standard.

Model no. 12220-35-01 ALPHA Advanced with forearm attachment

Model no. 12230-35-01 ALPHA Basic with forearm attachment

Model no. 12220-02-35-87 ALPHA Advanced with platform attachment

Model no. 12230-02-35-87 ALPHA Basic with platform attachment

Colour: Navy Blue

The Dolomite Alpha is a special rollator designed specifically for users who deal with the added complexity of walking with arthritic or stroke conditions. The unique, adjustable forearm supports or platform option allow for maximum stability. Alpha also offers two braking systems to accommodate the specific needs of the user.

## Features

- Frame construction maximizes stability and encourages proper posture and balance
- Forearm supports adjustments in angle, width and depth
- Individually adaptable hand grip can be adjusted in height and angle
- Unique advanced brake grip is extremely easy to use and accommodates limited hand dexterity
- Flip-up seat for greater reach access and gait training
- Dolomite Advanced Brake System (DABS II)

## Specifications

	Forearm/Platform
Distance between rear wheels	22"
Minimum handle height	37 ½"/38 ½"
Maximum handle height	47 ½"/48 ½"
Seat height	24 ½"
Overall width	26"
Overall length	27 ½"
Width between handles	14 - 17"
Rollator weight (lbs.)	28 ½"/29 ½"
Weight capacity (lbs.)	275



## Dolomite Rollator Warranty

All Dolomite Rollator frames have a lifetime warranty. Brake Handle and Housings are covered for a period of five years and all other components for a period of one year. Warranty does not cover normal wear.



## Dolomite Rollator Accessories

An extensive selection of unique, functional accessories to meet many different needs.



### Deluxe Basket

*These stylish baskets with built-in carrying handles and inside zippered pocket are ideal for storing your belongings. 10 kg./22 lb. capacity.*

Model no. 1560336, Grey  
Model no. 1564104, Red  
(Fits Jazz models after SN 13CHX000681)



### Padded Backrest

*Fits over Jazz 1522062 standard backrest to provide additional comfort while seated.*

Model no. 1551654 (Jazz After SN13CHX000681)

# Dolomite® Rollator Accessories



## Padded Back Support

*Made of soft velvet for more comfort when seated.*

Model no. 12514 (Legacy, Soprano, Maxi+, Symphony)



## Cane Holder

*Can be easily mounted on the front or rear of any of the rollators for use with any single point cane.*

Model no. 12606 (Legacy, Soprano, Maxi+, Symphony, Alpha)

Model no. 1512912 (Jazz After SN13CHX000681)



## Oxygen Tank Holder

*To meet respiratory equipment needs, this durable holder is suitable for "D" sized oxygen tanks.*

(Not available on Super Low models)

Model no. 12682-25 (Symphony, Legacy, Maxi+, Alpha)

Model no. 12683 (Soprano)



## Zippered Bag

*The water-resistant, zippered bag attaches easily to the frame beneath the seat providing extra security for personal belongings.*

Model no. 12523 (Legacy, Soprano, Maxi+, Symphony, Alpha)



## Food and Beverage Tray

*Ideal for in-home use, these affordable and easy-to-clean trays attach simply to the rollator seat. 15.75" x 12" (70311).*

*15.5" x 10.25" (1546264).*

Model no. 70311 (Symphony, Legacy, Maxi+, Alpha [std.], Soprano)

Model no. 1546264 (Jazz After SN13CHX000681)



## Alpha Platform

*Soft, comfortable platform attachment is ideal for increased support and gait training.*

Model no. 13018



## Curb Climber


*The curb climber makes it easier to negotiate over obstacles by lifting the front of the rollator over short heights. Can be mounted on right or left side.*

Model no. 12562-21 (Symphony, Legacy)

Model no. 12560-21 (Maxi+, Alpha)

Model no. 1563251 (Jazz after SN 13CHX000681)

# Dolomite® Rollator Accessories



13041-21

12543KIT

## Slow Down Brake

*Unique design allows user to control speed for increased safety and gait training without wearing on the wheel. Ideal for Parkinson's or where cognition is an issue.*

Model no. 13041-21 (Symphony, Maxi+, Alpha, Jazz (before SN13CHX000680) & Legacy Rollators)  
Model no. 12543KIT (Soprano)



## Hemiplegic (HP) Handle

*The foam handle can be angled to suit individual needs to promote a safe grip on the weak side. Brake handle can be rotated for comfort. Can be mounted on right or left.*

Model no. 12611-L-60, 12611-R-60 (Symphony)



## Anatomical Hand Grips

*The larger handle surface reduces fatigue through improved weight distribution and prevents pain and numbness in the hands through increased circulation. Ideal for arthritic conditions.*

Model no. 12509 (Symphony, Legacy, Maxi+)  
Model no. 1563282 (Jazz after SN 13CHX000681)



## Tall Handles

*These chrome-plated handles with durable foam grips improve rollator fit and comfort giving extra handle height for those who need it.*

Model no. 12621 (Legacy)



## One-Hand Brake

*If hand function is impaired on one side, the one hand brake allows for braking or locking of both brakes from the user's strong side. Ideal for arthritic or stroke clients.*

Model no. 12505-L, 12505-R (Legacy, Maxi+)  
Model no. 12514-L, 12514-R (Alpha Basic)  
Model no. 12515 (Alpha Advanced)  
Model no. 12507-R or 12507-L (Symphony)

## Accessory Chart

	Soprano	Symphony	Legacy / Legacy Deluxe	Jazz	Alpha	Maxi +
Slow Down Brake		•	•		•	•
One-Hand Brake		•	•		•	•
Tall Handles			*•			
Anatomical Hand Grips		•	•	•		•
HP Handle		•				
Oxygen Tank Holder		*•	*•		•	•
Food & Beverage Tray	•	•	•	•	•	•
Cane Holder	•	•	•	•	•	•
Bag for Basket		•	•		•	•
Padded Back Support/Back Rest			•	•		•
Padded Seat Cushion		•	•		•	
Deluxe Basket				•		

\* Not available in Super Low.



**Yes, you can.®**

**Invacare Canada LP**

[www.invacare.ca](http://www.invacare.ca)

570 Matheson Blvd. E.  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

All rights reserved.  
Trademarks are identified  
by the symbols <sup>TM</sup>, <sup>SM</sup> and <sup>®</sup>.  
All trademarks are owned  
by or licensed to Invacare  
Corporation unless  
otherwise noted.

Specifications are subject to  
change without notification.

© 2018  
Invacare Canada  
Form No. 05-430C  
Rev 6/18