

INVACARE® ETUDE HC HOMECARE BED





Etude HC Homecare Bed

The ETUDE HC Homecare bed, with its attractive aesthetics, flexible design, simple operation and easy handling is the ideal choice for community home care.

The ETUDE HC is the only advanced homecare bed to offer Shear Management Technology as a standard feature, reducing shear and friction providing superior comfort and safety. The SMT feature is provided through the auto-regression of the head section while repositioning in ETUDE HC through use of the 8 function hand control.

Features and Options













- 1. Built in Dual Height Brackets.
- 2. Shear Management Technology.
- 3. Transport and Storage Kit.
- 4. Support Handle.
- 5. Lifting Support Trapeze Kit.
- 6. 4" and 8" Extension Kit.

Modern and Flexible

ETUDE HC delivers a modern, streamlined and discreet design perfectly complimenting any interior. The ETUDE HC can be configured with attractive wood side panels and a choice of safety rails to accommodate specific transfer, safety and support requirements. ETUDE HC provides a deck height ranging from 7.8" to 26.5", the best in the industry. The simple to use 8 function hand control enables the user or caregiver to lock out specific functions for increased safety.

Easy Handling & Storage

ETUDE HC recognizes that portability, storage and set-up are important features for users and caregivers alike. The simplistic design allows for simple assembly without tools and transport can be managed safely and comfortably as each section of ETUDE HC weighs less than 52 lbs. Smooth roll 3" locking casters enable the bed to be easily manoeuvered within the home.

Highest Safety Levels

The ETUDE HC conforms to the IEC bed standard IEC60601-2-52.

This equates to the highest degree of safety for both the client and the caregiver. The 8 function pendant enables the caregiver to lock out any of the beds functions for enhanced patient safety. ETUDE HC's selection of side rail supports have all met the demands for strength and entrapment guidelines required in the new safety standards for beds. Visit www.invacare.ca for more information.

A





Length outside: 86.5" Length inside: 82"



Height adjustable: 7.8" - 26.5"



Mattress support dimensions: 36" - 80"



0 - 70° Backrest angle



0 - 18° Thigh angle



0 - 25° Legrest lift



11° Tilt option*



Product weight: 166 lbs.



Max. Bed weight capacity: 440 lbs. Max. User weight capacity: 360 lbs.

Frame colour and wood work







Beech

Warranty

5 years on frame, 3 years on actuators, 2 years on pendant and electronics.

Warranty against manufacturing defects and does not include normal wear, misuse or abuse. Custom modifications to ETUDE HC void manufacture Warranty.



Invacare Canada LP www.invacare.ca

570 Matheson Blvd. E.,

Mississauga, Ontario L4Z 4G4 Canada (800)

Invacare Corporation. All

to change without notice.

Form No. 15-912C

Rev: 4/19

Unit 8

668-5324

rights reserved.

Trademarks are identified by the symbols ™ and ®.

All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Specifications are subject

©2019

INVACARE®

ETUDE HC








Supporting your needs. The Etude HC homecare bed features a wide range of options and accessories providing you the opportunity to customize the positioning, transfer assistance and design of your bed. These options makes the Invacare Etude HC one of the most unique and versatile homecare beds available today.

















1. Verso Side Support Rails 15x59"

Model # EVSR-1823

The Verso Side Support Rail is a ¾ length rail that ascends the ETUDE HC to a new level of patient care. The Verso Side Support Rails provide a unique collapsable, folding solution, accommodating the needs of those at a higher risk of falls or requiring enhanced care.

2. Lifting Support Trapeze Kit

Model # ELSP-2093

The Lifting Support Trapeze Kit provides safe and convenient overhead repositioning assistance while in bed. This kit is exclusive to the ETUDE HC and simply and securely mounts into the receivers of the head section. It can be mounted completely tool-free on either the left or right hand side.

3. Side Support Rail 16x12 & 16x16

Model # ESR-2477 & ESR 2478

The 16x12 and 16x16 Side Support Rails deliver safe and simple transfers. These $\frac{1}{2}$ and $\frac{1}{2}$ rail inspired designs provide the necessary accommodations for safe positioning and transfer support, while ensuring user safety by meeting global bed standard entrapment guidelines.

4. Side Swivel Transfer Supports

ESSHL-2479 - LEFT & ESSHR-2480 - RIGHT

The ETUDE HC Side Swivel Transfer Supports are a unique solution for users in search of a low profile, discreet transfer aid. The swivel handles remain below the surface of the mattress and are easily swung out and locked into position to provide necessary transfer support on the right or left hand side.

5. Beech Accent Panels

Model # ESP-2716

The attractive Beech Accent Panels fit easily over the ETUDE HC side frame. The ESP-2716 adds a designer expression providing subtle contrasts enhancing the ETUDE HC's modern style, delivering a versatile and appealing addition to any décor. This option cannot be used in conjunction with the EVSR-1823 Verso Side Support

6. Transportation and Storage Adaptor Kit

Model # ETA-3779

Storage and delivery management of the ETUDE HC is made simple, safe and efficient with the use of the Transportation and Storage Adaptor kit. The ETA-3779 allows the ETUDE HC to be stored together in one compact mobile unit using the beds existing casters, enhancing the ease of storage and transport and vastly reducing the opportunity for component loss and damages.

7. Deck Extension Kit

Model # EEK-3993 - EEK-3994

The Deck Extension Kit adds an additional 4" or 8" of length to the ETUDE HC to accommodate applications where a deck length of greater than 80" is required.

8. Emergency Battery back-up kit

Model # EBBK-4186

Add an additional level of safety, comfort and piece of mind with the EBBK-4186. The 12V 1.3 AH rechargeable battery provides a temporary power supply enabling the user or caregiver to safely re-position the bed in emergency situations where 110V AC becomes unavailable. Up to four cycles can be achieved with a full charge.

Invacare Canada LP

www.invacare.ca

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

©2019

Invacare Corporation. All rights reserved.

Trademarks are identified by the symbols ™ and ®.

All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Specifications are subject to change without notice.

Form No. 15-916C Rev: 4/19

