INVACARE®

FULL-ELECTRIC BED





The **Invacare® Full-Electric Bed** offers the greatest convenience for the patient and caregiver. The easy-to-use pendant control provides motorized positioning of the upper body and knees, and also adjusts bed frame height for efficient use of accessories and for safe and easy transfer.

١

Features

High impact bed end panels are durable, impact and scratch resistant; washable for easy cleaning

Universal bed ends can be used as either the headboard or footboard on the bed

Rubber caster tread for better tracking with non-marring material

DC low voltage, compact motor system Quiet motor operation and fast cycle time on lowering height of bed for efficient bed set-up

Motor system with waterproof enclosures is rated IPX6. Washable with standard water pressure. The bed must not be power washed

Color-coded connections on junction box for easier maintenance


Double-insulated junction box helps provide extra protection against electrical shock

Ergonomic hand pendant comes with durable strain relief and heavy-duty cord. Waterproof casing for easy cleaning, conveniently attaches to bed sheet with clip or rail with hand pendant cradle (sold separately)

Interchangeability simplifies storage, reduces inventory complexity











- 1. Patented universal bed ends
- 2. Ergonomic hand pendant
- 3. Color coded connections
- 4. Emergency manual crank
- 5. Thermal overload protection

- 6. Double insulated electrical components/wiring











7" W 24" D







Model #	Dimensions	Sleep Surface	Weight Capacity	Shipping Weight	Carton Dimensions	Regulation Compliance	Limited Warranty
5410IVC	88" L 15"-23" H 36" W	36" W	Safe Working Load: 450 lb.* lax. Patient Weight 350 lb.		N/A	UL962 CSA Certified	5 Years Welds 2 Years Mechanical/ Electrical
5000IVC Headspring				45 lb.	43" L 8" W 36" D		5 Years Welds 2 Years Mechanical/ Electrical
5490IVC Footspring				81 lb.	47" L 9" W 36" D		5 Years Welds 2 Years Mechanical/ Electrical
5301IVC				55 lb.	39" L		



Invacare Corporation

www.invacare.com

One Invacare Way Elvria, Ohio

570 Matheson Blvd. E.,

Mississauga, Ontario L4Z 4G4 Canada (800) 668 - 5324

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 17-987C Rev. 03/17

44035-4190 (800) 333-6900 Canada

USA

Unit 8

@2017

Universal Bed End

^{*}Includes patient, mattress and all accessories.