



INVACARE® MICROAIR® MA900
THERAPEUTIC SUPPORT SURFACE



Yes, you can.®

microAIR®

Invacare® microAIR® MA900

Lateral Rotation Low Air Loss Mattress System

Medium - High Risk
Stage I-IV

Alternating Lateral Rotation and True Low Air Loss help provide Pressure Relief and Moisture Management, with built-in side bolsters for added safety



Highly vapor permeable and oversized pliable quilted nylon top cover provides low shear, friction and moisture protection



Cable management system keeps cords away from the ground to help prevent trip hazards





Nylon base cover help to reduce shear and friction during bed articulation



Features & Benefits

10 adjustable comfort settings depending on the patient's therapy requirement

Auto Firm Mode provides quick maximum inflation within seconds to aid with transfers and nursing procedures



Features **continuous 20° and 40° lateral rotation therapy**, which can help lower the risk of infection, pneumonia and other pulmonary complications in at-risk patients

True Low Air Loss Therapy provides extra comfort for moisture and temperature control

Customizable cycle times to meet patient's varying needs



Foam base helps prevent bottoming out in the event of a power failure

Soft nylon urethane backed cells provide easy cleaning and helps minimize the proliferation of bacterial growth

Inflatable built-in side bolsters with manual connectors can be removed as needed to aid with ease of ingress and egress

External exhaust keeps the mattress surface cooler and the pump quieter

Request a Quote or Place an Order Today!



Contact your Inside Sales Representative

or



Contact Customer Service at 800.333.6900

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2019

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ©. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form No. 19-028
190207



Invacare microAIR MA900 Lateral Rotation Low Air Loss Mattress System

| Model # | Dimensions | Weight Capacity | Flow Rate (Pump) | Therapy Type | Alarms | Limited Warranty |
|----------|--|--|------------------|---------------------------------------|----------------|--|
| MA900 | 80" L 36", 42", 48" W 10" H | 36" 600 lbs. 42", 48" 1000 lbs. | 1275 LPM | Lateral Rotation True Low Air Loss | Audio & Visual | 1 Year (Cover/Mattress) 1 Year (Power Unit) |
| MA900B42 | microAIR MA900 42" Lateral Rotation Low Air Loss Mattress System | | | | | |
| MA900B48 | microAIR MA900 48" Lateral Rotation Low Air Loss Mattress System | | | | | |
| MA900M | microAIR MA900 36" Mattress only | | | | | |
| MA900M42 | microAIR MA900 42" Mattress only | | | | | |
| MA900M48 | microAIR MA900 48" Mattress only | | | | | |
| MA900P | microAIR MA900 Pump only | | | | | |