


INVACARE®  
**MARINER™**  
REHAB SHOWER CHAIR

---



---

The **Invacare® Mariner™ Rehab Shower Commode Chair** offers rust-resistant aluminum frames and stainless steel hardware, which make them ideal for use in the shower. In addition, these new chairs fit over most standard and elongated bowls for over the toilet use as well. They feature flip-back padded arms that lock down for safety, swing-away front riggings with tool-free adjustable height footrests and a four-position padded seat. The Invacare Mariner Commode Chair offers various seat widths and wheel options for a custom fit along with additional features and benefits.

---

## Features

Aluminum frame and stainless steel hardware  
 Fits over most standard and elongated toilet bowls  
 Padded four-position seat  
 Flip-back padded arms that lock down for safety  
 Swing-away front riggings

Tool-less adjustable height footrests  
 Folds into a compact unit for storage and transport  
 Removable pan and pan holder  
 Safety strap included  
 23" treaded urethane tire and 5" caster models



1. Padded four-position seat
2. Swing-away front riggings with tool-less adjustable height footrests
3. Flip-back padded arms that lock down for safety



6895: 23" treaded urethane tires with 18" seat

**Invacare Corporation**  
 www.invacare.com

**USA**  
 One Invacare Way  
 Elyria, Ohio  
 44035-4190  
 (800) 333-6900

**Canada**  
 570 Matheson Blvd. E.,  
 Unit 8  
 Mississauga, Ontario  
 L4Z 4G4 Canada  
 (800) 668-5324

©2012  
 Invacare Corporation.  
 All rights reserved.  
 Trademarks are identified  
 by the symbols ™ and  
 ®. All trademarks are  
 owned by or licensed  
 to Invacare Corporation  
 unless otherwise noted.  
 Specifications are subject  
 to change without  
 notification.  
 Form No. 03-101  
 Rev. 09/22  
 111278



Model #	Dimensions HxWxD	Width Between Arms	Width Folded	Seat Dimensions	Seat-to-Floor Height	Seat-to-Arm Height	Back Dimensions
6895	39"x26.5"x32"	18.5"	12"	18.25"x18.25"	22.75"	8"	15.5"Hx18"W



	Footrest Extensions	Clearance Over Toilet HxWxD	Weight	Weight Capacity	Units Per Carton	Limited Warranty
6895	13"-18"	20.25"x16"x26.5"	41 lb.	300 lb.	1	3 Years

