



SMOOV

Suggested Canadian Price List
Effective April 1, 2022

PROVIDER NAME:	SHIPPING ADDRESS:
ACCOUNT NUMBER:	CITY:
PO NUMBER:	PROVINCE:
CONTACT NAME:	POSTAL CODE:
TAG NAME:	PHONE:

SMOOV Wheel System

MSRP

Each SMOOV system includes: 1 SMOOV drive unit with integrated battery, 1 Control unit w/ Bracket, 1 Battery charger, 1 USB-C charger cable, Rigid or Folding chair Adapter Axle / Bracket Set.

Fast, swift and low- maintenance, 10km/h - maximum speed. Quick assembly. Easy to dock and undock. Easy to transport, weight of just 16lbs.

The SMOOV Mobility App: Connect the SMOOV with your smartphone, make settings, keep battery charge, and range always in view.

- 1592897 SMOOV one O10, 6km / h
- 1592898 SMOOV one O10, 10km / h

\$6,495
\$6,995



WARRANTY: THE SMOOV IS COVERED FOR 24 MONTHS AGAINST MANUFACTURING DEFECT.
WEARABLE ITEMS ARE COVERED FOR A PERIOD SIX MONTHS AND BATTERIES ARE COVERED FOR ONE YEAR.

Wheelchair Brackets

* MUST SELECT ADAPTER AXLE OR BRACKET, MUST SPECIFY WHEELCHAIR DETAIL FOR CORRECT BRACKETS AND AXLES.

Note: Adapter Axle are only compatible to 0° - 2°
*ADAPTER AXLES - Folding Chairs

- 1593053 Adapter axle Folding frame 28.5 - 30.5 cm
- 1593054 Adapter axle Folding frame 30.5 - 32.5 cm
- 1593058 Adapter axle Folding frame 32.5 - 34.5 cm
- 1593055 Adapter axle Folding frame 34.5 - 36.5 cm
- 1593056 Adapter axle Folding frame 36.5 - 38.5 cm
- 1593057 Adapter axle Folding frame 38.5 - 40.5 cm
- 1593970 Adapter axle Folding frame 40.5 - 42.5 cm



\$305
\$357
\$305
\$357
\$335
\$357
\$357

OR

*BRACKET SETS - Rigid Frame

- 1593073 Bracket set axle 30mm
- 1593076 Bracket set axle 32mm
- 1593074 Bracket set axle 35mm
- 1593075 Bracket set axle 40mm
- 1593077 Bracket set axle 1 1/4"



\$229
\$229
\$226
\$226
\$229

NOTE: All specifications and dimensions are approximate
For the most current pricing information visit www.invacare.ca

Wheelchair Brackets Cont'd

MUST SPECIFY WHEELCHAIR DETAIL FOR CORRECT BRACKETS AND AXLES

*BRACKET SETS - CONTROL UNIT

<input type="checkbox"/>	1593029	Bracket set control unit 23mm
<input type="checkbox"/>	1593035	Bracket set control unit 25mm
<input type="checkbox"/>	1593036	Bracket set control unit profil 1
<input type="checkbox"/>	1593037	Bracket set control unit profil 2
<input type="checkbox"/>	1593038	Bracket set control unit 30mm
<input type="checkbox"/>	1593278	Bracket set control unit 32mm



\$339.00
\$344.00
\$339.00
\$339.00
\$400.00
\$400.00

Note: Brackets are no charge with every Smoov order.

Wheelchair Specifications (Required)

Mounting by Invacare (Optional):


Provider is responsible for shipping to and from Invacare. TS number will be supplied by Invacare Tech Services.
Cost of installation will vary based on chair configuration.

TS# _____


MANUFACTURER: _____

MODEL: _____


YEAR OF MANUFACTURE: _____



STF HEIGHT: *(Rear, sling bottom edge measurement)*



SEAT WIDTH



DISTANCE: *Stub axle to stub axle (between)*

Invacare Canada reserves the right to change prices and specifications without notice.

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324 (Canada) Fax: (905) 890-5244 (800) 668-5478 (Canada)
Invacare®, and Yes, you can., are trademarks of Invacare Corporation.

WWW.INVACARE.CA

NOTE: All specifications and dimensions are approximate
For the most current pricing information visit www.invacare.ca