

INVACARE®
STORM TORQUE®-LE
REAR WHEEL DRIVE POWER WHEELCHAIR



INVACARE® STORM® TORQUE-LE

The Invacare Storm Torque-LE is Invacare's rear-wheel drive rear-wheel drive power wheelchair. Equipped with Performance 4-Pole motors, a 6.4 mph top speed and outfitted with Maxx adjustable rehab seating, the Torque provides superior positioning, comfort and function.

Torque's optimal flexibility and function is enhanced with Invacare MK6i™ Electronics and Rubber Element standard rear spring suspension or optional Monroe® Spring Shock Suspension. Aluminum mag wheels and a wrap-around battery enclosure provide sporty and rugged good looks.

Features

Maxx ADJASBA seat is adjustable in width and depth.
 Rubber Element standard rear-wheel suspension helps absorb significant jolts and impacts for a smoother ride. Also available with optional Monroe® independent rear suspension.
 Equipped with wheelchair transport brackets (TRBKTS).
 Manufactured in North America.

Invacare Core Technologies

Invacare® High-Performance 4-pole motor package
 Invacare G-Trac™ Technology provides precise tracking for imprecise environments
 Invacare MK6i™ Electronics are the benchmark in innovative technology
 Can be interfaced with Motion Concepts power positioning systems

TORQUE-LE Features



- 1) Secure armrest attachment provides stability and enables easy removal of the flip-back arm.
- 2) Invacare® G-Trac™ Technology improve tracking making it possible for the power wheelchair user to drive straight on an intended path with fewer veer corrections, regardless of the terrain or obstacles that cause the chair to turn toward one side or the other.
- 3) The Motion Concepts Maxx Rehab and Modular Maxx Powered Positioning Systems provide maximum function and adjustment.
- 4) Invacare® MK6i™ Technology

*The Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of restraints available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.



The Torque SP can be outfitted with Invacare Matrix® seating and positioning products to provide optimal comfort and support. The extensive range of sizes and hardware options provide solutions and adjustability for a variety of individuals.

Invacare Canada LP
www.invacare.ca

570 Matheson Blvd. E.,
 Unit 8
 Mississauga, Ontario
 L4Z 4G4 Canada
 (800) 668-5324

16773, boulevard Hymus
 Kirkland, Quebec
 H9H 3L4
 (800) 363-0137

©2016
 Invacare Corporation.
 All rights reserved.
 Trademarks are identified by the symbols ™ and ®.
 All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
 Specifications are subject to change without

Form No. 15-919C
 Rev. 12/16



Model #	Overall Dimensions	Rehab Seat Dimensions	Rehab Seat-to-Floor Height	Drive Wheels	Ground Clearance	Turning Radius	Incline
3GTQSP-LE	43.7" L 26.2" W 40.4" H	12" - 24" x 12" - 25"	17.5" - 19.5"	14"	3"	33.1"	9°
	Total Weight	Weight Capacity	Motor	Max Speed	Battery	Battery Range	Limited Warranty
3GTQSP-LE	Up to 364 lb.	300 lb.	4 Pole	6.4 mph	Group 24 22 NF	Up To 19.36 Miles	5 Years Base Frame 3 Years Seat Frame 1 Year Electronics/ Motor/Gearbox