



## Thank you for using Motion Concepts eForms.

Motion Concepts is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.



### Adobe Acrobat Reader DC

\* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

#### Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

**SAVE**

#### Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

**SUBMIT**

**CLEAR FORM**

#### Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

**905-695-0583**

**PRINT FORM**



**TDX-SP2 HD BASE WITH  
MODULAR MAXX POWER POSITIONING SYSTEM**  
CDN PRICE LIST AND ORDER FORM  
Price Effective August 12, 2022



PRINT FORM

SAVE

CLEAR FORM

**REQUIRED INFORMATION**

Quote  Order

Date of Order: \_\_\_\_\_ Dealer Account #: \_\_\_\_\_

Dealer Name: \_\_\_\_\_

PO #: \_\_\_\_\_ Tag \_\_\_\_\_

Purchasing Contact: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

E-mail: \_\_\_\_\_

RTS/Therapist: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

City: \_\_\_\_\_ Province: \_\_\_\_\_


Postal Code: \_\_\_\_\_ Client Gender: M F

Special Client Conditions: \_\_\_\_\_

**REQUIRED MEASUREMENTS**

To ensure system is accurately configured please fill in all required

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Cushion Thickness:	



**Customer Service: 1.866.748-7943 | Fax: 1.905.695.0583 | www.motionconcepts.com**

**INVACARE POWER BASES**

<input type="checkbox"/> TDXSP2HD	TDX-SP2 HD Base .....	\$8,275
Select STF <sup>(1)</sup> :	HD-SLOW <input type="checkbox"/> 17.25" <sup>(2)</sup> .....	\$250
	HD-LOW <input type="checkbox"/> 17.75"	
	HD-MED <input type="checkbox"/> 18.75"	
	HD-TALL <input type="checkbox"/> 19.75"	

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size.
2. Must order Transport Brackets (TRBKTS). Only available on 17" seat depths and greater.

**USER WEIGHT LIMITS**

<input type="checkbox"/> U400	Weight Capacity up to 400 lbs .....	STD
-------------------------------	-------------------------------------	-----

**TRANSPORT TIE DOWN**

<input type="checkbox"/> TRBKTS	Wheelchair Transport Brackets <sup>(1)</sup> .....	STD
---------------------------------	--	-----

1. For unoccupied use only

**TIRE OPTIONS**

<input type="checkbox"/> B1431-3	14" x 3" Black Tire with Gel Foam Inserts .....	STD
<input type="checkbox"/> B1430-3T	14" x 3" Black Tire - Pneumatic Tire .....	N/C

**FORK OPTIONS**

<input type="checkbox"/> DSKF	Double-Sided Fork Package .....	STD
<input type="checkbox"/> FKPKG	Single-Sided Fork Package .....	\$375

**BATTERY TRAY TYPE**

<input type="checkbox"/> 24TRY	Group 24 Style Tray - Standard Base 25.5" .....	N/C
--------------------------------	---	-----

**FRAME FINISH**

<input type="checkbox"/> 163P		Black Ice Glossy .....	N/C
<input type="checkbox"/> 162P		Space Station Silver .....	N/C
<input type="checkbox"/> 160P		Invacare Blue .....	N/C
<input type="checkbox"/> 158P		Lights Out Black Matte .....	N/C
<input type="checkbox"/> 154P		Rockstar Red .....	N/C
<input type="checkbox"/> 155P		Tangy Orange .....	N/C
<input type="checkbox"/> 161P		Grape Jelly Bean .....	N/C
<input type="checkbox"/> 157P		Island Blue .....	N/C
<input type="checkbox"/> 156P		Monster Green .....	N/C
<input type="checkbox"/> 159P		White Out .....	N/C

**RIM INSERT COLORS - MUST PICK ONE**

<input type="checkbox"/> 163PR		Black Ice Glossy .....	N/C
<input type="checkbox"/> 162PR		Space Station Silver .....	N/C
<input type="checkbox"/> 160PR		Invacare Blue .....	N/C
<input type="checkbox"/> 158PR		Lights Out Black Matte .....	N/C
<input type="checkbox"/> 154PR		Rockstar Red .....	N/C
<input type="checkbox"/> 155PR		Tangy Orange .....	N/C
<input type="checkbox"/> 161PR		Grape Jelly Bean .....	N/C
<input type="checkbox"/> 157PR		Island Blue .....	N/C
<input type="checkbox"/> 156PR		Monster Green .....	N/C
<input type="checkbox"/> 159PR		White Out .....	N/C



**LIGHTS AND INDICATORS**

LIGHTS LED Lights and Indicators<sup>(1)</sup> ..... \$995

1. Must order REM216, REM400 or REM500

**CONTROLLERS**

EXPC Expandable Controller ..... STD

GTRAC-LX LiNX G-TRAC Module ..... \$1,200

**NON-EXPANDABLE DRIVE CONTROL**

REM110 LiNX Drive Only Remote/Joystick ..... STD

REM210 LED Non-Expandable Remote/Joystick ..... STD

Three drive function, control of single/dual actuator through drive control.

1. Lights are not available with non-expandable remote.

**EXPANDABLE DRIVE CONTROL**

REM211 LED Expandable Remote/Joystick<sup>(1)(2)</sup> ..... \$875  
Three drive functions, control of up to 6 actuators through the drive control. Not for use with specialty controls.

REM216 LED Expandable Remote/Joystick with Lights<sup>(2)(3)</sup> ..... \$875

Three drive functions, control of up to 6 actuators through the drive control and light operation. Not for use with specialty controls.

REM400 Color 3.5" Touch Screen Remote/Joystick<sup>(2)(4)(5)</sup> ..... \$1,100

SW400TGL Toggle Kit for REM400<sup>(6)</sup> ..... \$200

REM500 LiNX Color 3.5" Touch Screen Display Only - No Driver Control<sup>(2)(3)(4)</sup> ..... \$1,100

Select mounting position

R5MR Right  
 R5ML Left

PWH Harness Required for Expandable System ..... N/C

OUTPUT Output Module ..... \$645

1. Not available with lights and indicators
2. Must choose EXPC Expandable Controller and PWH Harness.
3. Must choose lighting package separately.
4. Can control lighting, must choose lighting package separately.
5. May also be used as display for specialty controls, chose specialty controls separately.
6. Adds 2 toggle switches to REM400.

**ATTENDANT DRIVE CONTROL CHOICES**

ACU Proportional Attendant Control ..... \$1,024

Select seating style for mounting type

TAM Tilt & Elevate Mount  
 RAM Recline Mount

Select mounting position

RIGHT-ACU Right  
 LEFT-ACU Left

**SPECIALTY PROPORTIONAL CONTROLS REQUIRES REM400 OR REM500**

CREM LiNX Compact Remote<sup>(1)</sup> ..... \$1,097

CREM-LF LiNX Compact Remote Low Force<sup>(1)</sup> ..... \$1,097

Select mounting position

MR Right  
 ML Left

1. Cannot be mounted on the same side as the REM400/REM500

**LINX\* DIGITAL CONTROLS SIP-N-PUFF**

PKG32666 Therafin Sip-N-Puff Breath Tube Kit ..... \$440

INPUT LiNX Input Module & Sip-N-Puff Interface<sup>(1)(2)</sup> (9 Pin Connection) ..... \$1,800

1. Only one input module per chair.

2. This is required with ASL Drive Controls.

**MOUNTING FOR DRIVER CONTROLS**

FM-ML Fixed Mounting Left ..... N/C

FM-MR Fixed Mounting Right ..... N/C

SAQR Motion Height Adjustable Swing-Away Quad Link (Right Hand) ..... \$344

SAQL Motion Height Adjustable Swing-Away Quad Link (Left Hand) ..... \$344

MRR Maxx Resolve Swing-Away Joystick Mount (Right Hand) ..... \$450

MRL Maxx Resolve Swing-Away Joystick Mount (Left Hand) ..... \$450

**LINX ELECTRONIC ACCESSORIES**

LAK LiNX Access Key<sup>(1)</sup> ..... \$50

MUSB Motion Dual USB Charger ..... N/C

1. Required for programming

**JOYSTICK TOPS**

PC101A Bodypoint U Shaped Handle 3" ..... \$120

PC102A Bodypoint U Shaped Handle 4" ..... \$120

PC107A Bodypoint Rubber Dome ..... \$120

1560 T Handle Flexible Joystick Extension ..... \$120

1561 Straight Handle Flexible Joystick Extension ..... \$120

1826 Chin Cup ..... \$120

1. To remove/disengage the joystick knob from the REM400 Remote pull straight up on the joystick (DO NOT TWIST) otherwise damage may occur and may void the warranty.

**INSTALLATION OF BATTERIES**

24BATTERY 24 Gel Battery (On Chair) ..... \$840

**BATTERY CHARGER OPTION**

110CHARGER 110 Volt Battery Charger ..... N/C

L8900 Omit Battery Charger ..... N/C

**Maxx Bariatric Power Positioning**  
on Invacare TDX-SP2 HD Base

NOTE: Systems Include Seat Pan, Back Pan, and Standard Armrests. They are designed for use with a 2" - 3" cushion

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**BARIATRIC REHAB SEAT**

MBRS Maxx Bariatric Rehab Seat . . . . . \$4,995

**BARIATRIC POWER POSITIONING SYSTEMS**

MBT Maxx Bariatric CG Tilt . . . . . \$7,995

MBTRE<sup>(1)</sup> Maxx Bariatric CG Tilt (45°) & Power Recline (150°) w/Extended Shear Reduction(ESR) . \$12,595

MBRESR Modular Maxx Power Recline (150°) w/Extended Shear Reduction(ESR) . . . . \$10,350

1. Tilt will be reduced to 45° if adding an Elevating Lift Seat Module

**ADD: HD ELEVATING LIFT SEAT MODULE**

ESM10 10" Elevating Lift Seat Module<sup>(1)(2)</sup> . . . . \$2,995  
Weight capacity up to 300lbs

1. Only Available with Tilt. Uses an integrated Tilt/Lift mechanism, tilt is 45°

2. Not available on seat sizes greater than 22" wide by 22" deep

**MODULES MUST SELECT SYSTEM TO ADD MODULES**

PRETILT 5° Fixed Pre-Tilt<sup>(1)</sup> . . . . . N/C

8PR 8° Precline Module<sup>(1)</sup> only available with Recline . . . . N/C

30PR 30° Precline Module<sup>(1)</sup> only available with Recline . . . \$500

1. Total degrees of recline is 150° ie. if 30° of precline is used, then 120° of recline will be available.

**ADULT SEAT WIDTH SETTINGS**

BW120 20" Wide (Adjustable from 20" to 24") . . . . . N/C

BW121 21" Wide (Adjustable from 20" to 24") . . . . . N/C

BW122 22" Wide (Adjustable from 20" to 24") . . . . . N/C

BW123 23" Wide (Adjustable from 20" to 24") . . . . . N/C

BW124 24" Wide (Adjustable from 20" to 24") . . . . . N/C

BW222 22" Wide (Adjustable from 22" to 26") . . . . . N/C

BW223 23" Wide (Adjustable from 22" to 26") . . . . . N/C

BW224 24" Wide (Adjustable from 22" to 26") . . . . . N/C

BW225 25" Wide (Adjustable from 22" to 26") . . . . . N/C

BW226 26" Wide (Adjustable from 22" to 26") . . . . . N/C

**ADULT ULTRA RAIL SEAT DEPTH SETTINGS**

SD117 17" Deep (Adjustable from 17" to 21") . . . . . N/C

SD118 18" Deep (Adjustable from 17" to 21") . . . . . N/C

SD119 19" Deep (Adjustable from 17" to 21") . . . . . N/C

SD120 20" Deep (Adjustable from 17" to 21") . . . . . N/C

SD121 21" Deep (Adjustable from 17" to 21") . . . . . N/C

SD219 19" Deep (Adjustable from 19" to 23") . . . . . N/C

SD220 20" Deep (Adjustable from 19" to 23") . . . . . N/C

SD221 21" Deep (Adjustable from 19" to 23") . . . . . N/C

SD222 22" Deep (Adjustable from 19" to 23") . . . . . N/C

SD223 23" Deep (Adjustable from 19" to 23") . . . . . N/C

**MATRIX SEATING OPTIONS**

Matrix Cushions available on page 10

Matrix Backs available on page 11

Matrix Headrests available on page 12



**CHEST STRAPS**

MFUM BP Monoflex, Underarm RLS M(4" x 18.25") \$195

MFUL BP Monoflex, Underarm RLS L(4.5" x 21.5") \$195

MFCM BP Monoflex, Center RLS M(3.25" x 18.25") \$195

MFCL BP Monoflex, Center RLS L(3.5" x 21.5") . . . \$195

**LAP BELTS**

BELT60 Push Button Style Lap Belt, Length 60" . . . . N/C

BELT71 Push Button Style Lap Belt, Length 71" . . . . N/C

PLB2 Bodypoint Padded Lap Belt 2 Point<sup>(1)</sup> . . . . \$250

PLB4 Bodypoint Padded Lap Belt 4 Point<sup>(1)</sup> . . . . \$315

1. 67" length only option

**MOTION CONCEPTS HEADREST OPTIONS**

MCSH Motion Concepts Standard<sup>(1)</sup> . . . . . N/C

MCOH Motion Concepts Auto Style<sup>(1)</sup> . . . . . \$235

FM Fixed Mounting (non-removable)<sup>(2)</sup> . . . . . N/C

OHR Omit: Headrest . . . . . N/C

1. Must select Fixed Mounting (FM) or Matrix Headrest Hardware on page 12.

2. Only compatible with Standard Rehab, Contour or High Back Models. All other upgradeable backs require removeable hardware.

**MATRIX HEADRESTS OPTIONS**

Matrix Elan Headrests and Hardware available on page 12



**FINISHED BACK HEIGHT<sup>(1)(2)</sup> FOR TILT ONLY/ELEVATE ONLY**

18"  19"  20"  21"  22"  23"  24"  25"  
FBR18 FBR19 FBR20 FBR21 FBR22 FBR23 FBR24 FBR25

1. 4" Gap from the seat pan to the bottom of back pan is recommended.  
2. Finished back height is measured from seat pan to top of backrest.

**FINISHED BACK HEIGHT<sup>(1)(2)</sup> FOR RECLINE OR TILT/RECLINE**

20"  21"  22"  23"  24"  25"  26"  27"  
FBR20 FBR21 FBR22 FBR23 FBR24 FBR25 FBR26 FBR27

1. 6" Gap from the seat pan to the bottom of back pan is recommended.  
2. Finished back height is measured from seat pan to top of backrest.

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### BACK CANES (NON RECLINE SYSTEMS)

<input type="checkbox"/> SBC18	18" Straight Back Canes	STD
<input type="checkbox"/> SBC20	20" Straight Back Canes	STD
<input type="checkbox"/> SBC22	22" Straight Back Canes	STD
<input type="checkbox"/> SBC24	24" Straight Back Canes	STD
Select Cane Angle	<input type="checkbox"/> 75° <input type="checkbox"/> 79° <input type="checkbox"/> 85° <input type="checkbox"/> 90° <input type="checkbox"/> 95° <input type="checkbox"/> 101° <input type="checkbox"/> 106° <input type="checkbox"/> 112° <input type="checkbox"/> 116°	
	BA75 BA79 BA85 BA90 BA95 BA101 BA106 BA112 BA116	
<input type="checkbox"/> ABC20	20" Angled (5°) Back Canes	N/C
<input type="checkbox"/> ABC22	22" Angled (5°) Back Canes	N/C
<input type="checkbox"/> ABC24	24" Angled (5°) Back Canes	N/C
Select Cane Angle	<input type="checkbox"/> 81° <input type="checkbox"/> 84° <input type="checkbox"/> 91° <input type="checkbox"/> 96° <input type="checkbox"/> 101° <input type="checkbox"/> 106° <input type="checkbox"/> 111° <input type="checkbox"/> 117° <input type="checkbox"/> 121°	
	BA81 BA84 BA91 BA96 BA101 BA106 BA111 BA117 BA121	
<input type="checkbox"/> EBC22	22" Endomorph Back Canes <sup>(1)</sup>	\$350
<input type="checkbox"/> EBC24	24" Endomorph Back Canes <sup>(1)</sup>	\$350
<input type="checkbox"/> EBCN1	1" Narrower	N/C
<input type="checkbox"/> EBCN2	2" Narrower	N/C
<input type="checkbox"/> EBCN4	4" Narrower	N/C

1. Comes standard with Cantilever Armrest Cane Mounted.

### CANE OPTIONS FOR RECLINE SYSTEMS

<input type="checkbox"/> BFS2	Recline Body Fit 2" Narrower <sup>(1)</sup>	\$995
<input type="checkbox"/> BFS4	Recline Body Fit 4" Narrower <sup>(1)</sup>	\$995
<input type="checkbox"/> AHH	Adjustable Hooking Handles <sup>(2)</sup> (pair)	\$450
<input type="checkbox"/> RPHS	Push Handles for Recline - Straight <sup>(3)</sup> (pair)	\$350
<input type="checkbox"/> RPHA	Push Handles for Recline - Angled <sup>(3)</sup> (pair)	\$350

1. Minimum finished back widths of 16".
2. Handle measures 6.5" high and has .75" height adjustment with angle adjustment.
3. 4" of total vertical mounting adjustment.

### BACK OPTIONS

<input type="checkbox"/> SRB	Standard Rehab Back & Cushion	N/C
<input type="checkbox"/> PIPF	Flat Planar Interface Plate <sup>(1)</sup>	N/C
<input type="checkbox"/> PIPR	Recessed Planar Interface Plate <sup>(1)(2)</sup>	N/C
<input type="checkbox"/> ORBC	Omit Rigid Back & Cushion <sup>(3)</sup>	N/C

1. Not available on tilt only or tilt/elevate systems, for use with Recline Module.
2. Included with all Recline Systems when no back is selected. May also be used as an interface for after market backs.
3. Not available with Recline/ESR.

### MATRIX BACKREST OPTIONS

Matrix Backs available on page 11



### ARMRESTS (NON RECLINE SYSTEMS)

<input type="checkbox"/> DP1	Adult Dual Post Adj. 9.75"-12.75"	STD
<input type="checkbox"/> DP2	Adult Dual Post Adj. 12.75"-15.75"	STD
<input type="checkbox"/> DP3	Adult Dual Post Adj. 15.75"-18.75"	STD
Select:	<input type="checkbox"/> FP Flat Pouch	STD
	<input type="checkbox"/> GB Glove Boxes	N/C
<input type="checkbox"/> OAL	ADD: OutBack Arm Option Left <sup>(1)</sup> each	\$195
<input type="checkbox"/> OAR	ADD: OutBack Arm Option Right <sup>(1)</sup> each	\$195
<input type="checkbox"/> OCTA1	Outback Cantilever Tilt Armrest 9"-12" <sup>(3)</sup>	\$488
<input type="checkbox"/> OCTA2	Outback Cantilever Titl Armrest 12"-15" <sup>(3)</sup>	\$488
<input type="checkbox"/> RMCAN	Ultra Rail Seat Mnt Flip Back Cantilever MAXX Tilt Arm (9-12.5") <sup>(2)</sup>	\$488
<input type="checkbox"/> DPM	ADD: Dual Post Module to RMCAN (pair)	\$250
<input type="checkbox"/> CAFEB	Cantilver Arms for Endomorph Back Canes (pair) (see page 4)	STD

1. Only available with Dual Post Armrests.
2. Mounts to Seat Rail. Slotted upper tube for infinite fore/aft adjustment of arm pad.
3. Not available with the Universal Cup Holder or Flip and Fold Lap Tray.

### RECLINE ARMRESTS

<input type="checkbox"/> RC1	Reclining Adjustable Height. 9.5"-13"	STD
<input type="checkbox"/> RC2	Reclining Adjustable Height. 13"-16"	STD
<input type="checkbox"/> OAL	ADD: OutBack Arm Option Left <sup>(1)</sup> each	\$195
<input type="checkbox"/> OAR	ADD: OutBack Arm Option Right <sup>(1)</sup> each	\$195
<input type="checkbox"/> OCRA1	Outback Cantilver Recline Armrest 9"-12" <sup>(2)</sup>	\$488
<input type="checkbox"/> OCRA2	Outback Cantilver Recline Armrest 12"-15" <sup>(2)</sup>	\$488
<input type="checkbox"/> MCAN1	Maxx Style Cane Mounted Height, Width and Angle Adjustable Cantilever Flip Back Armrests for Recline Systems 9.5"-13" <sup>(2)</sup>	\$488
<input type="checkbox"/> MCAN2	Maxx Style Cane Mounted Height, Width and Angle Adjustable Cantilever Flip Back Armrests for Recline Systems 12.5"-16" <sup>(2)</sup>	\$488
<input type="checkbox"/> CANQ1	Cane Mounted Height & Angle Adjustable Cantilever Flip Back Armrest 9.5"-12" <sup>(3)</sup>	\$488
<input type="checkbox"/> CANQ2	Cane Mounted Height & Angle Adjustable Cantilever Flip Back Armrest 10"-13" <sup>(3)</sup>	\$488
<input type="checkbox"/> CANQ3	Cane Mounted Height & Angle Adjustable Cantilever Flip Back Armrest 13"-16" <sup>(3)</sup>	\$488

1. Only available with Reclining Armrests.
2. Not available with the Universal Cup Holder or Flip and Fold Lap Tray.
3. Mounts to Seat Rail. Slotted upper tube for infinite fore/aft adjustment of arm pad.

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**ARMREST TUBE LENGTH - PLEASE SPECIFY**

- ATFL Left Full . . . . . STD
- ATDL Left Desk . . . . . STD
- ATFR Right Full . . . . . STD
- ATDR Right Desk . . . . . STD

**ARMPADS**

- Standard**
- SAPFL Left Full (2.25" x 14") . . . . . STD
  - SAPDL Left Desk (2.25" x 10") . . . . . STD
  - SAPFR Right Full (2.25" x 14") . . . . . STD
  - SAPDR Right Desk (2.25" x 10") . . . . . STD

- Waterfall Arm Pads<sup>(1)</sup>**
- WAPFL Left Full (13") . . . . . N/C
  - WAPDL Left Desk (9") . . . . . N/C
  - WAPFR Right Full (13") . . . . . N/C
  - WAPDR Right Desk (9") . . . . . N/C

- Modular Arm Pad - Left<sup>(2)</sup>**
- MAPTFL Left Tray Full (14") . . . . . \$95
  - MAPTDL Left Tray Desk (10") . . . . . \$95

Select Pad Insert below:

- ISSFL Iskin - Integrated soft skin foam
- SCVFL Startex Covered Visco foam

- Modular Arm Pad - Right<sup>(2)</sup>**
- MAPTFR Right Tray Full (14") . . . . . \$95
  - MAPTDR Right Tray Desk (10") . . . . . \$95

Select Pad Insert below:

- ISSFR Iskin - Integrated soft skin foam
- SCVFR Startex Covered Visco foam

- Flat Multi Position Arm Pads**
- MPAFL Left Full (14") . . . . . \$150
  - MPADL Left Desk (10") . . . . . \$150
  - MPAFR Right Full (14") . . . . . \$150
  - MPADR Right Desk (10") . . . . . \$150

- Gel Pads**
- GAPNL Narrow Left (2" x 12") . . . . . \$88
  - GAPNR Narrow Right (2" x 12") . . . . . \$88
  - GAPWL Wide Left (3.5" x 12") . . . . . \$96
  - GAPWR Wide Right (3.5" x 12") . . . . . \$96
  - GAPLL Long Left (3.5" x 14") . . . . . \$103
  - GAPLR Long Right (3.5" x 14") . . . . . \$103

- Arm Troughs**
- EATL Ergonomic Arm Trough - Left . . . . . \$121
  - EATR Ergonomic Arm Trough - Right . . . . . \$121
  - FHPL Flat Hand Pad - Left<sup>(3)</sup> . . . . . \$91
  - FHPR Flat Hand Pad - Right<sup>(3)</sup> . . . . . \$91

- Arm Pad Storage/Cup Holder<sup>(4)</sup>**
- APSL Left . . . . . \$295
  - APSR Right . . . . . \$295

**ARMPADS - CONTINUED**

- OMIT Omit Armpads(7/8" diameter tube) . . . . . N/C

1. Full length pad not compatible with the desk length armrest tube if using Ultra Rail mounted Tilt Arms.
2. Requires interchangeable pad insert.
3. Must order with ergonomic arm trough.
4. Not available on the joystick side.

**ADDITIONAL ARM ACCESSORIES**

- UCH Universal Cup Holder<sup>(1)</sup> . . . . . \$35
- MACESL Multi Axis Upper Extremity Support<sup>(2)</sup>(each)  
Left Side . . . . . \$225
- MACESR Multi Axis Upper Extremity Support<sup>(2)</sup> (each)  
Right Side . . . . . \$225
- OAS2L 1" Offset Armrest Spacer Left(each) . . . . \$125
- OAS2R 1" Offset Armrest Spacer Right(each) . . . \$125
- OAS1L 1/2" Offset Armrest Spacer Left(each) . . . \$125
- OAS1R 1/2" Offset Armrest Spacer Right(each) . . \$125
- MPAPA Multi Purpose Arm Pad  
Adapter Block<sup>(3)</sup> (each) Quantity ( ) . . . \$35

1. Not available on the joystick side.
2. Allows rotation of the armpad, only compatible with Gel Arm pad, Flat Arm pad or Ergonomic Arm Troughs.
3. Not available with Standard Arm Pads or Waterfall Arm Pads.

**ELBOW BLOCKS**

- Elbow Blocks available on page 14

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**POWER CENTER MOUNT FOOT PLATFORM**


MUST BE COMPLETED	Measurement in Inches
Knee to Heel:	
Cushion Thickness: (Subtract)	
Cushion Compression:(Add)	
Seat Pan to Footplate:(Total)	

- LNX **LNX Power Center Mount Foot Platform<sup>(1)(2)</sup> 13" - 19"** Seat pan to footplate . . . \$3,200  
With Quick Release Height Adjustment
- IFP **Comes with a rubber coated footplate** (11.5" W x 10") . . . . . N/C
- ICP **Comes standard with individual calf pads.** . . . . . N/C
- OPCP **One Piece Malleable Calf Panel** . . . . . \$250

1. 3" ground clearance required.  
2. Comes standard with moulded i-skin calf pads.


**MANUAL CENTER MOUNT FRONT RIGGING**

- MSFCMFP **Maxx Style Fixed Center Mount Foot Platform<sup>(1)(2)</sup>** . . . . . \$395  
With Quick Release Height Adjustment
  - FBIFP **Comes with a standard footplate with** (17" w x 13"d) . . . . . N/C
  - ICP **Individual calf pads** . . . . . \$200
- Select seat pan to footplate:
- SPTF1 **9" - 13"** . . . . . N/C
  - SPTF14 **13" - 19"** . . . . . N/C
- Set at:
- SA70 **Set at 70°** . . . . . N/C
  - SA90 **Set at 90°** . . . . . N/C
  - SA97 **Set at 97°** . . . . . N/C



1. Depth and height adjustable.  
2. 3" ground clearance required.


**FOOT PLATFORM OPTIONS FOR FIXED/LNX CMFP**



- MFFP **Medium Platform<sup>(1)</sup> (12"w x 10.75"d)** . . . . . \$180
- LFFP **Large Platform<sup>(1)</sup> (14.5"w x 13"d)** . . . . . \$250
- XLFFP **Extra Large Platform<sup>(1)</sup> (17"w x 13"d)** . . . . . \$250
- XXLFFP **Extra Extra Large Platform<sup>(1)</sup> (21"w x 13"d)** . \$250

1. May need to adjust leg position to avoid interference with front castors

**INDIVIDUAL FOOT PLATES FOR FIXED/LNX CMFP**



- IFSL **Small<sup>(1)</sup> - (5" x 7.5" deep) - Left** . . . . . \$165
- IFSR **Small<sup>(1)</sup> - (5" x 7.5" deep) - Right** . . . . . \$165
- IFML **Medium<sup>(1)</sup> - (5.5" x 9.5" deep) - Left** . . . . . \$165
- IFMR **Medium<sup>(1)</sup> - (5.5" x 9.5" deep) - Right** . . . . . \$165
- IFLL **Large<sup>(1)</sup> - (6" x 11.5" deep) - Left** . . . . . \$165
- IFLR **Large<sup>(1)</sup> - (6" x 11.5" deep) - Right** . . . . . \$165
- MFOL **Multi-Axis Optional Footplate Left (6"W x 10"D)** . \$225
- MFOR **Multi-Axis Optional Footplate Right (6"W x 10"D)** \$225

1. Foot plates may interfere with casters in certain configurations

**MAXX STYLE POWER ELEVATING LEGRESTS**

MSPPL **MAXX Swing Away Power Elevating<sup>(1)(2)</sup> ..... \$2,490**

Select Power Legrest Function:

INDA **Individual Movement ..... N/C**  
 COMA **Combined Movement ..... N/C**  
 SYNA **Synchronized/Recline Movement ..... N/C**

Select Receivers:

RS **Straight Receivers ..... N/C**  
 RF **Flared 1" Outward Receivers ..... N/C**

Select Seat Pan to Footplate:

SPTF5 **Length: 10"-14" ..... N/C**  
 SPTF2 **Length: 13"-17" ..... N/C**  
 SPTF7 **Length: 16"-20" ..... N/C**

FCPMS **Flat Calf Pad Option (pair) ..... \$275**

CCPMS **Curved Calf Pad Option (inside with 6") (pair) .. \$275**

1. Must choose footplate on page 8
2. Comes standard with moulded i-skin calf pads.

**SWING AWAY LEGREST**

HD70 **Heavy Duty 70° Swing Away Footrest(pair)<sup>(1)(2)</sup> ..... \$495**

Select Receivers:

RS **Straight Receivers ..... N/C**  
 RF **Flared 1" Outward Receivers ..... N/C**

Select Seat Pan to Footplate:

SPTF11 **Length: 6"-10" ..... N/C**  
 SPTF8 **Length: 11"-15" ..... N/C**  
 SPTF2 **Length: 13"-17" ..... N/C**  
 SPTF7 **Length: 16"-20" ..... N/C**

HDFCP **Flat Calf Pad Option (pair) ..... \$275**

HD70CCP **Curved Calf Pad Option (inside with 6") (pair) .. \$275**

DCS **Dynamic Coil Springs (Pair) for Heavy Duty 70° Swingaway Footrests ..... \$190**

1. Must choose footplate on page 8
2. Heavy duty durable construction (4.5lbs ea.).

**MAXX STYLE MANUAL ELEVATING LEGRESTS**

MMPPSA **MAXX Swing Away Manual Elevating<sup>(1)(2)</sup> ..... \$680**

Select Receivers:

RS **Straight Receivers ..... N/C**  
 RF **Flared 1" Outward Receivers ..... N/C**

Select Seat Pan to Footplate:

SPTF5 **Length: 10"-14" ..... N/C**  
 SPTF2 **Length: 13"-17" ..... N/C**  
 SPTF7 **Length: 16"-20" ..... N/C**


FCPMS **Flat Calf Pad Option (pair) ..... \$275**

CCPMS **Curved Calf Pad Option (inside with 6") (pair) .. \$275**

1. Must choose footplate on page 8
2. Comes standard with moulded i-skin calf pads




**FOOTPLATE FOR SWING AWAY/ELEVATING LEGS**



1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

**FOOTPLATE/CUT OUTS SIZE SELECTION**



1. Certain footplate sizes may not fit all seat widths depending on legrests selected

**ADJUSTABLE ANGLE FOOTPLATES**

<input type="checkbox"/> AFSL	Left Small Footplate (4.25"W x 8"D) .....	\$165
<input type="checkbox"/> AFSR	Right Small Footplate (4.25"W x 8"D) .....	\$165
<input type="checkbox"/> AFML	Left Medium Footplate (5.25"W x 8"D) .....	\$165
<input type="checkbox"/> AFMR	Right Medium Footplate (5.25"W x 8"D) .....	\$165
<input type="checkbox"/> AFLL	Left Large Footplate (6.25"W x 8"D) .....	\$165
<input type="checkbox"/> AFLR	Right Large Footplate (6.25"W x 8"D) .....	\$165
<input type="checkbox"/> AFXL	Left Extra Large Footplate (9"W x 11"D) .....	\$165
<input type="checkbox"/> AFXR	Right Extra Large Footplate (9"W x 11"D) .....	\$165

**ADJUSTABLE ANGLE FOOTPLATES/CUT OUTS**

<input type="checkbox"/> AFCSL	Left Small Footplate (7.25"W x 9"D) .....	\$165
<input type="checkbox"/> AFCSR	Right Small Footplate (7.25"W x 9"D) .....	\$165
<input type="checkbox"/> AFCML	Left Medium Footplate (7.5"W x 11"D) .....	\$165
<input type="checkbox"/> AFCMR	Right Medium Footplate (7.5"W x 11"D) .....	\$165
<input type="checkbox"/> AFCLL	Left Large Footplate (10"W x 14"D) .....	\$165
<input type="checkbox"/> AFCLR	Right Large Footplate (10"W x 14"D) .....	\$165

**MULTI-AXIS ADJUSTABLE ANGLE FOOTPLATES**

<input type="checkbox"/> MFSL	Left Small Footplate (4.25"W x 8"D) .....	\$225
<input type="checkbox"/> MFSR	Right Small Footplate (4.25"W x 8"D) .....	\$225
<input type="checkbox"/> MFML	Left Medium Footplate (5.25"W x 8"D) .....	\$225
<input type="checkbox"/> MFMR	Right Medium Footplate (5.25"W x 8"D) .....	\$225
<input type="checkbox"/> MFLL	Left Large Footplate (6.25"W x 8"D) .....	\$225
<input type="checkbox"/> MFLR	Right Large Footplate (6.25"W x 8"D) .....	\$225
<input type="checkbox"/> MFXL	Left Extra Large Footplate (9"W x 11"D) .....	\$225
<input type="checkbox"/> MFXR	Right Extra Large Footplate (9"W x 11"D) .....	\$225

**MULTI-AXIS ADJUSTABLE ANGLE FOOTPLATES/ CUT OUTS**

<input type="checkbox"/> MFCSL	Left Small Footplate (7.25"W x 9"D) .....	\$225
<input type="checkbox"/> MFCSR	Right Small Footplate (7.25"W x 9"D) .....	\$225
<input type="checkbox"/> MFCML	Left Medium Footplate (7.5"W x 11"D) .....	\$225
<input type="checkbox"/> MFCMR	Right Medium Footplate (7.5"W x 11"D) .....	\$225
<input type="checkbox"/> MFCLL	Left Large Footplate (10"W x 14"D) .....	\$225
<input type="checkbox"/> MFCLR	Right Large Footplate (10"W x 14"D) .....	\$225

**HEEL LOOPS, CUPS AND TOE STRAPS**

<input type="checkbox"/> HLL	Left Heel Loop (each) .....	\$20
<input type="checkbox"/> HLR	Right Heel Loop (each) .....	\$20
<input type="checkbox"/> HCL	Left Cup Loop (each) .....	\$20
<input type="checkbox"/> HCR	Right Cup Loop (each) .....	\$20
<input type="checkbox"/> TSL	Left Toe Strap (each) .....	\$20
<input type="checkbox"/> TSR	Right Toe Strap (each) .....	\$20

**CALF STRAPS AND CALF PANELS**

<input type="checkbox"/> CSL	Body Point Airmesh Calf Strap for Chairs 18"-22" wide with 3.5" H Pad <sup>(1)</sup> .....	\$200
<input type="checkbox"/> CPL	Body Point Airmesh Calf Panel for Chairs 18"-22" wide with 38.5" W x 9" H Pad <sup>(1)</sup> . . .	\$315

**THORACIC LATERALS AND HIP/KNEE SUPPORTS**

- Matrx Laterals please see page 12*
- Maxx Laterals and Hip/Knee supports please see page 13*

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**ELECTRONIC SWITCH MOUNTING**

- ESML Left Side (Armpad)
- ESMR Right Side (Armpad)
- ESMJ Joystick (Mounts to Joystick Bracket)
- ESMD Do not mount

**SINGLE-DUAL FUNCTION ELECTRONICS**

- SFCB Single/Dual Function Electronics<sup>(1)</sup> . . . . . STD

Select Switch:

- SPBO Single Push Button
- STO Single Toggle
- DPBO Dual Push Button
- SOS Omit Switch
- JOY Through Joystick
- SC Through Specialty Control

Optional:

- RDE Reduced Drive for Elevating Seat<sup>(2)</sup> . . . . . \$400
- SPBSA Additional Single Push Button . . . . . \$97
- STA Additional Single Toggle Switch . . . . . \$97
- DPBA Additional Dual Push Button . . . . . \$97

1. These electronics include drive lock-out and tilt limit.  
2. If reduced drive is not selected the system will be in drive lock-out when elevated.

**MULTI-FUNCTION ELECTRONICS**

- MFCB Multi-Function through Switch/Toggle<sup>(1)</sup> . . . . . \$1,950

Select Switch:

- QPBO Quad Push Button
- FWTO Four-Way Toggle
- EWRO 8-Way Rocker<sup>(2)</sup>

- MFCBD Multi-Function through Drive Control<sup>(1)</sup> . . . . . \$1,950

Select Control:

- JOY Through Joystick
- SC Through Specialty Control

Optional:

- RDY Include Reduced Drive for Elevating Seat<sup>(3)</sup>
- RDN Do not Induce Reduced Drive for Elevating Seat<sup>(3)</sup>
- QPBA Additional Quad Button . . . . . \$97
- FWTA Additional Four-Way Toggle Switch . . . . . \$97
- EWRA Additional 8-Way Rocker . . . . . \$150

1. These electronics include seating attendant control, drive lock-out, tilt/recline limit and actuator speed control.  
2. Dual function is not available with 8-Way Rocker.  
3. If reduced drive is not selected system will be in drive lock-out when elevated.

**ELECTRONICS OPTIONS**

- M616S M616 Converter<sup>(1)</sup> . . . . . \$125
- FTSQ Feather Touch Switch Single Function Quarter Size . . . . . \$295
- FTMD Feather Touch Switch Multi-Function Dual Dime Size . . . . . \$295
- FTMQ Feather Touch Switch Multi-Function Quad Dime Size . . . . . \$295
- ASL320R ASL Pulse Mechanical Switch Red . . . . . \$125
- ASL320G ASL Pulse Mechanical Switch Green . . . . . \$125
- ASL320BK ASL Pulse Mechanical Switch Black . . . . . \$125
- SESG Stealth Egg Switch - Single Function Green . . \$135
- SESB Stealth Egg Switch - Single Function Black . . \$135
- PSOS Power Shut Off Switch - Toggle<sup>(1)</sup> . . . . . \$175
- LJS Lazarus Jump Starter<sup>(1)</sup> . . . . . \$150
- APPS Accessory Port Power Supply<sup>(1)</sup> (3 amp capacity) . . . . . \$195
- TPSMH 12 Piece Switch Mounting Hardware Kit (additional kits) . . . . . \$39

1. Not available with Single-Dual Function Electronics.

**VENT TRAYS, O2 HOLDERS, ACCESSORIES**

- For O2 Holders and Additional Accessories  
Please see page 14

**INSTALLATION**

- MCI Motion Concepts Installation . . . . . Included

**PLEASE SPECIFY PAGES INCLUDED WITH ORDER**

- Pages 1-9 Excludes Matrix and Maxx Seating Accessories Thoraic Supports, Hip Supports, Vent Trays and Accessories
- Pages 1-14 Complete Order Form

**Total Retail Value \$** \_\_\_\_\_

**SUBMIT**

**PRINT FORM**      **SAVE**      **CLEAR FORM**



**CUSHION COVER REFERENCE GUIDE**

Meshtex	Startex	Infection Control			
Spacetex	Onyx				
Cushion	Meshtex	Startex	Infection	Spacetex	Onyx
Libra/Adjust	N/A	STD	\$50	N/C	N/A
PS/PSP	N/A	STD-Reversible	STD-Reversible	N/A	N/A
Vi	STD	\$103	\$103	N/A	\$103

**FABRIC OPTIONS FOR CUSHIONS**

- SFLF Spacetex Fabric Option for Libra Cushion ..... N/C
- ICFL Infection Control Fabric Upgrade for Libra Cushion ..... \$50
- OCUMA Onyx Fabric Upgrade for Matrix Vi Cushion ..... \$103
- SFVS Smooth side Startex Fabric Upgrade for Matrix Vi Cushion ..... \$103
- SFVF Fabric side Startex Fabric Upgrade for Matrix Vi Cushion ..... \$103

**MATRIX SEAT CUSHION OPTIONS**

Matrix Libra Adjust Seat Cushion - 20" wide ..... \$985

- LC2016-ADJ  LC2018-ADJ  LC2020-ADJ

Matrix Libra Adjust Seat Cushion 22" depth ..... \$1,265

- LC2022-ADJ

Matrix Libra Adjust Seat Cushion 21" - 24" wide ..... \$1,265

- LC2118-ADJ  LC2120-ADJ  LC2122-ADJ  
 LC2218-ADJ  LC2220-ADJ  LC2222-ADJ  
 LC2318-ADJ  LC2320-ADJ  LC2322-ADJ  
 LC2418-ADJ  LC2420-ADJ  LC2422-ADJ

Matrix Libra Seat Cushion - 20" wide ..... \$745

- LC2016  LC2018  LC2020

Matrix Libra Seat Cushion 22" depth ..... \$1,025

- LC2022

Matrix Libra Seat Cushion 21" - 24" wide ..... \$1,025

- LC2118  LC2120  LC2122  LC2218  LC2220  
 LC2320  LC2322  LC2418  LC2420  
 LC2422

Matrix PS Seat Cushion 20" wide ..... \$569

- PS2016  PS2018  PS2020

Matrix PS HD Seat Cushion 22" - 26" wide ..... \$975

- PS2218  PS2220  PS2222  PS2418  PS2420  
 PS2422  PS2618  PS2620  PS2622

Matrix PSVF Seat Cushion - 20" wide ..... \$589

- PSVF2016  PSVF2018  PSVF2020

Matrix PSP Seat Cushion - 20" wide ..... \$645

- PSP2016  PSP2018  PSP2020

Matrix Vi Seat Cushion 20" - 22" wide ..... \$681

- MA2016-VI  MA2018-VI  MA2020-VI  MA2218-VI  
 MA2220-VI

Matrix Vi HD Seat Cushion up to 26" wide ..... \$840

(For cushions 21" width and depth up to 26")

Specify size: \_\_\_\_\_



Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**BACK COVER REFERENCE GUIDE**

Back	Meshtex	Infection	Onyx
High	STD	N/A	\$103
Elite/ Elite HD/Elite Deep	STD	\$103	N/A
Elite TR/ Elite TR HD	STD	\$103	N/A
PB/PB HD/PB Deep/PB HD Deep	STD	\$103	N/A

**UPGRADABLE BACK OPTIONS**

Matrix PB Posture Back - 20" wide (5" Contour) . . . . . \$734

- MPB2012    MPB2016    MPB2020

Matrix PB Posture Back Deep - 20" wide (7" Contour) . . . . . \$892

- PBD2012    PBD2016    PBD2020

Matrix PB Heavy Duty Back 20" - 24" wide (5" Contour) . \$1,051

- HDB2016    HDB2020    HDB2216    HDB2220  
HDB2416    HDB2420

Matrix PB Heavy Duty Deep Back 20"-22" wide (5" Contour) \$1,209

- DHDB2016    DHDB2020    DHDB2216    DHDB2220  
DHDB2416    DHDB2420

**UPGRADABLE BACK OPTIONS**

High Back 20" - 24" wide . . . . . \$595

- HB2016    HB2018    HB2020    HB2116  
HB2118    HB2120    HB2216    HB2218  
HB2220    HB2316    HB2318    HB2320  
HB2416    HB2418    HB2420

Matrix Elite Back 20" wide (3" Contour) . . . . . \$766

- PBE2014    PBE2016    PBE2018    PBE2020

Matrix Elite HD Back 21" - 24" wide (3" Contour) . . . . \$1,095

- PBE2116-HD    PBE2118-HD    PBE2120-HD  
PBE2216-HD    PBE2218-HD    PBE2220-HD  
PBE2316-HD    PBE2318-HD    PBE2320-HD  
PBE2416-HD    PBE2418-HD    PBE2420-HD

Matrix Elite HD Back 26" wide (3" Contour) . . . . . \$1,666

- PBE2616-HD    PBE2620-HD

Matrix Elite Deep Back 20" wide (3" Contour) . . . . . \$918

- EDB2014    EDB2016    EDB2018    EDB2020

Matrix Elite Deep HD Back 21" - 24" wide (6" Contour) . \$1,095

- EDB2116-HD    EDB2118-HD    EDB2120-HD  
EDB2216-HD    EDB2218-HD    EDB2220-HD  
EDB2316-HD    EDB2318-HD    EDB2320-HD  
EDB2416-HD    EDB2418-HD    EDB2420-HD

Matrix Elite TR Back 20" wide (3" Contour) . . . . . \$866

- PBE2016-TR    PBE2018-TR    PBE2020-TR

Matrix Elite TR HD Back 21" - 24" wide (3" Contour) . . \$1,151

- PBE2120-TRHD    PBE2220-TRHD    PBE2320-TRHD  
PBE2420-TRHD

**OTHER UPGRADABLE BACK OPTIONS**

- EBSS    Super Soft Foam for Elite Back Upgrade . . . . \$103  
OCUHB    Onyx Fabric Upgrade for High Backs. . . . . \$103  
WSCU    Infection Control Cover for Elite/PB Backs . . . \$103  
*purchased with back*  
WSC    Infection Control Cover for Elite/PB Backs . . . \$150  
*purchased with without back*  
CHIEB    Chest Harness Interface Kit for Elite Back<sup>(1)(2)</sup> . \$103  
CHIEB    Chest Harness Interface Kit for PB Back<sup>(1)(2)</sup> . . . \$52

1. Does not include the chest strap.  
2. This Interface Kit is only Required for 4 point harness



Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**MATRIX ELAN HEADREST PAD<sup>(1)(2)</sup>**

<input type="checkbox"/> ESP6	Matrx Elan Standard Pad 6" .....	\$310
<input type="checkbox"/> ESP10	Matrx Elan Standard Pad 10" .....	\$310
<input type="checkbox"/> ESP14	Matrx Elan Standard Pad 14" .....	\$310
<input type="checkbox"/> EOP9	Matrx Elan Occipital Pad 9" .....	\$341
<input type="checkbox"/> EOP12	Matrx Elan Occipital Pad 12" .....	\$341
<input type="checkbox"/> E4PTSM	Matrx Elan 4-Point Pad 9"W x 7"H .....	\$448
<input type="checkbox"/> E4POINT	Matrx Elan 4-Point Pad 11"W x 10"H .....	\$448
<input type="checkbox"/> ESP6-IC	Matrx Elan Standard Pad 6" with Infection Control .....	\$335
<input type="checkbox"/> ESP10-IC	Matrx Elan Standard Pad 10" with Infection Control .....	\$335
<input type="checkbox"/> ESP14-IC	Matrx Elan Standard Pad 14" with Infection InfectionControl .....	\$335
<input type="checkbox"/> EOP9-IC	Matrx Elan Occipital Pad 9" with Infection Infection Control .....	\$366
<input type="checkbox"/> EOP12-IC	Matrx Elan Occipital Pad 12"with Infection Control .....	\$366
<input type="checkbox"/> E4PTSM-IC	Matrx Elan 4-Point Pad 9"W x 7"H with Infection Control .....	\$473
<input type="checkbox"/> E4POINT-IC	Matrx Elan 4-Point Pad 11"W x 10"H with Infection Control .....	\$473

1. Not compatible with sling back (ASBA-2)  
2. If mounting on PB Elite must order removable hardware

**HEADREST MOUNTING HARDWARE OPTIONS**

<input type="checkbox"/> MEHWRE	Elan Multi-Axis Removable HDW.....	\$285
---------------------------------	------------------------------------	-------

1. Only compatible w/Standard Rehab, Contour & High Back Models. All other head rests REQUIRE an upgrade to Removable Hardware.

**MATRIX FIXED LATERALS - FOR ELITE AND ELITE TR BACKS**

<input type="checkbox"/> PBSSL	Matrx Small Lateral Fixed-Left(3x5.5") .....	\$113
<input type="checkbox"/> PBSSR	Matrx Small Lateral Fixed-Right(3x5.5"). . . .	\$113
<input type="checkbox"/> PBLSL	Matrx Standard Lateral Fixed-Left(7.5x5") .	\$116
<input type="checkbox"/> PBLSR	Matrx Standard Lateral Fixed-Right (7.5x5")	\$116
<input type="checkbox"/> PBOLL	Matrx Offset Lateral Fixed Left . . . . .	\$125
<input type="checkbox"/> PBOLR	Matrx Offset Lateral Fixed Right . . . . .	\$125

**MATRIX SWING AWAY LATERALS - FORELITE AND ELITE TR BACKS**

Matrx Lateral Trunk Support Left <sup>(1)</sup>		
<input type="checkbox"/> EPADXL	Size XS (3.25x4.25") .....	\$115
<input type="checkbox"/> EPADSL	Size S (3.5x5.25") .....	\$115
<input type="checkbox"/> EPADM L	Size M (3.75x5.75") .....	\$115
<input type="checkbox"/> EPADLL	Size L (4.25x6.5") .....	\$115
Matrx Lateral Trunk Support Right <sup>(1)</sup>		
<input type="checkbox"/> EPADXR	Size XS (3.25x4.25") .....	\$115
<input type="checkbox"/> EPADSR	Size S (3.5x5.25") .....	\$115
<input type="checkbox"/> EPADMR	Size M (3.75x5.75") .....	\$115
<input type="checkbox"/> EPADLR	Size L (4.25x6.5") .....	\$115

1. Must choose swing away hardware below

**MATRIX SWING AWAY HARDWARE - FORELITE AND ELITE TR BACKS**

<input type="checkbox"/> SAEHBL	Matrx Elite Swing Away Hardware Left . . . . .	\$250
<input type="checkbox"/> SAHEBR	Matrx Elite Swing Away Hardware Right . . . . .	\$250
<input type="checkbox"/> OSAHL	Matrx Offset Elite Swing Away Hardware Left . . . . .	\$250
<input type="checkbox"/> OSAHR	Matrx Offset Elite Swing Away HardwareRight . . . . .	\$250

**MAXX LATERALS TRUNK SUPPORTS<sup>(1)</sup>**

**Left Fixed Mount Lateral Trunk Support<sup>(1)</sup>**

Select Left Trunk Support Size

- LATXL Size XS (3.25" x 4.25") . . . . . \$100
- LATSL Size S (3.5" x 5.25") . . . . . \$100
- LATML Size M (3.75" x 5.75") . . . . . \$100
- LATLL Size L (4.25" x 6.5") . . . . . \$100

**Right Fixed Mount Lateral Trunk Support<sup>(1)</sup>**

Select Right Trunk Support Size

- LATXR Size XS (3.25" x 4.25") . . . . . \$100
- LATSR Size S (3.5" x 5.25") . . . . . \$100
- LATMR Size M (3.75" x 5.75") . . . . . \$100
- LATLR Size L (4.25" x 6.5") . . . . . \$100

1. These laterals are not available with Matrix PB, PB Elite Deep and PB Deep Backs.

**MAXX STYLE SWING-AWAY TRUNK SUPPORT HARDWARE<sup>(1)</sup>**

- SALATL Left MAXX Style Swing Away HDW . . . . . \$195

Select Left Trunk Support Size

- LATXL Size XS (3.25" x 4.25") . . . . . \$100
- LATSL Size S (3.5" x 5.25") . . . . . \$100
- LATML Size M (3.75" x 5.75") . . . . . \$100
- LATLL Size L (4.25" x 6.5") . . . . . \$100

- SALATR Right MAXX Style Swing Away HDW<sup>(1)</sup> . . . \$195

Select Right Trunk Support Size

- LATXR Size XS (3.25" x 4.25") . . . . . \$100
- LATSR Size S (3.5" x 5.25") . . . . . \$100
- LATMR Size M (3.75" x 5.75") . . . . . \$100
- LATLR Size L (4.25" x 6.5") . . . . . \$100

- ATLL Additional Telescoping Link - Left<sup>(2)</sup> . . . \$75

- ATLR Additional Telescoping Link - Right<sup>(2)</sup> . . \$75

1. These laterals are not available with Matrix PB, PB Elite Deep and PB Deep Backs  
2. Up to 2.75" medial offset

**FIXED HIP/KNEE SUPPORTS**

**Lateral Hip Support with Fixed Mounting Hardware Left<sup>(1)</sup>**

Select Pad Size and mounting position at Hip:

- HSPXLH Size XS (5.5" W x 4" H) mounted at Hip . \$120
- HSPSMLH Size S/M (7"W x 4"H) mounted at Hip . . \$120
- HSPSLH Size S (7"W x 6"H) mounted at Hip . . . . \$120
- HSPMLH Size M (9" W x 4"H) mounted at Hip . . . \$120
- HSPLLH Size L (13" W x 4"H) mounted at Hip . . . \$120

Select Pad Size and mounting position at Knee:

- HSPXLK Size XS (5.5"W x 4" H) mounted at Knee . . . \$120
- HSPSMLK Size S/M (7"W x 4" H) mounted at Knee . . . \$120
- HSPSLK Size S (7"W x 6" H) mounted at Knee . . . . . \$120
- HSPMLK Size M (9"W x 4"H) mounted at Knee . . . . \$120
- HSPLLK Size L (13"W x 4"H) mounted at Knee . . . . . \$120

**FIXED HIP/KNEE SUPPORTS- CONTINUED**

**Lateral Hip Support with Fixed Mounting Hardware Right<sup>(1)</sup>**

Select Pad Size and mounting position:

- HSPXRH Size XS (5.5"W x 4"H) mounted at Hip . . . . \$120
- HSPSMRH Size S/M (7"W x 4"H) mounted at Hip . . . . \$120
- HSPSRH Size S (7"W x 6"H) mounted at Hip . . . . . \$120
- HSPMRH Size M (9"W x 4"H) mounted at Hip . . . . . \$120
- HSPLRH Size L (13"W x 4"H) mounted at Hip . . . . . \$120

Select Pad Size and mounting position at Knee:

- HSPXRK Size XS(5.5"W x 4"H) mounted at Knee . . . \$120
- HSPSMRK Size S/M (7"W x 4"H) mounted at Knee . . . \$120
- HSPSRK Size S (7"W x 6"H) mounted at Knee . . . . . \$120
- HSPMRK Size M (9"W x 4"H) mounted at Knee . . . . . \$120
- HSPLRK Size L (13"W x 4"H) mounted at Knee . . . . . \$120

1. Must choose swing away hardware below if required.

**REMOVABLE HIP/KNEE SUPPORT HARDWARE**

**Lift Off Removable Hardware for Hip Support**

- LOHL Left . . . . . \$195
- LOHR Right . . . . . \$195

**MAXX Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports**

Select mounting position:

- MHSHLH Left mounted at Hip . . . . . \$215
- MHSHLK Left mounted at Knee . . . . . \$215
- MHSHRH Right mounted at Hip . . . . . \$215
- MHSHRK Right mounted at Knee . . . . . \$215

**Swing Away Removable Hardware for Hip Supports**

Select mounting position:

- SHSHLH Left mounted at Hip . . . . . \$215
- SHSHLK Left mounted at Knee . . . . . \$215
- SHSHRH Right mounted at Hip . . . . . \$215
- SHSHRK Right mounted at Knee . . . . . \$215

**Extended Maxx Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports<sup>(1)</sup>**

Select mounting position:

- EMQHLH Left mounted at Hip . . . . . \$215
- EMQHLLK Left mounted at Knee . . . . . \$215
- EMQHRH Right mounted at Hip . . . . . \$215
- EMQHRK Right mounted at Knee . . . . . \$215

1. Recommended when finished seat cushion thickness will be greater than 3.5" high.

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**ELBOW BLOCKS**

**Left Elbow Block<sup>(1)</sup>**

Select Pad Size and Both Extensions:

- EBXL Size XS (3.25" x 4.25") . . . . . \$195
- EBSL Size S (3.5" x 5.25") . . . . . \$195
- EBML Size M (3.75" x 5.75") . . . . . \$195
- EBLL Size L (4.25" x 6.5") . . . . . \$195
- V3 Left Verticle Extension3"
- V5 Left Verticle Extension 5"
- H5 Left Horizontal Extension 5"
- H7 Left Horizontal Extension 7"

**Right Elbow Block<sup>(1)</sup>**

Select Pad Size and Both Extensions:

- EBXL Size XS (3.25" x 4.25") . . . . . \$195
- EBSL Size S (3.5" x 5.25") . . . . . \$195
- EBML Size M (3.75" x 5.75") . . . . . \$195
- EBLL Size L (4.25" x 6.5") . . . . . \$195
- V3 Right Verticle Extension3"
- V5 Right Verticle Extension 5"
- H5 Right Horizontal Extension 5"
- H7 Right Horizontal Extension 7"

1. Prevents arm from slipping.

**ADDITIONAL ARM ACCESSORIES**

- UCH Universal Cup Holder<sup>(1)</sup> . . . . . \$35
- MACESL Multi Axis Upper Extremity Support<sup>(2)</sup>(each)  
Left Side . . . . . \$225
- MACESR Multi Axis Upper Extremity Support<sup>(2)</sup> (each)  
Right Side . . . . . \$225
- OAS2L 1" Offset Armest Spacer Left(each) . . . . \$125
- OAS2R 1" Offset Armest Spacer Right(each) . . . \$125
- OAS1L 1/2" Offset Armest Spacer Left(each) . . . \$125
- OAS1R 1/2" Offset Armest Spacer Right(each) . . \$125
- MPAPA Multi Purpose Arm Pad  
Adapter Block<sup>(3)</sup> (each) Quantity ( ) . . . . \$35

1. Not available on the joystick side.  
2. Allows rotation of the armpad, only compatible with Gel Arm pad, Flat Arm pad or Ergonomic Arm Troughs.  
3. Not available with Standard Arm Pads or Waterfall Arm Pads.

**O2 HOLDERS<sup>(1)</sup>**

- NAO2H Non-Articulating O<sup>2</sup> Gas Holder . . . . . \$395
- AOHG Articulating O<sup>2</sup> Gas Holder  
Cylinder Style<sup>(2)(3)</sup> . . . . . \$1,250
- AOHL Articulating Liquid O<sup>2</sup> Holder  
(6.5" x 4.5")<sup>(2)(3)</sup> . . . . . \$1,250

1. The user weight limit is 325 lbs. The maximum weight capacity of the O2 holder is 20 lbs.  
2. Required for use on Recline Systems  
3. Articulating O2 Holders are not standard configuration on Elevating or Precline Systems.

**ADDITIONAL OPTIONS & ACCESSORIES**

- THS Transfer Handles 4" - 6.5" . . . . . \$350
- THL Transfer Handles 9" - 11.5" . . . . . \$350
- FFLTL Flip and Fold Away 1/2 Lap Tray  
Left Mounting(14"W x 10"D)<sup>(1)(2)</sup> . . . . . \$795
- FFLTR Flip and Fold Away 1/2 Lap Tray  
Right Mounting(14"W x 10"D)<sup>(1)(2)</sup> . . . . . \$795

1. Not available on joystick side and only available with a flat arm pad or modular arm pad.  
2. The maximum weight capacity of the tray is 10 lbs.

**SUBMIT**

PRINT FORM

SAVE

CLEAR FORM