INVACARE® IV WHEELCHAIR







The **Invacare® Tracer® IV Wheelchair** features superior durability, rollability and streamlined looks. It is designed for the active individual who needs the strength of a reinforced heavy-duty frame. Its many options include components interchangeable within the Tracer series. The Tracer IV wheelchair is long-lasting and low-maintenance.

Features

Dual axle position allows variation of seat-to-floor heights

Durable, low-maintenance, triple chrome-plated carbon steel frame is long-lasting

Reinforced frame includes gussetted sides, caster journals and dual crossbraces

Corded urethane rear tires (350 lb.) and pneumatic flat free tires

(450 lb.) mounted on "no flex" wheels offer superior performance

Low-maintenance, dual-sealed precision bearings

Heavy-duty inner liners keep seat and back from stretching

Upholstery is dual-embossed

HCPCS Code: K0007













- 1. Dual axle positions
- 2. Corded urethane tires
- 3. Dual-embossed upholstery
- 4. Bariatric sizes available
- 5. Seating available
- 6. Builder options



















Model# T420, T422, T424, T4X22, T4X24

Overall Width 29'

Overall Height 36"

Overall Depth w/o riggings 30"

Folded Width 13"

Seat Dimensions 18"-24"x18", 20"

Seat-to-Floor Height 17.5", 19.5"

Back Height 16"



8"







Weight







Limited

T420, T422, T424, T4X22, T4X24

Wheel Caster Height

Height 24"

Weight w/o Rigging 42 lb.

Capacity

350 lb.,450 lb.

Shipping Weight 53 lb.

Carton Dimensions 32"x13"x 37"

Warranty 5 Years



Invacare Corporation www.invacare.com

570 Matheson Blvd. E.,

Mississauga, Ontario L4Z 4G4 Canada

(800) 668 - 5324

All trademarks are

unless otherwise

owned by or licensed

noted. Specifications are subject to change without notification. Form No. 11-228 Rev. 04/14 131070

to Invacare Corporation

USA One Invacare Way Elyria, Ohio

44035-4190 (800) 333-6900

Canada

Unit 8

©2014 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.