



NAME _____ PHONE NO _____
ACCOUNT # _____ P.O. _____ DATE _____
ADDRESS _____

FRAME & SEAT

CHAIR TYPE

☐ TRSX5 Tracer SX5 \$950

1 FRAME TYPE

☐ ADULT Adult 19 1/2" No Charge
☐ HEMI Hemi 17 1/2" No Charge

FRAME WIDTH & DEPTH

☐ WD46 14" Wide X 16" Deep \$145
☐ WD48 14" Wide X 18" Deep \$145
☐ WD66 16" Wide X 16" Deep No Charge
☐ WD68 16" Wide X 18" Deep \$145
☐ WD86 18" Wide X 16" Deep No Charge
☐ WD88 18" Wide X 18" Deep \$145
2 ☐ WD06 20" Wide X 16" Deep No Charge
2 ☐ WD26 22" Wide X 16" Deep \$145

ARM & BACK

3 ARM TYPE

☐ 18FB Full Length Fixed Height Flip-Back 10" \$86
☐ 19FB Full Length Adjustable Height Flip-Back 10" - 14" \$199
☐ 28FB Desk Length Fixed Height Flip-Back 10" \$86
☐ 29FB Desk Length Adjustable Height Flip-Back 10" - 14" \$199

BACK TYPE

☐ BH16 Back Height 16" Standard
☐ BH18 Back Height 18" No Charge

COLOURS

UPHOLSTERY COLOURS

☐ U240 Black Nylon Upholstery No Charge
☐ U550 Midnight Blue Vinyl Upholstery No Charge
☐ U67 Black Vinyl Standard
☐ U86 Ruby Red Vinyl No Charge

CASTERS

CASTER OPTIONS

☐ AT902 6" Composite Urethane 3 Spoke \$120
☐ 1228 8" X 1 3/4" Semi Pneumatic Casters Light Grey \$150
☐ AT903 8" X 1 1/4" Composite Urethane Standard
☐ 1252 8" X 1 1/4" Pneumatic Casters \$116
☐ 4252 8" x 1 1/4" Pneumatic Casters w/FF Inserts \$228
☐ 1255 8" X 1" Solid Rubber Dark Grey No Charge

WHEELS

WHEEL OPTIONS

☐ U1763BK 20" Composite Urethane Tires No Charge
☐ U1764BK 22" Composite Urethane Tires No Charge
☐ 2244BK 24" Composite Pneumatic Tires \$231
☐ 4244BK 24" Composite Pneumatic w/FF Inserts \$353
☐ U2244BK 24" Composite Urethane Treaded Tires No Charge
☐ U2222C 24" Composite Urethane Corded Tires Standard

FOOTNOTES

- 1 Based on 24" wheels and 8" casters, see seat-to-floor matrix.
2 20" and 22" widths have 300 lb wt limit, standard with 1228 casters, 24" rear wheels and CHDRM handrims
3 Arms are flipback non-removable.

4 WHEEL OPTIONS (ONE ARM DRIVE)

<input type="checkbox"/> 99457R	One Arm Drive w/2244BK/Flat Free RH	\$782
<input type="checkbox"/> 99458L	One Arm Drive w/2244BK/Flat Free LH	\$782
<input type="checkbox"/> 9456R	One Arm Drive w/U2222C RH	\$672
<input type="checkbox"/> 9456L	One Arm Drive w/U2222C LH	\$672

HANDRIM OPTIONS

5 <input type="checkbox"/> AHDRM	Aluminum Handrims	No Charge
6 <input type="checkbox"/> COM	Composite Handrims	Standard

WHEEL LOCKS

<input type="checkbox"/> 9690B	Pull to Lock Wheel Locks (Standard on Hemi, Bolt Through)	Standard
<input type="checkbox"/> 9691B	Push to Lock Wheel Locks (Standard on Adult, Bolt Through)	Standard

POSITIONING STRAPS
SEAT POSITIONING STRAPS

<input type="checkbox"/> 1515	Auto Style Buckle Seat Positioning Strap	\$84
<input type="checkbox"/> 1321BK	Hook & Loop Seat Positioning Strap	\$66

ACCESSORIES
BACK ACCESSORIES

<input type="checkbox"/> L8601	Omit Back Upholstery	Less net (\$20)
<input type="checkbox"/> 1160	Carry Pocket	\$83
7 <input type="checkbox"/> STC	Stencil Back Upholstery	\$59

ARM ACCESSORIES

<input type="checkbox"/> 1316	Black Plastic Arm Pads	No Charge
-------------------------------	------------------------	-----------

SEAT ACCESSORIES

<input type="checkbox"/> L8651	Omit Seat Upholstery	Less net (\$15)
9 <input type="checkbox"/> EXT2	Seat Extension Kit 2"	\$145

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/> 1360	Adjustable Rear Anti-Tipper For 8" Casters	\$109
8 <input type="checkbox"/> 1360A	Adjustable Rear Anti-Tipper For 6" Casters	\$123
10 <input type="checkbox"/> 9357L	Amputee Attachment	\$131
<input type="checkbox"/> 1329BK	Toggle Wheel Lock Extension (8")	No Charge
<input type="checkbox"/> 1335	Non-Folding Device	\$63
8,10 <input type="checkbox"/> 1577QR	Quick Release Axles Dual Height	\$142
10 <input type="checkbox"/> 2218	IV Rod & Bracket - Telescopic	\$185
11 <input type="checkbox"/> MUPH	Matching Upholstery	No Charge

FRONT RIGGING

12 <input type="checkbox"/> T93HC	Hemi Footrests	\$159
12 <input type="checkbox"/> T93HA	Hemi Footrests Aluminum Footplates	\$159
2,13 <input type="checkbox"/> T94HC	Hemi Elevating Legrests	\$302
2,13 <input type="checkbox"/> T94HA	Hemi Elevating Legrests Aluminum Footplates	\$309
2,13 <input type="checkbox"/> AHL4A	Hemi Smartleg Articulating Legrests Aluminum Footplates	\$399
2,13 <input type="checkbox"/> AHL4C	Hemi Smartleg Articulating Legrests Composite Footplates	\$392
3,14 <input type="checkbox"/> AHS4A	Hemi Smartleg Articulating Legrests 14" Wide Chairs	\$399
14 <input type="checkbox"/> T93HS	Hemi Swingaway Footrests Aluminum Footplates 14" wide	\$158
3,14 <input type="checkbox"/> T94HS	Hemi Swingaway Elevating Legrests Aluminum Footplates 14" wide	\$302
12 <input type="checkbox"/> T93HE	Hemi Footrests, w/o Heel Loop, Composite Footplates	\$119
12 <input type="checkbox"/> T94HE	Hemi Legrest, Composite Footplate, Black Plastic Calf Pad	\$239
12 <input type="checkbox"/> T93HAB	Hemi Bolt-Thru Footrest, Aluminum Footplates	\$158
12 <input type="checkbox"/> T93HCB	Hemi Bolt-Thru Footrest, Composite Footplates	\$158
14 <input type="checkbox"/> T93HSB	Hemi Bolt-Thru Footrest, Aluminum Footplates 14" wide	\$158
12 <input type="checkbox"/> T94HAB	Hemi Bolt-Thru Legrest, Aluminum Footplates	\$302
12 <input type="checkbox"/> T94HCB	Hemi Bolt-Thru Legrest, Composite Footplates	\$302
14 <input type="checkbox"/> T94HSB	Hemi Bolt-Thru Legrest, Aluminum Footplate 14" wide	\$302

FOOTNOTES

- 4 Not available with amputee attachment.
5 Not available with U2222C rear wheels.
6 Available with 24" rear wheels only
7 3/4" block letters
8 Not available on 20" and 22" widths.
9 Not available on 18" deep chair.
10 Not available with one-arm drive
11 Applies to arm pads and calf pads.
12 Not available on 14" width.
13 All calf pads are black.
14 Available only on 14" width.
15 Standard on swingaway footrests.

FOOTPLATE TYPE

<input type="checkbox"/> 1350	Extra Large Footplates	\$106
<input type="checkbox"/> AT5543	Adjustable Angle Flip-Up Footplates	\$132

FOOTREST/LEGREST ACCESSORIES

<input type="checkbox"/> 1820	Non-Removable Device	\$41
15 <input type="checkbox"/> 1346	Heel Loops	\$63
<input type="checkbox"/> L8515	Omit Calf Pads	Less net (\$10)

STOCK OPTIONS

10 <input type="checkbox"/> 1497	Oxygen Holder Package Assembly	\$198	
10 <input type="checkbox"/> 1496	IV/Oxygen Holder Package Assembly	\$198	FOOTNOTES
17 <input type="checkbox"/> 1514	Overhead Anti-Theft Device	\$136	16 Not available with 14" and 22" widths.



The Tracer SX5 Wheelchair is also available in several stock configurations. Please refer to order form 06-479C.

Note: All specifications and dimensions are approximate.

Specifications	
Seat Widths	14", 16", 18", 20" or 22"
Overall Width	Seat width plus 8 1/2"
Seat Depth	16" or 18"
Overall Depth	30 1/4" (18"x16" without front riggings)
Back Heights	16" or 18"
Seat-to-Floor	17 1/2" or 19 1/2"
Weight Limit	250 lbs. (300 lbs. on 20" or 22" widths)
Chair Weight (approx.)	33 3/4 lbs. (18"x16" without front riggings)
Shipping Weight (approx.)	44 1/4 lbs. (18"x16" without front riggings)
Please Note: Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.	

6-INCH FRONT CASTERS			
SEAT-TO-FLOOR HEIGHT (IN INCHES)	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
15 1/2	TOP	20 inch	TOP
16 1/2	MIDDLE	22 inch	TOP
17 1/2	BOTTOM	24 inch	TOP
17 1/2	BOTTOM	20 inch	BOTTOM

8-INCH FRONT CASTERS			
SEAT-TO-FLOOR HEIGHT (IN INCHES)	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
17 1/2	TOP	20 inch	BOTTOM
17 1/2	TOP	24 inch	TOP
18 1/2	MIDDLE	22 inch	BOTTOM
19 1/2	BOTTOM	24 inch	BOTTOM



Invacare Canada reserves the right to change prices and specifications without notice.

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324 (Canada) Fax: (905) 890-5244 (800) 668-5478 (Canada)
 nvacare®, and Yes, you can., are trademarks of Invacare Corporation.
 WWW.INVACARE.CA