INVACARE® TRACER® SX5 SERIES WHEELCHAIRS





The **Invacare® Tracer® SX5 Wheelchair's** lightweight frame, weighing less than 36 lb., makes this chair perfect for rental and long- or short-term use. This wheelchair offers multiple seat widths, arm styles, sectional back and seat-to-floor heights by utilizing the dual axle position frame base. The 20" and 22" widths include a semi-pneumatic caster to accommodate a weight capacity of 300 lb. Combining the design and technology of the Invacare® 9000 series wheelchairs make it possible to interchange components with the Tracer SX5 wheelchair.

The **Invacare® Tracer® SX5 Recliner Wheelchair** offers the durability of a lightweight manual wheelchair frame with the added versatility and comfort of a full reclining back. The Tracer SX5 Recliner wheelchair offers dynamic recline range from 90° to 180°. This chair is available with multiple seat widths, arm styles, front riggings, a full range of accessories and a 300 lb. weight capacity on 20" and 22" widths.

Invacare® Tracer® SX5 Wheelchair Features

Durable, low-maintenance, carbon steel frame is long-lasting Dual axle positions allow variation of seat-to-floor heights

14-gauge crossbraces add strength and durability

Urethane rear tires mounted on "no flex" wheels offer superior performance

Sectional back available in 16" and 18" back heights Low-maintenance, dual-sealed precision bearings

Heavy-duty inner liners keep seat and back from stretching Meets RESNA standards for manual wheelchairs in product

performance

Invacare® Tracer® SX5 Recliner Wheelchair Features

Dynamic range of recline from 90° to 180°

Durable, low-maintenance, carbon steel frame

Urethane rear tires, mounted on "no flex" wheels, offer superior performance

Low-maintenance, dual-sealed precision bearings

Heavy-duty inner liners keep seat and back from stretching Upholstery is dual-embossed

Meets RESNA standards for manual wheelchairs in product performance

HCPCS Code: K0003/K0006 & E1226 (recliner)











- 1. Dual axle positions
- 2 Urethane rear tires

TRSX50, TRSX52

TRSX56, TRSX58 TRSX5RC6, TRSX5RC8

- 3. Headrest
- 4. Seating options available
- 5. Builder options available



Invacare Corporation

www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668 - 5324

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 11-225 Rev. 04/14 131066



















Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Dimensions	Seat-to- Floor Height	Back Height
TRSX50, TRSX52 TRSX56, TRSX58	22.5"-30.5"	N/A	34"	N/A	14"-22"x16"	17.5", 19.5"	16", 18"
TRSX5RC6, TRSX5RC8	N/A	N/A	30.25"	N/A	14"-22"×16"	17.5", 19.5"	34"





Wheel





Weight



Shipping



Carton

Dimensions

32"x13"x 37"

Limited Warranty		
5 Years	-	
5 Years		

Caster	
Height	

Height 6", 8" 20", 22", 24" 8" 20", 22", 24"

Weight w/o Rigging 34 lb

53 lb

Capacity 250 lb./300 lb.

Weight 43 lh 250 lb./300 lb.

62 lb

36"x13"x 44" 5 Years

Yes, you can.