



NAME _____ PHONE NO _____

ACCOUNT # _____ P.O. _____ DATE _____

ADDRESS _____

STOCK WHEELCHAIR MODELS - Silver Vein Finish, Midnight Blue Vinyl Upholstery

Please select a front rigging option.

- TRSX56FBP 16" X 16" Frame with Flip-Back Fixed Height, Desk Length Arms \$819
TRXS5/WD66/ADULT/28FB/BH16/1255/U2222C/U240
- TRSX58FBP 18" X 16" Frame with Flip-Back Fixed Height, Desk Length Arms \$819
TRXS5/WD86/ADULT/28FB/BH16/1255/U2222C/U240
- TRSX58FBFP 18" X 16" Frame with Flip-Back Fixed Height, Full Length Arms \$819
TRXS5/WD86/ADULT/18FB/BH16/1255/U2222C/U240
- TRSX50FBP 20" X 16" Frame with Flip-Back Fixed Height, Desk Length Arms \$819
TRXS5/WD06/ADULT/28FB/BH16/1228/U2222C//U240
- TRSX50FBFP 20" X 16" Frame with Flip-Back Fixed Height, Full Length Arms \$819
TRXS5/WD06/18FB/BH16/ADULT/U2222C/1228/U240
- TRSX52FBP 22" X 16" Frame with Flip-Back Fixed Height, Desk Length Arms \$969
TRXS5/WD26/ADULT/28FB/BH16/1228/U2222C/U240

STOCK FRONT RIGGING OPTIONS - Wet Black Finish

- T93HAP Hemi Footrests, Aluminum Footplates with heel loops \$159
- T93HCP Hemi Footrests, Composite Footplates with heel loops \$159
- T94HAP Hemi Elevating Legrests; Aluminum Footplates & Padded Calf Pads \$309
- T94HCP Hemi Elevating Legrests; Composite Footplates & Padded Calf Pads \$309
- AHL4A Hemi Smartleg Articulating Legrests; Aluminum Footplates & Padded Calf Pads \$399
- T93HEP Economy Hemi Footrests; Composite Footplates without Heel Loops \$119
- T94HEP Economy Hemi Elevating Legrests; Composite Footplates & Calf Pads \$239

OTHER FRONT RIGGINGS - Silver Vein Finish (2-3 week lead time)

- T93HA_PTO Hemi Footrests, Aluminum Footplates with heel loops \$159
- T93HC_PTO Hemi Footrests, Composite Footplates with heel loops \$159
- T94HA_PTO Hemi Elevating Legrests; Aluminum Footplates & Padded Calf Pads \$309
- T94HC_PTO Hemi Elevating Legrests; Composite Footplates & Padded Calf Pads \$309
- T93HE_PTO Economy Hemi Footrests; Composite Footplates without Heel Loops \$119
- T94HE_PTO Economy Hemi Elevating Legrests; Composite Footplates & Calf Pads \$239

STOCK OPTIONS & ACCESSORIES

ANTI-TIPPERS

- 1360 Adjustable Rear Anti-Tippers \$109

MISCELLANEOUS

- 1329BKST 8" Wheel Lock Extensions \$63
- 1335ST Non-Folding Device \$111
- 1514 Overhead Anti-Theft Device \$136
- 1515ST Push Button Style Seat Positioning Strap \$127
- 1496 IV/Oxygen Cylinder Holder Package Assembly \$198
- 1497 Oxygen Cylinder Holder Package Assembly \$198
- 1820ST Non-Removable Device for Front Riggings \$52

Footnotes:
 1 Not available on 20" and 22" widths
 2 Not available on 22" widths

All Stock Tracer SX5 Wheelchairs Include: Standard	
Dual Axle Positions (17.5" and 19.5" Seat-to-Floor)	
16" Back Height	
8" x 1" Solid Rubber Casters (8 x 1.75" on 20" and 22" widths)	
24" Composite Wheels with Urethane Tires	
Composite Handrims (Chrome Plated on 20" and 22" widths)	
Black Nylon Upholstery (U240)	
Padded Arm pads	
Wheel Locks	

Specifications	
Seat Widths	16", 18", 20", or 22"
Overall Width	Seat Width plus 8.5"
Seat Depth	16"
Overall Depth	31.875" (18 x 16" without front riggings)
Back Heights	16"
Seat-to-Floor	17.5" or 19.5"
Weight Limit	250 lbs (300 lbs on 20" and 22" widths)
Chair Weight	34.75 lbs (18 x 16" without front riggings)
Shipping Weight	44.25 lbs

The Tracer SX5 Wheelchair can also be custom ordered. Please refer to order form 05-362C

Note: All specifications and dimensions are approximate.

Invacare Canada reserves the right to change prices and specifications without notice.
HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)
Invacare®, and Yes, you can., are trademarks of Invacare Corporation.
WWW.INVACARE.CA