



Invacare® Tracer® SX5 Recliner
Suggested Canadian Price List
Effective April 1, 2022

NAME _____ PHONE NO _____

ACCOUNT # _____ P.O. _____ DATE _____

ADDRESS _____

STOCK WHEELCHAIR MODELS - Silver Vein Finish, Midnight Blue Vinyl Upholstery

<input type="checkbox"/> TRSX5RC6P	16"x16" Recliner, Fixed Height Conventional Desk Arms	\$1,379
	TRSX5RC/WD66/28/R/AT903/U2222C/U550	
<input type="checkbox"/> TRSX5RC8P	18"x16" Recliner, Fixed Height Conventional Desk Arms	\$1,379
	TRSX5RC/WD86/28/R/AT903/U2222C/U550	

STOCK FRONT RIGGINGS - Black Finish

<input type="checkbox"/> T93HAP	Hemi Swingaway Footrests Aluminum Footplates	\$159
<input type="checkbox"/> T93HCP	Hemi Swingaway Footrests Composite Footplates	\$159
<input type="checkbox"/> T94HAP	Hemi Swingaway Elevating Legrests Aluminum Footplates	\$309
<input type="checkbox"/> T94HCP	Hemi Swingaway Elevating Legrests Composite Footplates	\$309
<input type="checkbox"/> AHL4AP	Hemi Smartleg Articulating Legrest Aluminum Footplates	\$399
<input type="checkbox"/> T94HEP	Hemi Legrests, Composite Footplate, Black Plastic Calf Pad	\$239
<input type="checkbox"/> T93HEP	Hemi Footrests without Heel Loops, Composite Footplates	\$119

Specifications

Seat Widths	14", 16", 18", 20" or 22"
Overall Width	22 1/2", 24 1/2", 26 1/2", 28 1/2", 30 1/2"
Seat Depth	16" and 18"
Overall Depth	34" (18"x16" without front riggings)
Back Heights	23" (34" with headrest)
Seat-to-Floor	19 1/2"
Weight Limit	250 lbs. (300 lbs. on 20" or 22" widths)
Chair Weight (approx.)	53 lbs. (18"x16" without front riggings)
Shipping Weight (approx.)	61 lbs. (18"x16" without front riggings)

Please Note: Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

FRAME & SEAT

CHAIR TYPE

<input type="checkbox"/> TRSX5RC	Tracer SX5	\$1,296
----------------------------------	------------	---------

FRAME WIDTH & DEPTH

<input type="checkbox"/> WD46	14" Wide X 16" Deep	No Charge
<input type="checkbox"/> WD48	14" Wide X 18" Deep	\$145
<input type="checkbox"/> WD66	16" Wide X 16" Deep	No Charge
<input type="checkbox"/> WD68	16" Wide X 18" Deep	\$145
<input type="checkbox"/> WD86	18" Wide X 16" Deep	No Charge
<input type="checkbox"/> WD88	18" Wide X 18" Deep	\$145
<input type="checkbox"/> WD06	20" Wide X 16" Deep	No Charge
<input type="checkbox"/> WD08	20" Wide X 18" Deep	\$145
<input type="checkbox"/> WD26	22" Wide X 16" Deep	\$145
<input type="checkbox"/> WD28	22" Wide X 18" Deep	\$271

Standard on SX5 Recliner

- Rear anti-tippers
- Patient only operated push-to-lock wheel lock

Weight Limitation 250 lbs., 300 lbs. 20" & 22" width

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

ARM & BACK

2 ARM TYPE - Removable, Conventional, Non Flip-Back

<input type="checkbox"/> 18	Full Length Fixed Height 10"	\$86
<input type="checkbox"/> 19	Full Length Adj. Height 10" - 14"	\$199
<input type="checkbox"/> 28	Desk Length Fixed Height 10"	\$86
<input type="checkbox"/> 29	Desk Length Adj. Height 10" - 14"	\$199

3 BACK TYPE

<input type="checkbox"/> R	Recliner Frame	\$556
----------------------------	----------------	-------

COLOURS

UPHOLSTERY COLOURS

<input type="checkbox"/> U240	Black Nylon Upholstery	No Charge
<input type="checkbox"/> U550	Midnight Blue Vinyl Upholstery	Standard
<input type="checkbox"/> U67	Black Vinyl	No Charge
<input type="checkbox"/> U86	Ruby Red Vinyl	No Charge

CASTERS

CASTER OPTIONS

<input type="checkbox"/> AT903	8" X 1 1/4" Composite Urethane	
<input type="checkbox"/> 1228	8" X 1 3/4" Semi Pneumatic Casters Lt Grey	
<input type="checkbox"/> 1252	8" X 1 1/4" Pneumatic Casters	
<input type="checkbox"/> 4252	8" x 1 1/4" Pneumatic. Casters w/FF Inserts	

Seat-to-Floor 8" Casters			
Seat-to-Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
17 1/2"	Top	20"	Recliner Axle Mounting Bracket
18 1/3"	Middle	22"	Recliner Axle Mounting Bracket
19 1/2"	Bottom	24"	Recliner Axle Mounting Bracket

FOOTNOTES

- Standard**
- 1 20" and 22" widths have 300 lbs. weight limit requiring 1228 casters and 24" rear wheels.
 - 2 Removable, conventional arm. Non flip-back.
 - 3 90° to 180° of recline. 24" back height with 10" headrest extension.

POSITIONING STRAPS

SEAT POSITIONING STRAPS

<input type="checkbox"/> 1515	Auto Style Buckle Seat Positioning Strap	\$84
<input type="checkbox"/> 1321BK	Hook & Loop Seat Positioning Strap	\$66

WHEELS

WHEEL OPTIONS

<input type="checkbox"/> U1763BK	20" Composite Urethane Tires	No Charge
<input type="checkbox"/> U1764BK	22" Composite Urethane Tires	No Charge
<input type="checkbox"/> 2244BK	24" Composite Pneumatic Tires	\$231
<input type="checkbox"/> 4244BK	24" Composite Pneumatic w/FF Inserts	\$353
<input type="checkbox"/> U2244BK	24" Composite Urethane Treaded Tires	No Charge
4 <input type="checkbox"/> U2222C	24" Composite Urethane Corded Tires	Standard

HANDRIM OPTIONS

5 <input type="checkbox"/> COM	Composite Handrims	Standard
6 <input type="checkbox"/> AHDRM	Aluminum Handrims	No Charge

WHEEL LOCK OPTIONS

<input type="checkbox"/> 9690B	Pull to Lock Wheel Lock (Bolt Thru) - Std Adult	No Charge
<input type="checkbox"/> 9691B	Push to Lock Wheel Lock (Bolt Thru) - Std Hemi	No Charge

ACCESSORIES

BACK ACCESSORIES

<input type="checkbox"/> L8601	Omit Back Upholstery	Less net (\$20)
7 <input type="checkbox"/> STC	Stencil Back Upholstery	\$59

ARM ACCESSORIES

<input type="checkbox"/> 1316	Black Plastic Arm Pads	No Charge
-------------------------------	------------------------	-----------

SEAT ACCESSORIES

<input type="checkbox"/> L8651	Omit Seat Upholstery	Less net (\$15)
--------------------------------	----------------------	-----------------

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/> 1329BK	Toggle Wheel Lock Extension (8")	\$61
8 <input type="checkbox"/> MUPH	Matching Upholstery	No Charge
<input type="checkbox"/> 1322SXR	IV Rod & Bracket, Telescopic, Recliner	\$178

FRONT RIGGING OPTIONS

FRONT RIGGING

9 <input type="checkbox"/> T93HC	Hemi Footrests	\$159
9 <input type="checkbox"/> T93HA	Hemi Footrests Aluminum Footplates	\$159
9,10 <input type="checkbox"/> T94HC	Hemi Elevating Legrests	\$302
9,10 <input type="checkbox"/> T94HA	Hemi Elevating Legrests Aluminum Footplates	\$309
9,10 <input type="checkbox"/> AHL4A	Hemi Smartleg Articulating Legrests Aluminum	\$399
9,10 <input type="checkbox"/> AHL4C	Hemi Smart Articulating Legrests Composite Footplates	\$392
10,11 <input type="checkbox"/> AHS4A	Hemi Smartleg Articulating Legrests-14" Wide Chairs	\$399
13 <input type="checkbox"/> T93HS	Hemi Swingaway Footrests Aluminum Footplates 14" wide	\$158
12,13 <input type="checkbox"/> T94HS	Hemi Swingaway Elevating Legrests Aluminum Footplates 14" Wide	\$302
9 <input type="checkbox"/> T93HE	Hemi Footrests, w/o Heel Loop, Composite Footplates	\$119
9 <input type="checkbox"/> T94HE	Hemi Artic. Legrests, Composite Footplates Black Plastic Calf Pad	\$239
9 <input type="checkbox"/> T93HAB	Hemi Bolt-Thru First, Aluminum Footplates	\$158
9 <input type="checkbox"/> T93HCB	Hemi Bolt-Thru First, Composite Footplates	\$158
9 <input type="checkbox"/> T93HSB	Hemi Bolt-Thru First, Composite Footplates 14" wide	\$158
9 <input type="checkbox"/> T94HAB	Hemi Bolt-Thru Legrest, Aluminum Footplates	\$302
9 <input type="checkbox"/> T94HCB	Hemi Bolt-Thru Legrest, Composite Footplates	\$302
9 <input type="checkbox"/> T94HSB	Hemi Bolt-Thru Legrest, Aluminum Footplate 14" wide	\$302

FOOTPLATE TYPE

9 <input type="checkbox"/> 1350	Extra Large Footplates	\$106
<input type="checkbox"/> AT5543	Adjustable Angle Flip-Up Footplates	\$132

FOOTNOTES

- 4 Available with 20" and 22" wide chairs only.
- 5 Available with 24" rear wheels only.
- 6 Aluminum handrims not available on U2222C wheels.
- 7 3/4" block letters, 12 or more chairs, no charge.
- 8 Applies to arm pads and calf pads.
- 9 Not available on 14" width.
- 10 All calf pads are black
- 11 Only available with 14" width.
- 12 Standard on swingaway footrests T93HA and T93HC
- 13 Not available with 14" and 22" widths.

FOOTREST/LEGREST ACCESSORIES

<input type="checkbox"/> 1820	Non-Removable Device	\$41
12 <input type="checkbox"/> 1346	Heel Loops	\$63
<input type="checkbox"/> L8515	Omit Calf Pads	Less net (\$10)

STOCK OPTIONS

13 <input type="checkbox"/> 1514	Overhead Anti-Theft Device	\$136
----------------------------------	----------------------------	-------

Invacare Canada reserves the right to change prices and specifications without notice.

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324 (Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)

Invacare®, and Yes, you can., are trademarks of Invacare Corporation.

WWW.INVACARE.CA