





**alber® Power Mobility Products**  
**Suggested Canadian Price List**  
**Effective April 27th, 2020**

Item Number	Description	MSRP
<b>e-motion M15 - Power Assist Wheels for Manual Wheelchairs</b>		
	<p><i>e-motion combines active mobility with therapeutic benefit. The electric motors integrated in the wheel hubs effectively support the propelling movement of the wheelchair driver. The user remains active and trains his or her vital functions while muscles and joints are relieved of stress. Very little effort is required to manoeuvre, propel and brake on various terrains as e-motions are perfectly adjustable to the needs of the user: all driving parameters are freely programmable! <b>ADP Approved! ADP # WPS1DO005</b></i></p> <p><b>Kit Includes:</b> 2 programmable drive wheels available in 22"/24" diameter, programmable remote control, battery pack of 25.2V/6 AH, lithium-ions, leak-proof and maintenance free battery, automatic quick charger and brackets for mounting to wheelchair.</p> <p><b>In Photo:</b> e-motion with optional anti-tippers with jack-up function.</p>	
1490998-22	22" e-motion M15 Power Assist Wheels with ECS (Ergonomic Control System)	\$8,515
1490998-24	24" e-motion M15 Power Assist Wheels with ECS (Ergonomic Control System)	\$8,515
<b>Mounting brackets included with e-motion unit purchase.</b>		No Charge
<b>Quotes Required: Please supply manufacturer, chair model, width, depth, seat to floor height and camber.</b>		
<b>e-Motion M15 Accessories, Modifications, Replacement Service Parts</b>		
<b>Mounting Bracket</b>	See e-Motion approved bracket listing for North American chairs on Chart A, page 3.	See Chart A
1489214	Anti-tippers with jack-up function (pair) - <i>Mounting hardware is not included. Please call for quote.</i>	\$650
1456301-1	Push rim cover e-Motion, 22" (each) (gummy/rubber style)	\$134
1456273	Push rim cover e-Motion, 24" (each) (gummy/rubber style)	\$134
1461521	Rubber coated push rim e-motion 24" (each)	\$175
1464497	24" x 1 Poly Urethane Tire Grey (each) (solid tire, comes in 24" only)	\$185
1490975	Spoke guard 22" (each)	\$178
1490973-1	Spoke guard 24" (each)	\$178
1490821	Bracket for remote control (each)	\$364
1490492-1	Spare battery, 6.0 AH Lithium-ions (each)	\$700
1456265	Axle Conversion for manual wheels compatible to Alber products (each)	\$396
1467780	Brake widening set (pair) (required, depending on the wheelchair type, to adapt the parking brake to the e-Motion wheels)	\$56
1490707	Replacement Programmable Remote control	\$1,039
<b>viamobil V25 - Attendant Push and Brake Aid</b>		
	<p><i>The viamobil V25 is the new high-performance pushing and braking aid for wheelchairs. It's extremely high versatility (it will fit to practically any wheelchair) and it's power (high user weight capacity and excellent climbing ability) is setting standards. The utmost attention was paid to the safety and user-friendliness of the viamobil to reduce caregiver effort.</i></p> <p><b>Kit Includes:</b> Drive unit with hub drive motor, battery-pack with automatic charger, control unit and brackets for mounting to wheelchair.</p> <p><b>In Photo:</b> viamobil with optional anti-tippers and additional push handles.</p>	
V25	alber viamobil V25	\$5,355
<b>Mounting brackets included with viamobil unit purchase.</b>		No Charge
<b>Quotes Required: Please supply manufacturer, chair model, width, depth, seat to floor height and camber.</b>		
<b>viamobil Accessories, Replacement Service parts</b>		
1456959	Additional push handles for V15 model (pair)	\$392
1521030	Additional push handles for V25 model (pair)	\$392
1457357	Anti-tippers wide complete exceeding 16.4W (pair) for V15	\$514
1457354	Anti-tippers narrow complete up to 16.4W (pair) for V15	\$514
1564967	Battery Pack for Viamobil V25, Complete	\$760

**scalamobil S35/38 - Mobile Wheelchair Stairclimber Device**



*Keeping safety and comfort for the wheelchair user and the attendant in mind, the scalamobil will work on all types of staircases: wide, narrow, tall, short and spiral. The stairs are optimally protected thanks to the patented climbing mechanism. Two different models available depending on step height requirements.*

*Caregiver effort is reduced thanks to the new ErgoBalance handles and the ability to individually adjust the width and height of the handlebars. It is compatible with most wheelchair models thanks to special brackets and can be docked onto the wheelchair in a matter of seconds without any tools being required.*

**Kit Includes:** Extendable handlebar, electronic speed control, single-step or multi-step mode, SoftStep function for a smooth ride, battery, automatic battery charger, and brackets for mounting to wheelchair.

<b>S35</b>	scalamobil S35 - Mobile stair climber (For up to 7.87" high steps)	\$8,529
<b>S38EXT</b>	scalamobil S38 - Mobile stair climber with extended stair height (For up to 9.84" high steps)	\$9,179

**Mounting brackets included with scalamobil unit purchase.**

No Charge

**Quotes Required: Please supply manufacturer, chair model, width, depth, seat to floor height and camber.**

**scalamobil Accessories, Replacement Service parts**

<b>1457162*</b>	Anti-tippers wide version, 15 mm (pair)	\$576
<b>1457163*</b>	Anti-tippers narrow version, 12 mm (pair)	\$576

**scalaport**

*A portable stair climbing system that helps you transport most manual wheelchairs on stairs of all kinds.*

<b>X6S25</b>	scalaport X6, universal rack with headrest	<b>For S25 model</b>	\$1,689
<b>X6S30*</b>	scalaport X6, universal rack with headrest adaptor	<b>For S25 and S30 models</b>	\$1,809
<b>X7S30</b>	scalaport X7 for scalamobil S30	<b>For S30 model</b>	\$1,819
<b>X7S35</b>	scalaport X7 for scalamobil S35/S38	<b>For S35/S38 model</b>	\$2,090
<b>1490239</b>	Headrest with clips (new version)	<b>For S25/S27/S28/S30/S31/X7 models</b>	\$308

**scalachair X3**

*A good alternative to a normal wheelchair, particularly in small rooms or where the staircase is narrow. The bracket for the scalamobil is integrated into the scalachair X3. The chair can be easily detached for use as a indoor push chair without the need of additional hardware.*

<b>X3FOLD35</b>	Scalachair X3, Folding	<b>For S35 model</b>	\$2,075
<b>X3RIGID35</b>	Scalachair X3, Rigid	<b>For S35 model</b>	\$1,969
<b>X3FOLD38</b>	Scalachair X3, Folding	<b>For S38 model</b>	\$2,075
<b>X3RIGID38</b>	Scalachair X3, Rigid	<b>For S38 model</b>	\$1,969

Product	Product Weight	User Weight	Total Weight Capacity Includes user, wheelchair and wheelchair add-on drive
e-motion M15	46.2 lbs	286 lbs	396 lbs
viamobil V25	29.7 lbs	352 lbs	462 lbs
scalamobil S35/S38	55.1 lbs	308 lbs	396 lbs
e-fix 25	58.3 lbs	264 lbs	374 lbs
e-fix 26	71.5 lbs	352 lbs	462 lbs

**Invacare Canada reserves the right to change prices and specifications without notice.**

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)  
Invacare®, and Yes, you can., are trademarks of Invacare Corporation.  
WWW.INVACARE.CA